



Patient Last Name: _____
Fitter Last Name: _____
Fitter Title: _____
Date: _____

Patient First Name: _____
Fitter First Name: _____
(example: PT/OT/PTA)




Dayspring™

Order form (1/2)

General Measurements

Referring to the diagram provided at right, complete the corresponding fields below.
All measurements should be taken in centimeters (cm) to the nearest half unit.

MEASUREMENT	(A)	DISTANCE FROM WRIST TO ELBOW:		cm
MEASUREMENT	(B)	DISTANCE FROM WRIST TO AXILLA:		cm
MEASUREMENT	(C)	DISTANCE FROM WRIST TO KNUCKLES:		cm
MEASUREMENT	(D)	CIRCUMFERENCE OF KNUCKLES:		cm



Proceed to Page 2.



Patient Last Name: _____

Fitter Last Name: _____

Fitter Title: _____

Date: _____

Patient First Name: _____

Fitter First Name: _____

(example: PT/OT/PTA)



Dayspring™

Order form (2/2)

Circumference Measurements

The size of the garment depends on the length of the patient’s arm (B), measured previously on page one. Mark and complete the corresponding section below.

Short:
(B) is below 44cm

DISTANCE FROM WRIST

0 — 4 — 11 — 18 — 25 — 31 — 38

cm cm cm cm cm cm cm

CIRCUMFERENCE

cm

cm

cm

cm

cm

cm

cm

WRIST CREASE

SECTION 4S

Medium:
(B) is between 44cm and 48cm

DISTANCE FROM WRIST

0 — 4 — 12 — 20 — 27 — 35 — 43

cm cm cm cm cm cm cm

CIRCUMFERENCE

cm

cm

cm

cm

cm

cm

cm

WRIST CREASE

SECTION 4M

Long:
(B) is above 48cm

DISTANCE FROM WRIST

0 — 4 — 12 — 19 — 26 — 33 — 41 — 48

cm cm cm cm cm cm cm

CIRCUMFERENCE

cm

cm

cm

cm

cm

cm

cm

cm

WRIST CREASE

SECTION 4L