




Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____

lymphedivas[®]

medically correct *fashion* for lymphedema



	MEASUREMENTS	SMALL	MEDIUM	LARGE
FINGER		4.4 – 6 cm	5.5 – 7.5 cm	6.6 – 9 cm
PALM		15 – 18.7 cm	19 – 21.9 cm	22 – 26 cm
WRIST		14.6 – 19 cm	16.5 – 21.6 cm	19 – 23.5 cm
FOREARM		17.8 – 23.2 cm	23.5 – 28.7 cm	28.9 – 34.3 cm
BICEP		21.6 – 39.4 cm	27 – 45 cm	32.4 – 50.5 cm
LENGTH		SHORT 37 - 44 cm	LONG 43 - 56 cm	

SLEEVE	SIZE	CLASS	LENGTH
	<input type="checkbox"/> SMALL	<input type="checkbox"/> CLASS 1 (20–30 mmHg)	<input type="checkbox"/> SHORT Fits arms 37–44 cm
	<input type="checkbox"/> MEDIUM	<input type="checkbox"/> CLASS 2 (30–40 mmHg)	<input type="checkbox"/> LONG Fits arms 43–56 cm
	<input type="checkbox"/> LARGE		

GAUNTLET OR GLOVE	GARMENT	HAND	SIZE	CLASS
	<input type="checkbox"/> GAUNTLET	<input type="checkbox"/> LEFT	<input type="checkbox"/> SMALL	<input type="checkbox"/> CLASS 1 (20–30 mmHg)
	<input type="checkbox"/> GLOVE	<input type="checkbox"/> RIGHT	<input type="checkbox"/> MEDIUM	<input type="checkbox"/> CLASS 2 (30–40 mmHg)
			<input type="checkbox"/> LARGE	

USE WITH LYMPHEDIVAS ONLY