



Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____



MEDI NON-CUSTOM FLAT KNIT HAND GLOVE MEASUREMENT FORM

Measure in centimeters

Left Right

Measurement Points (Two measurements are required)

Circumference **point 'a'** is found across the metacarpus (hand).

Circumference **point 'c'** is found distally on the arm where the wrist meets the hand. One should find the smallest circumference at this point (wrist).


Circumference 'a' = _____

Circumference 'c' = _____

Item	Color(s) Available
MEDI ESPRIT GLOVE ONLY Class II Only	MEDI ESPRIT GLOVE - "CARAMEL ONLY"

* MEDI ESPRIT GLOVE has full compression in the fingers; with seams

SIZE CHART



Flat Knit - Medi Esprit

	I	II	III	IV	V
A	17.5-18.5 cm 6¾"- 7¼"	18.5 - 20 cm 7¼"- 7¾"	20-22 cm 7¾"- 8¾"	22-24 cm 8¾"- 9½"	24- 27cm 9½"-10¾"
C	15-17 cm 5¾"- 6¾"	17-19 cm 6¾"- 7½"	19-21 cm 7½"- 8¼"	21-23 cm. 8¼"- 9¼"	23-25 cm 9¼"- 10¼"