



Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____



MEDI NON-CUSTOM CIRCULAR KNIT HAND GLOVE & GAUNTLET MEASUREMENT FORM

Measure in centimeters

Left Right

Measurement Points (Two measurements are required)

Circumference **point 'a'** is found across the metacarpus (hand).

Circumference **point 'c'** is found distally on the arm where the wrist meets the hand. One should find the smallest circumference at this point (wrist).

Circumference 'a' = _____

Circumference 'c' = _____

Item	Class I	Class II	Color(s) Available
HARMONY GAUNTLET	20-30 mmHg <input type="checkbox"/>	30-40 mmHg <input type="checkbox"/>	HARMONY GAUNTLET Sand <input type="checkbox"/> Caramel <input type="checkbox"/> Black <input type="checkbox"/>
HARMONY GLOVE	20-30 mmHg <input type="checkbox"/>	30-40 mmHg <input type="checkbox"/>	HARMONY GLOVE Sand <input type="checkbox"/> Caramel <input type="checkbox"/> Black <input type="checkbox"/>
SEAMLESS HARMONY GAUNTLET	20-30 mmHg <input type="checkbox"/>	30-40 mmHg <input type="checkbox"/>	SEAMLESS HARMONY GAUNTLET Sand <input type="checkbox"/> Caramel <input type="checkbox"/> Black <input type="checkbox"/>
SEAMLESS HARMONY GLOVE	20-30 mmHg <input type="checkbox"/>	30-40 mmHg <input type="checkbox"/>	SEAMLESS HARMONY GLOVE Sand <input type="checkbox"/> Caramel <input type="checkbox"/> Black <input type="checkbox"/>

* MEDI HARMONY GLOVE has 50% compression in the fingers; no seams

lymphology measuring charts

mediven® harmony sleeves, sleeve with gauntlet, gauntlet and gloves

circumference (cm)	1	2	3	4	5	6	7	8
"G" extra-wide	28-31	30-34	32.5-36.5	34.5-39.5	36.5-41.5	39-44	41-46.5	44.5-50.5
"G" standard	25-28	27.5-31	29.5-33.5	31.5-36	33-37.5	35.5-40.5	37.5-42.5	40-45.5
"E" extra-wide	24-26.5	26-29.5	28-31.5	29.5-34	31-35.5	33-37.5	35-40	38-43
"E" standard	23-25	24-27	25.5-29	27-31	28.5-32.5	30-34.5	32-36.5	34.5-39.5
"C1"	15-16.5	16-18.5	17.5-20	19-22	20.5-24	22.5-26	24.5-28.5	27-31
"C"	13.5-14.5	14-16	15.5-18	17-20	18.5-21.5	20-23.5	22-25.5	24.5-28.5
"A"	15-18	15.5-19.5	16.5-20.5	18.5-23	19.5-24	21.5-26	23.5-28	25.5-30

gloves & gauntlets available in sizes 2-7

armsleeves available in sizes 1-8