



Account Number: 120167

Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example PT/OT/PTA)
 Date: _____

**MEDI NON-CUSTOM CIRCULAR KNIT SLEEVE & GAUNTLET
 COMBO MEASUREMENT FORM**

Measure in centimeters

Left Right

Length Measurement from axilla to wrist (c-g) = _____



Measurement Points

Circumference **point 'g'** is taken around the upper arm at the axilla. For ease in determining this point, place a folded piece of paper up into the arm pit area. Point 'g' on the arm will be even with the top of the paper. **Measure with moderate tension.**

Circumference **point 'e'** is found in the middle of the elbow. **Measure without tension.**

Circumference **point 'c'** is the point of greatest compression and therefore the most important point. This point is found distally on the arm where the wrist meets the hand. One should find the smallest circumference at this point. **Measure without tension.**

Circumference 'g' = _____ **SILICONE BORDER ONLY**

Circumference 'e' = _____


Circumference 'c1' = _____

Circumference 'c' = _____

Circumference 'a' = _____

Item	Class I	Class II	Colors Available	Note
HARMONY SLEEVE with ATTACHED GAUNTLET ONLY	20-30 mmHg <input type="checkbox"/>	30-40 mmHg <input type="checkbox"/>	HARMONY SLEEVE w/GAUNTLET Sand <input type="checkbox"/> Caramel <input type="checkbox"/> Black <input type="checkbox"/>	SILICONE BORDER ONLY * MEDI HARMONY COMBO DOES NOT COME WITH GLOVE; GAUNTLET ONLY

SIZE CHART



lymphology measuring charts

mediven® harmony sleeves, sleeve with gauntlet, gauntlet and gloves

circumference (cm)	1	2	3	4	5	6	7	8
"G" extra-wide	28-31	30-34	32.5-36.5	34.5-39.5	36.5-41.5	39-44	41-46.5	44.5-50.5
"G" standard	25-28	27.5-31	29.5-33.5	31.5-36	33-37.5	35.5-40.5	37.5-42.5	40-45.5
"E" extra-wide	24-26.5	26-29.5	28-31.5	29.5-34	31-35.5	33-37.5	35-40	38-43
"E" standard	23-25	24-27	25.5-29	27-31	28.5-32.5	30-34.5	32-36.5	34.5-39.5
"C1"	15-16.5	16-18.5	17.5-20	19-22	20.5-24	22.5-26	24.5-28.5	27-31
"C"	13.5-14.5	14-16	15.5-18	17-20	18.5-21.5	20-23.5	22-25.5	24.5-28.5
"A"	15-18	15.5-19.5	16.5-20.5	18.5-23	19.5-24	21.5-26	23.5-28	25.5-30

gloves & gauntlets available in sizes 2-7

armsleeves available in sizes 1-8