



Account Number: 120167

Patient Last Name: _____ Patient First Name: _____
Fitter Last Name: _____ Fitter First Name: _____
Fitter Title: _____ (example PT/OT/PTA)
Date: _____



MEDI NON-CUSTOM FLAT KNIT HAND GLOVE MEASUREMENT FORM

Measure in centimeters

Left [] Right []

Measurement Points (Two measurements are required)

Circumference point 'a' is found across the metacarpus (hand).

Circumference point 'c' is found distally on the arm where the wrist meets the hand. One should find the smallest circumference at this point (wrist).

Circumference 'a' = _____


Circumference 'c' = _____

Table with 2 columns: Item, Color(s) Available. Item: MEDI ESPRIT GLOVE ONLY Class II Only. Color(s) Available: MEDI ESPRIT GLOVE - "CARAMEL ONLY"

* MEDI ESPRIT GLOVE has full compression in the fingers; with seams

SIZE CHART

Flat Knit - Medi Esprit



Size chart table with columns I, II, III, IV, V and rows A, C. Values include cm and inch measurements.