



Account Number: 120167

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Fitter Title: \_\_\_\_\_ (example PT/OT/PTA)

Date: \_\_\_\_\_

**MEDI NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASUREMENT FORM AND SIZING CHART**



Measure in Centimeters Left  Right

Thigh Length \_\_\_\_\_

Calf Length \_\_\_\_\_

A-g STANDARD LENGTH IS > 71 CM  
A-g PETITE LENGTH < 71 CM

A-d STANDARD LENGTH IS > 38 CM  
A-d PETITE LENGTH IS < 38 CM

**For thigh length**, the length measurement is taken from the floor (inside of heel) along the inside of the leg to the inseam/groin, point 'g'; this length should end at the gluteal fold.


**For calf length**, the length measurement is taken from the floor (inside of heel) along the inside of the leg to approximately 3 1/2 centimeters below the crease of the knee, point 'd'

**Measurement Points**

Circumference **point 'g'** is found at the inseam/groin

Circumference **point 'c'** is the largest circumference of the calf

Circumference **point 'b'** is the point of greatest compression and therefore the most important point. This point is the smallest circumference of the ankle found closest to and above the ankle bone



**Circumference 'g'** = \_\_\_\_\_

**Measurement Point for Waist Attachment**

The following measurement will be necessary if you are ordering a thigh length stocking to be sewn to a waist attachment

**Circumference 'c'** = \_\_\_\_\_

Hip Circumference 't' = \_\_\_\_\_

Waist attachment option not available in Medi Forte stockings.

**Circumference 'b'** = \_\_\_\_\_

| Item       | Class I                             | Class II                            | Class III                           | Options**  |
|------------|-------------------------------------|-------------------------------------|-------------------------------------|--|
| Medi Plus  | 20-30 mmHg <input type="checkbox"/> | 30-40 mmHg <input type="checkbox"/> | 40-50 mmHg <input type="checkbox"/> | Silicone <input type="checkbox"/> No Silicone <input type="checkbox"/> Open Toe <input type="checkbox"/><br>Closed Toe* <input type="checkbox"/> |
| Medi Forte | N/A                                 | 30-40 mmHg <input type="checkbox"/> | 40-50 mmHg <input type="checkbox"/> | Silicone <input type="checkbox"/> No Silicone <input type="checkbox"/> Open Toe <input type="checkbox"/><br>N/A                                  |

\*Closed toe option not available in Medi Plus, Class III or Medi Forte stockings.

\*\*Medi Plus and Medi Forte stockings only come in beige.

|  | I                          | II                        | III                         | IV                         | V                             | VI                         | VII                             |
|--|----------------------------|---------------------------|-----------------------------|----------------------------|-------------------------------|----------------------------|---------------------------------|
| <b>g</b><br>panty & thigh  | 38-57 cm<br>15"-22½"       | 41-62 cm<br>16"-24¼"      | 44-67 cm.<br>17¼"-26¼"      | 47-72 cm.<br>18½"-<br>28¼" | 50-77 cm.<br>19¾"-30¼"        | 53-82 cm.<br>21"-32¼"      | 56-87 cm.<br>22¼"-<br>34¼"      |
| <b>g</b><br>thigh w/ silicone<br>only                            | 42-49 cm<br>16½"-19¼"      | 45-53 cm<br>17¾"-<br>20¾" | 48-57 cm.<br>19"-22¼"       | 52-61 cm.<br>20½"- 24"     | 56-65 cm.<br>22"- 25¾"        | 60-70 cm.<br>23½"-<br>27½" | 64-74 cm.<br>25"- 29"           |
| <b>c</b><br>calf & calf w/<br>silicone                           | 26-34 cm<br>10¼"-13¼"      | 28-37 cm<br>11"-14½"      | 30-40 cm.<br>11¾"-15¾"      | 32-43 cm.<br>12½"- 17"     | 34-46 cm.<br>13¼"-18"         | 36-49 cm.<br>14¼"-<br>19¼" | 38-52 cm.<br>15"- 20½"          |
| <b>c</b><br>*calf extra wide &<br>calf extra wide w/<br>silicone | 34-38 cm<br>13¼"-15"       | 37-42 cm<br>14½"-<br>16½" | 40-46 cm.<br>15¾"-18"       | 43-50 cm.<br>17"- 19½"     | 46-54 cm.<br>18"-21"          | 49-58 cm.<br>19¼"-<br>22¾" | 52-62 cm.<br>20½"-<br>24½"      |
| <b>b</b><br>ankle  | 16.5-18.5<br>cm<br>6½"-7¼" | 19-21 cm<br>7½"-8¼"       | 21.5-23.5<br>cm.<br>8½"-9¼" | 24-26 cm.<br>9½"- 10¼"     | 26.5-28.5<br>cm.<br>10½"-11¼" | 29-32 cm.<br>11½"-<br>12½" | 32.5-35<br>cm.<br>12¾"-<br>13¾" |

\*Medi Plus extra wide calf-length stockings are **only available in classes I and II.**

\*\*Medi Forte extra wide calf-length stockings are **only available in classes II and III.**

Luna Medical, Inc. · Specialists in Venous & Lymphatic Insufficiencies

1816 W. Belmont Avenue · Suite 1 · Chicago, IL 60657 · Phone (800) 380-4339 · Fax (888) 696-0299 · www.lunamedical.com