




SIENNAWRAP FOOT SIZE/MEASURING CHART

Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____



SIENNAWRAP FOOT				
	Small	Medium	Large	Extra Large
B	17-24cm	21-28cm	26-33cm	31-38cm
A	19-23cm	22-27cm	26-31cm	30-35cm
Lengths	Small	Medium	Large	Extra Large
Y	17cm	18cm	19cm	20cm

Measure on the medial side of the foot from the back of the heel to past the bony prominence of the metacarpal. Must be a flat measurement on the ground.