



**Patient Last Name:** \_\_\_\_\_ **Patient First Name:** \_\_\_\_\_  
**Fitter Last Name:** \_\_\_\_\_ **Fitter First Name:** \_\_\_\_\_  
**Fitter Title:** \_\_\_\_\_ **(example: PT/OT/PTA)**  
**Date:** \_\_\_\_\_

# sigvaris

## Chipflow Calf & Foot

### Product Information

Product includes one custom Chipflow Calf, one Compreboot Standard Foot wrap, one Chip Foam Liner, and one pair of non-compressive Basic Liners\*.

<input type="checkbox"/> Right Calf	Color: Black	Quantity:
<input type="checkbox"/> Right Foot	Color: Black:	Quantity:
<input type="checkbox"/> Left Calf	Color: Black	Quantity:
<input type="checkbox"/> Left Foot	Color: Black:	Quantity:

\*For the largest size, a pair of non-compressive Cotton Liners is included instead of the Basic Liners.

### Important

Exact measurements are critical for this garment to ensure a proper fit. If you'd like to learn more about measuring and fitting, attend our MCE Basic Fitter Training. To sign up, visit [sigvariseducation.com](http://sigvariseducation.com) or contact your local territory manager for more information.


Alternatively, call Customer Care Solution Center at **800-322-7744**, or e-mail [us\\_orders@sigvaris.com](mailto:us_orders@sigvaris.com), to receive a remote consultation/training.

### Supplies Needed

- Measuring instructions and forms.
- SIGVARIS GROUP Measuring tape and body pen (or eyeliner pencil).

### Circumference

	Left	Right
C1	_____	_____
C	_____	_____
B	_____	_____
A	_____	_____
G	_____	_____
I	_____	_____
J	_____	_____
K	_____	_____



**K:** Measure **medial** length from heel to 1st metatarsal head

### Chipflow Calf Size Chart (A-C1) (Available in Black Only)

**Note: Product contains Measurement Limitations with Min & Max sizing for ankle and calf**

	Min	Max
C1	31cm	69cm
C	31cm	69cm
B	26cm	59cm
A	19cm	51cm
<b>Black</b>		
G Length	27cm	42cm

#### Measurement Disclaimer

Max calf < 2x ankle