




Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____

sigvaris

Chipvest Custom



Product Information

Product includes one Chipvest Custom.

Custom Chipvest Item #: **2239-VS**

- Full (Bilateral)
- Right Side (Unilateral)
- Left Side (Unilateral)

Measuring Instructions

Step One:

Measure the length in centimeters from the suprasternal notch to the waist.

Step Two:

Measure circumferences in centimeters.

Custom Size Measurements

Length	Centimeters
L	
Circumference	Centimeters
Axilla	
Largest Chest	
Xyphoid Process	
Waist	
Mid-Hip	

Caution: Drawstring contains natural rubber latex.

Important

Exact measurements are critical for this garment to ensure a proper fit. If you'd like to learn more about measuring and fitting, attend our MCE Basic Fitter Training. To sign up, visit sigvariseducation.com or contact your local territory manager for more information.

Alternatively, call Customer Care Solution Center at **800-322-7744**, or e-mail us_orders@sigvaris.com, to receive a remote consultation/training.

Supplies Needed

- Cell phone with camera. Photos of the torso with measurement markings must be emailed to: us_orders@sigvaris.com
- Measuring instructions and forms.
- SIGVARIS GROUP Measuring tape and body pen (or eyeliner pencil).
- Signed Custom Order Terms & Conditions Form (include with order).