



Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example: PT/OT/PTA)  
 Date: \_\_\_\_\_

# sigvaris

## Comprefit Plus Calf & Foot

### Product Information

Product includes one Comprefit Standard Calf wrap, one flat foam liner, one Compreboot Standard Foot wrap, and 1 pair of Basic Liners. The Compreboot Standard Foot wrap included matches in size of the Comprefit Standard Calf wrap selected. **Note:** If Compreboot needed is a different size, select the Comprefit Plus Calf with separate foot piece.

<input type="checkbox"/> Right Calf	Size:	Item #:
Color: Black	Length:	Quantity:
<input type="checkbox"/> Left Calf	Size:	Item #:
Color: Black	Length:	Quantity:

### Companion Products (sold separately)

<input type="checkbox"/> Comprefit Strap Extender	Item #:
Color: <input type="checkbox"/> Black	Quantity:

**Note:** Liner and additional Compreboot options available, each sold separately.

### Compreboot Standard Foot Size Chart (I-K)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm
<b>Black</b>	1801-BTR	1801-BTL	1803-BTR	1803-BTL	1805-BTR	1805-BTL

### Comprefit Strap Extender Chart

	1 Strap	10 Straps
<b>Black</b>	1100-SE1	1100-SEB

### Circumference


Left    Right

C1 \_\_\_\_\_  
 C \_\_\_\_\_  
 B \_\_\_\_\_  
 A \_\_\_\_\_  
 I \_\_\_\_\_  
 J \_\_\_\_\_

### G Length

Regular: 30-35.9cm  
 Tall: 36-42cm

G \_\_\_\_\_  
 K \_\_\_\_\_



**K:** Measure **medial** length from heel to 1st metatarsal head

### Comprefit Standard Calf Regular Size Chart (A-C1)

	Small	Medium	Large	X Large	XX Large
	<b>Regular</b>				
C1	24-34cm	29-39cm	34-44cm	43-53cm	50-60cm
C	24-34cm	29-39cm	34-44cm	43-53cm	50-60cm
B	19-29cm	24-34cm	28-38cm	36-46cm	39-50cm
A	15-24cm	16-25cm	20-31cm	27-37cm	28-38cm
<b>Black</b>	1101-BPR	1102-BPR	1103-BPR	1104-BPR	1105-BPR

### Comprefit Standard Calf Tall Size Chart (A-C1)

	Small Tall	Medium Tall	Large Tall	X Large Tall	XX Large Tall
	<b>Tall</b>				
C1	24-34cm	29-39cm	34-44cm	43-53cm	50-60cm
C	21-31cm	26-36cm	30-40cm	39-49cm	45-55cm
B	16-26cm	20-30cm	25-35cm	31-41cm	35-45cm
A	15-24cm	16-25cm	20-31cm	27-37cm	28-38cm
<b>Black</b>	1101-BPT	1102-BPT	1103-BPT	1104-BPT	1105-BPT