



Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____

sigvaris

Comprefit Standard Calf & Foot

Product Information

Product includes one Comprefit Standard Calf wrap, one Compreboot Standard Foot wrap, and 1 pair of Basic Liners. The Compreboot Standard Foot wrap included matches the size of the Comprefit Standard Calf wrap selected. **Note:** If Compreboot needed is a different size, select the Comprefit Standard Calf with separate foot piece.

Right Calf Size: _____ Item #: _____
 Color: Black Beige Length: _____ Quantity: _____

Left Calf Size: _____ Item #: _____
 Color: Black Beige Length: _____ Quantity: _____

Companion Products (sold separately)

Comprefit Strap Extender Item #: _____
 Color: Black Beige Quantity: _____

Note: Liner and additional Compreboot options available, each sold separately.

Compreboot Standard Foot Size Chart (I-K)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm
Black	1801-BTR	1801-BTL	1803-BTR	1803-BTL	1805-BTR	1805-BTL
Beige	1811-BTR	1811-BTL	1813-BTR	1813-BTL	1815-BTR	1815-BTL

Comprefit Strap Extender Chart

	1 Strap	10 Straps
Black	1100-SE1	1100-SEB
Beige	1110-SE1	1110-SEB


Circumference

Left Right
 C1 _____
 C _____
 B _____
 A _____
 I _____
 J _____

G Length

Regular: 30-35.9cm
 Tall: 36-42cm

G _____
 K _____



K: Measure **medial** length from heel to 1st metatarsal head

Comprefit Standard Calf Regular Size Chart (A-C1)

	Small	Medium	Large	X Large	XX Large
	C1	29-39cm	34-44cm	39-49cm	48-58cm
C	29-39cm	34-44cm	39-49cm	48-58cm	55-65cm
B	24-34cm	29-39cm	33-43cm	41-51cm	44-55cm
A	20-29cm	21-30cm	25-36cm	32-42cm	33-43cm
Black	1101-BKR	1102-BKR	1103-BKR	1104-BKR	1105-BKR
Beige	1111-BKR	1112-BKR	1113-BKR	1114-BKR	1115-BKR

Comprefit Standard Calf Tall Size Chart (A-C1)

	Small Tall	Medium Tall	Large Tall	X Large Tall	XX Large Tall
	C1	29-39cm	34-44cm	39-49cm	48-58cm
C	26-36cm	31-41cm	35-45cm	44-54cm	50-60cm
B	21-31cm	25-35cm	30-40cm	36-46cm	40-50cm
A	20-29cm	21-30cm	25-36cm	32-42cm	33-43cm
Black	1101-BKT	1102-BKT	1103-BKT	1104-BKT	1105-BKT
Beige	1111-BKT	1112-BKT	1113-BKT	1114-BKT	1115-BKT