



Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____

sigvaris

Comprefit Standard Calf

Product Information

Product includes one Comprefit Standard Calf and one pair of non-compressive Basic Liners.

<input type="checkbox"/> Right Calf	Size: _____	Item #: _____
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length: _____	Quantity: _____
<input type="checkbox"/> Left Calf	Size: _____	Item #: _____
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length: _____	Quantity: _____

Companion Products (sold separately)

<input type="checkbox"/> Comprefit Strap Extender	Size: _____	Item #: _____
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Quantity: _____	

Note: Liner and Compreboot options available, each sold separately.

Comprefit Strap Extender Chart

	1 Strap	10 Straps
Black	1100-SE1	1100-SEB
Beige	1110-SE1	1110-SEB


Circumference

Left Right

C1 _____
 C _____
 B _____
 A _____

G Length
 Regular: 30–35.9cm
 Tall: 36–42cm

G _____



Comprefit Standard Calf Regular Size Chart (A–C1)

	Comprefit Standard Calf Regular Size Chart (A–C1)					
	Small	Medium	Large	X Large	XX Large	
Regular	C1	29–39cm	34–44cm	39–49cm	48–58cm	55–65cm
	C	29–39cm	34–44cm	39–49cm	48–58cm	55–65cm
	B	24–34cm	29–39cm	33–43cm	41–51cm	44–55cm
	A	20–29cm	21–30cm	25–36cm	32–42cm	33–43cm
	Black	1101-NF-BKR	1102-NF-BKR	1103-NF-BKR	1104-NF-BKR	1105-NF-BKR
Beige	1111-NF-BKR	1112-NF-BKR	1113-NF-BKR	1114-NF-BKR	1115-NF-BKR	

Comprefit Standard Calf Tall Size Chart (A–C1)

	Comprefit Standard Calf Tall Size Chart (A–C1)					
	Small Tall	Medium Tall	Large Tall	X Large Tall	XX Large Tall	
Tall	C1	29–39cm	34–44cm	39–49cm	48–58cm	55–65cm
	C	26–36cm	31–41cm	35–45cm	44–54cm	50–60cm
	B	21–31cm	25–35cm	30–40cm	36–46cm	40–50cm
	A	20–29cm	21–30cm	25–36cm	32–42cm	33–43cm
	Black	1101-NF-BKT	1102-NF-BKT	1103-NF-BKT	1104-NF-BKT	1105-NF-BKT
Beige	1111-NF-BKT	1112-NF-BKT	1113-NF-BKT	1114-NF-BKT	1115-NF-BKT	