



Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example: PT/OT/PTA)  
 Date: \_\_\_\_\_

# sigvaris

## Compreflex Transition Calf

### Product Information

Product includes one Compreflex Transition Calf wrap and one pair of Transition Liners. Liner included matches the garment size. Example: The medium garment includes a medium liner.

<input type="checkbox"/> Right Calf	Size: _____	Item #: _____
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length: _____	Quantity: _____
<input type="checkbox"/> Left Calf	Size: _____	Item #: _____
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length: _____	Quantity: _____

### Companion Products (sold separately)

<input type="checkbox"/> Compreflex Strap Extender	Item #: _____
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Quantity: _____


**Note:** Liner and Compreboot options available, each sold separately.

### Circumference

	Left	Right
C1	_____	_____
C	_____	_____
B	_____	_____
A	_____	_____

**G Length**  
 Regular: 30–33.9cm  
 Tall: 34–42cm

G \_\_\_\_\_



### Compreflex Transition Calf Size Chart (A–C1)

	Small	Medium	Large	X Large	XX Large
C1	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
B	24–34cm	29–39cm	34–44cm	39–49cm	44–55cm
A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm

Black					
Regular	1401-UC-BKR	1402-UC-BKR	1403-UC-BKR	1404-UC-BKR	1405-UC-BKR
Tall	1401-UC-BKT	1402-UC-BKT	1403-UC-BKT	1404-UC-BKT	1405-UC-BKT
Beige					
Regular	1411-UC-BKR	1412-UC-BKR	1413-UC-BKR	1414-UC-BKR	1415-UC-BKR
Tall	1411-UC-BKT	1412-UC-BKT	1413-UC-BKT	1414-UC-BKT	1415-UC-BKT

### Compreflex Strap Extender Chart

	1 Strap	10 Straps
Black	1401-SE1	1401-SEB
Beige	1411-SE1	1411-SEB