




# sigvaris

## Comprevest

**Patient Last Name:** \_\_\_\_\_ **Patient First Name:** \_\_\_\_\_  
**Fitter Last Name:** \_\_\_\_\_ **Fitter First Name:** \_\_\_\_\_  
**Fitter Title:** \_\_\_\_\_ **(example: PT/OT/PTA)**  
**Date:** \_\_\_\_\_

### Step One


Measure the circumference in inches at the fullest part of bust.



Bust Size

### Step Two

Measure the circumference in inches directly under the bust.



Band Size

### Step Three

Determine the cup by subtracting the band from the bust.

Cup Size

2" or less = **A/B**  
 2" - 4" = **C/D**  
 4" - 6" = **DD/DDD**

#### Product Information

Product includes one Comprevest.

Comprevest      Cup Size: \_\_\_\_\_      Item #: \_\_\_\_\_  
 Color:  Black    White      Band Size: \_\_\_\_\_      Quantity: \_\_\_\_\_

#### Companion Products (sold separately)

Chip Pad Bra Half      Item #: \_\_\_\_\_  
 Color:  Black    White      Quantity: \_\_\_\_\_

Chip Pad Bra Full      Item #: \_\_\_\_\_  
 Color:  Black    White      Quantity: \_\_\_\_\_

Drain Tube Pocket      Item #: \_\_\_\_\_  
 Color:  Black      Quantity: \_\_\_\_\_

**Caution:** Elastic band contains natural rubber latex.

#### Comprevest Size Chart

	Small	Medium	Large	X Large	XX Large	
<b>Band Size</b>	28-30in	32-34in	36-38in	40-42in	44-46in	
<b>Cup Size</b>						
<b>Black</b>	A/B	2201-BRR-A	2202-BRR-A	2203-BRR-A	2204-BRR-A	2205-BRR-A
	C/D	2201-BRR-C	2202-BRR-C	2203-BRR-C	2204-BRR-C	2205-BRR-C
	DD/DDD	2201-BRR-D	2202-BRR-D	2203-BRR-D	2204-BRR-D	2205-BRR-D
<b>White</b>	A/B	2281-BRR-A	2282-BRR-A	2283-BRR-A	2284-BRR-A	2285-BRR-A
	C/D	2281-BRR-C	2282-BRR-C	2283-BRR-C	2284-BRR-C	2285-BRR-C
	DD/DDD	2281-BRR-D	2282-BRR-D	2283-BRR-D	2284-BRR-D	2285-BRR-D