



**sigvaris**  
Medaboot

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example: PT/OT/PTA)  
 Date: \_\_\_\_\_

**Circumference**

Left Right

I \_\_\_\_\_

J \_\_\_\_\_


L \_\_\_\_\_

**Length**

Left Right

K \_\_\_\_\_

M \_\_\_\_\_



**K:** Measure **medial** length from heel to 1st metatarsal head

**Product Information**

Product includes one Medaboot.

Right Foot      Color: Black      Quantity:

Left Foot      Color: Black      Quantity:

**Important**

Exact measurements are critical for this garment to ensure a proper fit. If you'd like to learn more about measuring and fitting, attend our MCE Basic Fitter Training. To sign up, visit [sigvariseducation.com](http://sigvariseducation.com) or contact your local territory manager for more information.

Alternatively, call Customer Care Solution Center at **800-322-7744**, or e-mail [us\\_orders@sigvaris.com](mailto:us_orders@sigvaris.com), to receive a remote consultation/training.

**Supplies Needed**

- Cell phone with camera. Photos of the foot with measurement markings must be emailed to: [us\\_orders@sigvaris.com](mailto:us_orders@sigvaris.com)
- Measuring instructions and forms.
- SIGVARIS GROUP Measuring tape and body pen (or eyeliner pencil).
- Signed Custom Order Terms & Conditions Form (include with order).

**MEASURING INSTRUCTIONS**

All measurements should be recorded in centimeters. Apply slight tension to hold the tape measure in place.

**MEASURE FOOT**

1. Measure foot medially from heel to 1st metatarsal head (or desired boot length) and record length on line K.
2. Measure top of foot from the 3rd metatarsal head to ankle bend and record length on line M.
3. Encircle the ankle bend and heel with a tape measure and record the circumference on line I.
4. Encircle the foot across the metatarsal heads and record the circumference on line J.
5. Measure the ankle circumference and record on line L.