




SIGVARIS

NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASURING FORM & CHART

Patients Name:	Date:
Fitters Name:	Fitters Title: (example PT/OT/PTA)

SEND TO:
Luna Medical, Inc.
 1057 W Grand Ave, Ste. 1
 Chicago, IL 60642
lunamedical.com
 P: 800.380.4339
 F: 888.696.0299

LOWER EXTREMITY



CALF STYLE

- 1 Ankle Circumference at narrowest part of ankle, above ankle bone **(b)** _____
- 2 Calf Circumference at fullest part of the calf **(c)** _____
- 3 Calf Length from floor up to the bend of the knee (outside of calf) **(A-D)** _____

Series	Features	Best Suited For
MIDTOWN MICROFIBER & MICROFIBER SHADES	Advanced microfibre technology for ultimate comfort	Patients who prefer fashionable options
87% Nylon 13% Spandex	Glides on and off with ease Excellent breathability and moisture management	Management of venous insufficiency symptoms
(contents may vary by style)	Soft and durable sock for everyday wear	Long-distance travel
Latex-Free		

Donning Accessories	
<input type="checkbox"/> SIGVARIS ARION Sim-Slide (Open Toe Donning Device)	<input type="checkbox"/> SIGVARIS ARION Magnide On/Off (Closed Toe Donning Device)
<input type="checkbox"/> Latex-Free Donning Gloves	<input type="checkbox"/> SIMON (Donning Device)
<input type="checkbox"/> Rubber Donning Gloves	<input type="checkbox"/> SIGVARIS Doff N' Donner (Donning Device)
<input type="checkbox"/> Cotton Underliner	<input type="checkbox"/> S.O.S. (Slip On SIGVARIS Donning Device)

Products	Compression	Style		Colors
820 MIDTOWN MICROFIBER* (for men only)	<input type="checkbox"/> 15-20 mmHg <input type="checkbox"/> 20-30 mmHg <input type="checkbox"/> 30-40 mmHg	Closed Toe: <input type="checkbox"/> Calf <input type="checkbox"/> Calf w/grip-top	<input type="checkbox"/> Thigh w/grip-top <input type="checkbox"/> Extended foot sizes	<input type="checkbox"/> Navy (10) <input type="checkbox"/> Tan Khaki (32) <input type="checkbox"/> Steel Grey (93) <input type="checkbox"/> Black (99)
830 MICROFIBER SHADES	<input type="checkbox"/> 20-30 mmHg	Closed Toe: <input type="checkbox"/> Calf	<input type="checkbox"/> Dark Navy Mini-Stripe (03) <input type="checkbox"/> Graphite Mini-Stripe (03) <input type="checkbox"/> Dark Navy Stripe (13) <input type="checkbox"/> Pink Stripe (15) <input type="checkbox"/> Graphite Stripe (17)	<input type="checkbox"/> Onyx Stripe (19) <input type="checkbox"/> Graphite Heather (27) <input type="checkbox"/> Pink Argyle (45) <input type="checkbox"/> Purple Argyle (46) <input type="checkbox"/> Graphite Argyle (45) <input type="checkbox"/> Royal Blue Argyle (45)

*For all available styles, colors, and compression, please see the SIGVARIS product catalog.

For indications and contraindications, please see prescribing information at sigvarisusa.com

SIGVARIS is a registered trademark of SIGVARIS, CH-9014 Svt.Gallen/Switzerland, in many countries worldwide. © Copyright 2018 SIGVARIS, Inc. Rev.12.18