




# SIGVARIS

## NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASURING FORM & CHART

Patients Name:	Date:
Fitters Name:	Fitters Title: (example PT/OT/PTA)

**SEND TO:**  
**Luna Medical, Inc.**  
 1057 W Grand Ave, Ste. 1  
 Chicago, IL 60642  
**lunamedical.com**  
 P: 800.380.4339  
 F: 888.696.0299

### LOWER EXTREMITY



#### CALF STYLE

- 1 Ankle Circumference at narrowest part of ankle, above ankle bone **(b)** \_\_\_\_\_
- 2 Calf Circumference at fullest part of the calf **(c)** \_\_\_\_\_
- 3 Calf Length from floor up to the bend of the knee (outside of calf) **(A-D)** \_\_\_\_\_

Donning Accessories	
<input type="checkbox"/> SIGVARIS ARION Sim-Slide (Open Toe Donning Device)	<input type="checkbox"/> SIGVARIS ARION Magnide On/Off (Closed Toe Donning Device)
<input type="checkbox"/> Latex-Free Donning Gloves	<input type="checkbox"/> SIMON (Donning Device)
<input type="checkbox"/> Rubber Donning Gloves	<input type="checkbox"/> SIGVARIS Doff N' Donner (Donning Device)
<input type="checkbox"/> Cotton Underliner	<input type="checkbox"/> S.O.S. (Slip On SIGVARIS Donning Device)

Series	Features	Best Suited For
SEA ISLAND COTTON 43% Sea Island Cotton 39% Nylon 11% Spandex 7% Supima Cotton Latex-Free	Knitted from Sea Island Cotton, the world's rarest cotton, for a fashionable look and a soft, luxurious feel  Only cotton touches the skin to ensure comfort  Breathable and comfortable for sensitive skin	Patients who prefer natural fibres for everyday wear  Management of venous insufficiency symptoms  Long-distance travel
ALL SEASON MERINO WOOL 53% Merino Wool 39% Nylon 9% Spandex Latex-Free	Soft, comfortable, and itch-free Made from the highest quality extra-fine Merino Wool  Natural properties of wool make these socks odour-resistant  Thermoregulating for year-round wear and keeps feet dry	Patients who prefer fashionable and premium fabric options  Management of venous insufficiency symptoms  Long-distance travel
CUSHIONED COTTON 64% Supima Cotton 28% Nylon 8% Spandex Latex-Free	Thick, padded sole, heel and toe to protect feet  Flat, low-friction toe seam Soft, stretchable knit-in band prevents pinching  Ideal for everyday wear and for sensitive skin  Recommended for patients with diabetes	People who stand for prolonged periods of time  Work shoes and boots  Casual shoes and sports shoes Diabetic friendly

Products	Compression	Style	Colors
220 SEA ISLAND COTTON	<input type="checkbox"/> 20-30 mmHg	<b>Closed Toe:</b> <input type="checkbox"/> Calf	<input type="checkbox"/> Navy (10) <input type="checkbox"/> Black (99) <input type="checkbox"/> Brown (11)
240 ALL SEASON MERINO WOOL	<input type="checkbox"/> 20-30 mmHg	<b>Closed Toe:</b> <input type="checkbox"/> Calf	<input type="checkbox"/> Navy (10) <input type="checkbox"/> Black (99) <input type="checkbox"/> Brown (11)
360 CUSHIONED COTTON	<input type="checkbox"/> 20-30 mmHg	<b>Closed Toe:</b> <input type="checkbox"/> Calf	<input type="checkbox"/> White (00) <input type="checkbox"/> Black (99)

For indications and contraindications, please see prescribing information at [sigvarisusa.com](http://sigvarisusa.com)