




SIGVARIS

NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASURING FORM & CHART

SEND TO:
Luna Medical, Inc.
 1057 W Grand Ave, Ste. 1
 Chicago, IL 60642
lunamedical.com
 P: 800.380.4339
 F: 888.696.0299

Patients Name:	Date:
Fitters Name:	Fitters Title: (example PT/OT/PTA)

UPPER EXTREMITY



ARMSLEEVE WITHOUT GAUNTLET

- Wrist Circumference at narrowest part of wrist (above wrist bone)(**c**) _____
- Axilla (underarm) Circumference just underneath the armpit. (**f**) _____
- Arm length from wrist to axilla (measure outside of arm only)(**C-G**) _____

ARMSLEEVE WITH GAUNTLET

- (continue from steps 1-3)
- Palm Circumference at the widest part of hand (thumb excluded)(**a**) _____

Series	Features	Best Suited For
SECURE ARM 78% Nylon 22% Spandex (contents may vary by style) Not made with natural Latex rubber	High modulus in-lay yarns provide containment and work to maintain limb profile Secure dotted grip-top band is breathable and appropriate for sensitive skin Double-covered yarns ensure no bare spandex next to the skin Unique sizing system fits a wide range of regular and full sizes—up to 20 inch arms UPF 50+ sun protection	Ideal for patients with chronic edema Lymphedema Post-operative edema Other edema
GAUNTLET & GLOVE 60% Nylon 40% Spandex (contents may vary by style) Not made with natural Latex rubber	Seamless garment for optimal comfort Circular knit softness and increased range of motion Anatomically shaped with relief zones for thumb comfort Breathable fabric is color coordinated to matching armsleeves Wrist area has less compression to overlap with armsleeves Fingers can be adjusted to length and resist rolling down	Ideal for patients with chronic edema Mild to moderate lymphedema Post-operative edema Other edema

Donning Accessories	
<input type="checkbox"/> SIGVARIS ARION Sim-Slide (Open Toe Donning Device)	<input type="checkbox"/> SIGVARIS ARION Magnide On/Off (Closed Toe Donning Device)
<input type="checkbox"/> Latex-Free Donning Gloves	<input type="checkbox"/> SIMON (Donning Device)
<input type="checkbox"/> Rubber Donning Gloves	<input type="checkbox"/> SIGVARIS Doff N' Donner (Donning Device)
<input type="checkbox"/> Cotton Underliner	<input type="checkbox"/> S.O.S. (Slip On SIGVARIS Donning Device)

Products	Compression	Style	Colors
560 SECURE ARM*	<input type="checkbox"/> 15-20 mmHg <input type="checkbox"/> 20-30 mmHg <input type="checkbox"/> 30-40 mmHg	Armsleeve: <input type="checkbox"/> Armsleeve w/ grip top Sold Separately: <input type="checkbox"/> Glove* <input type="checkbox"/> Gauntlet*	<input type="checkbox"/> Black (99) <input type="checkbox"/> Beige (77) <input type="checkbox"/> Coco (84) <input type="checkbox"/> Dusty Rose (65)

*For all available styles, colors, and compression, please see the SIGVARIS product catalog.

For indications and contraindications, please see prescribing information at sigvarisusa.com