




SIGVARIS

NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASURING FORM & CHART

Patients Name:	Date:
Fitters Name:	Fitters Title: (example PT/OT/PTA)

SEND TO:
Luna Medical, Inc.
 1057 W Grand Ave, Ste. 1
 Chicago, IL 60642
lunamedical.com
 P: 800.380.4339
 F: 888.696.0299

LOWER EXTREMITY



CALF STYLE

- 1 Ankle Circumference at narrowest part of ankle, above ankle bone **(b)** _____
- 2 Calf Circumference at fullest part of the calf **(c)** _____
- 3 Calf Length from floor up to the bend of the knee (outside of calf) **(A-D)** _____

PANTYHOSE AND THIGH-HIGH STYLES (continue from steps 1-3)

- 4 Thigh Circumference at its widest circumference **(g)** _____
- 5 Thigh Length from floor up to the gluteal fold **(A-G)** _____
- 6 Hip Circumference at widest part of the hip **(f)** _____

Donning Accessories	
<input type="checkbox"/> SIGVARIS ARION Sim-Slide (Open Toe Donning Device)	<input type="checkbox"/> SIGVARIS ARION Magnide On/Off (Closed Toe Donning Device)
<input type="checkbox"/> Latex-Free Donning Gloves	<input type="checkbox"/> SIMON (Donning Device)
<input type="checkbox"/> Rubber Donning Gloves	<input type="checkbox"/> SIGVARIS Doff N' Donner (Donning Device)
<input type="checkbox"/> Cotton Underliner	<input type="checkbox"/> S.O.S. (Slip On SIGVARIS Donning Device)

Series	Features	Best Suited For
SELECT COMFORT 64% Nylon 36% Spandex (contents may vary by style) Latex-Free	From petite to plus, SELECT COMFORT offers a complete line of sizes designed to fit more than 99% of the population Durable construction combined with a smooth, opaque look Comfortable, stretchable material constructed from double covered yarn easily glides on and off Plus-size pantyhose is a great alternative for maternity Cotton-lined gusset on pantyhose style	Post-vein treatment (stripping, ligation, phlebectomy, sclero, EVLT, RF ablation)
COTTON 50% Nylon 25% Spandex 25% Supima Cotton (contents may vary by style) Latex-Free	Only Supima Cotton touches the skin to ensure comfort Breathable and durable Classic ribbed look Ideal for everyday wear and for sensitive skin	Containment of venous edema and early stages of lymphedema All-day cotton comfort in all types of climates Ideal for dress or casual shoes Ideal for long-distance travel

Products	Compression	Style	Colors
860 SELECT COMFORT*	<input type="checkbox"/> 20-30 mmHg <input type="checkbox"/> 30-40 mmHg	Closed Toe: <input type="checkbox"/> Calf <input type="checkbox"/> Petite calf <input type="checkbox"/> Calf plus <input type="checkbox"/> Calf w/grip-top <input type="checkbox"/> Thigh w/grip-top <input type="checkbox"/> Petite thigh-high w/grip-top <input type="checkbox"/> Extra-large thigh w/grip-top <input type="checkbox"/> Pantyhose <input type="checkbox"/> Petite pantyhose <input type="checkbox"/> Plus-size Pantyhose Open Toe: <input type="checkbox"/> Calf <input type="checkbox"/> Calf plus <input type="checkbox"/> Calf w/grip-top <input type="checkbox"/> Thigh w/grip-top <input type="checkbox"/> Extra-large thigh w/grip-top <input type="checkbox"/> Thigh-high w/ waist attachment	<input type="checkbox"/> Natural (33) <input type="checkbox"/> Suntan (36) <input type="checkbox"/> Crispa (66) <input type="checkbox"/> Black (99)
230 COTTON*	<input type="checkbox"/> 20-30 mmHg <input type="checkbox"/> 30-40 mmHg	Closed Toe: <input type="checkbox"/> Calf	<input type="checkbox"/> White (00) <input type="checkbox"/> Navy (10) <input type="checkbox"/> Black Mist (14) <input type="checkbox"/> Crispa (66) <input type="checkbox"/> Chocolate (88) <input type="checkbox"/> Black (99)

*For all available styles, colors, and compression, please see the SIGVARIS product catalog.