




# SIGVARIS

## NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASURING FORM & CHART

Patients Name:	Date:
Fitters Name:	Fitters Title: (example PT/OT/PTA)

**SEND TO:**  
**Luna Medical, Inc.**  
 1057 W Grand Ave, Ste. 1  
 Chicago, IL 60642  
**lunamedical.com**  
 P: 800.380.4339  
 F: 888.696.0299

### LOWER EXTREMITY



#### CALF STYLE

- 1 Ankle Circumference at narrowest part of ankle, above ankle bone (**b**) \_\_\_\_\_
- 2 Calf Circumference at fullest part of the calf (**c**) \_\_\_\_\_
- 3 Calf Length from floor up to the bend of the knee (outside of calf) (**A-D**) \_\_\_\_\_

#### PANTYHOSE AND THIGH-HIGH STYLES (continue from steps 1-3)

- 4 Thigh Circumference at its widest circumference (**g**) \_\_\_\_\_
- 5 Thigh Length from floor up to the gluteal fold (**A-G**) \_\_\_\_\_
- 6 Hip Circumference at widest part of the hip (**f**) \_\_\_\_\_

#### Donning Accessories

- |   |  |
|---|--|
| <input type="checkbox"/> SIGVARIS ARION Sim-Slide (Open Toe Donning Device) | <input type="checkbox"/> SIGVARIS ARION Magnide On/Off (Closed Toe Donning Device) |
| <input type="checkbox"/> Latex-Free Donning Gloves                          | <input type="checkbox"/> SIMON (Donning Device)                                    |
| <input type="checkbox"/> Rubber Donning Gloves                              | <input type="checkbox"/> SIGVARIS Doff N' Donner (Donning Device)                  |
| <input type="checkbox"/> Cotton Underliner                                  | <input type="checkbox"/> S.O.S. (Slip On SIGVARIS Donning Device)                  |

Series	Features	Best Suited For
SOFT OPAQUE	Comfortable tights with lasting microfibre softness	Patients who prefer fashionable options
79% Nylon 21% Spandex	Advanced moisture management that keeps legs cool and dry	Pre- and post-vein treatment with minimal edema
(content can vary slightly by style)	Comfortable wide top band prevents knee-highs from rolling down	Pregnancy
Latex-Free	Cotton-lined gusset on pantyhose style	Long-distance travel
EVERSHEER	Fashionable, sheer medical graduated compression hosiery	Patients who prefer fashionable options
64% Nylon 36% Spandex	Constructed from double-covered inlay yarns for softness against the skin	Pre- and post-vein treatment with minimal edema
(content may vary by style)	Excellent breathability and moisture management for added comfort	Ideal for business or dress
Latex-Free	Cotton-lined gusset on pantyhose style	

Products	Compression	Style	Colors
780 EVERSHEER* (for women only)	<input type="checkbox"/> 15-20 mmHg <input type="checkbox"/> 20-30 mmHg <input type="checkbox"/> 30-40 mmHg	<b>Closed Toe:</b> <input type="checkbox"/> Calf <input type="checkbox"/> Thigh-high w/grip-top <input type="checkbox"/> Pantyhose <input type="checkbox"/> Maternity Pantyhose <b>Open Toe:</b> <input type="checkbox"/> Calf <input type="checkbox"/> Thigh-high w/grip-top <input type="checkbox"/> Pantyhose <input type="checkbox"/> Maternity Pantyhose	<input type="checkbox"/> Dark Navy (08) <input type="checkbox"/> Natural (33) <input type="checkbox"/> Warm Sand (34) <input type="checkbox"/> Toasted Almond (37) <input type="checkbox"/> Suntan (36) <input type="checkbox"/> Café (73) <input type="checkbox"/> Mocha (85) <input type="checkbox"/> Nightshade (34) <input type="checkbox"/> Black (99)
840 SOFT OPAQUE* (for women only)	<input type="checkbox"/> 15-20 mmHg <input type="checkbox"/> 20-30 mmHg <input type="checkbox"/> 30-40 mmHg	<b>Closed Toe:</b> <input type="checkbox"/> Calf <input type="checkbox"/> Calf w/grip-top <input type="checkbox"/> Thigh-high w/grip-top <input type="checkbox"/> Pantyhose <input type="checkbox"/> Maternity Pantyhose <b>Open Toe:</b> <input type="checkbox"/> Calf <input type="checkbox"/> Thigh-high w/grip-top <input type="checkbox"/> Pantyhose	<input type="checkbox"/> Graphite (91) <input type="checkbox"/> Black (99)

\*For all available styles, colors, and compression, please see the SIGVARIS product catalog.

For indications and contraindications, please see prescribing information at [sigvarisusa.com](http://sigvarisusa.com)

SIGVARIS is a registered trademark of SIGVARIS, CH-9014 Svt.Gallen/Switzerland, in many countries worldwide. © Copyright 2018 SIGVARIS, Inc. Rev.12.18