




# SIGVARIS

## NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASURING FORM & CHART

|                |                                       |
|----------------|---------------------------------------|
| Patients Name: | Date:                                 |
| Fitters Name:  | Fitters Title:<br>(example PT/OT/PTA) |

**SEND TO:**  
**Luna Medical, Inc.**  
 1057 W Grand Ave, Ste. 1  
 Chicago, IL 60642  
**lunamedical.com**  
 P: 800.380.4339  
 F: 888.696.0299

### LOWER EXTREMITY



#### CALF STYLE

- 1 Ankle Circumference at narrowest part of ankle, above ankle bone **(b)** \_\_\_\_\_
- 2 Calf Circumference at fullest part of the calf **(c)** \_\_\_\_\_
- 3 Calf Length from floor up to the bend of the knee (outside of calf) **(A-D)** \_\_\_\_\_

#### THIGH-HIGH STYLE

(continue from steps 1-3)

- 4 Thigh Circumference at its widest circumference **(g)** \_\_\_\_\_
- 5 Thigh Length from floor up to the gluteal fold **(A-G)** \_\_\_\_\_
- 6 Hip Circumference at widest part of the hip **(f)** \_\_\_\_\_

| Series                             | Features  | Best Suited For   |
|------------------------------------|---|---|
| ALLURE                             | Stylish diamond-dot pattern for fashionable daily wear            | Management of early stages of chronic venous disease                |
| 64% Nylon<br>36% Spandex           | Durable construction with reinforced heel                         | Pre- and post-vein treatment with minimal edema                     |
| (contents may vary by style)       | Excellent breathability and moisture management for added comfort | Patients who prefer fashionable options Ideal for business or dress |
| Not made with natural Latex rubber | Cotton-lined gusset on pantyhose style                            |   |
|                                    | Sensinnov® band grip-top on thigh style                           |   |

#### Donning Accessories

|   |  |
|---|--|
| <input type="checkbox"/> SIGVARIS ARION Sim-Slide (Open Toe Donning Device) | <input type="checkbox"/> SIGVARIS ARION Magnide On/Off (Closed Toe Donning Device) |
| <input type="checkbox"/> Latex-Free Donning Gloves                          | <input type="checkbox"/> SIMON (Donning Device)                                    |
| <input type="checkbox"/> Rubber Donning Gloves                              | <input type="checkbox"/> SIGVARIS Doff N' Donner (Donning Device)                  |
| <input type="checkbox"/> Cotton Underliner                                  | <input type="checkbox"/> S.O.S. (Slip On SIGVARIS Donning Device)                  |

| Products    | Compression  | Style   | Colors   |
|-------------|--|---|--|
| 710 ALLURE* | <input type="checkbox"/> 15-20 mmHg<br><input type="checkbox"/> 20-30 mmHg | <b>Closed Toe:</b><br><input type="checkbox"/> Pantyhose<br><input type="checkbox"/> Thigh w/grip-top | <input type="checkbox"/> Graphite (91) <input type="checkbox"/> Black (99) |

\*For all available styles, colors, and compression, please see the SIGVARIS product catalog.

For indications and contraindications, please see prescribing information at [sigvarisusa.com](http://sigvarisusa.com)

Sigvaris and SIGVARIS GROUP are registered trademarks of SIGVARIS, CH-9014 Svt.Gallen/Switzerland, in many countries worldwide. © Copyright 2019 SIGVARIS, Inc.