



Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____

sigvaris


Chipsleeve Standard Arm

Product Information			
Product includes one Chipsleeve Standard Arm, one pair of Cotton Liners, and one black Oversleeve.			
<input type="checkbox"/> Right Arm	Size:	Item #:	
Color: Black	Length:	Quantity:	
<input type="checkbox"/> Left Arm	Size:	Item #:	
Color: Black	Length:	Quantity:	

Circumference

	Left	Right
C	_____	_____
BI	_____	_____
B	_____	_____
A	_____	_____

D Length	
Short:	38–42.9cm
Regular:	43–47.9cm
Long:	48–53cm
<i>Measure length along dorsal aspect</i>	
D	_____



Oversleeve Arm Size Chart					
Color	Small	Medium	Large	X Large	XX Large
Black	1801-ARR	1802-ARR	1803-ARR	1804-ARR	1805-ARR

Chipsleeve Standard Arm Size Chart (A–C) (Available in Black Only)					
	Small	Medium	Large	X Large	XX Large
C	23–32cm	28–37cm	33–43cm	39–49cm	44–54cm
BI	22–30cm	26–34cm	30–39cm	35–44cm	39–49cm
B	20–27cm	24–31cm	28–35cm	32–39cm	36–43cm
A	14–17cm	15–18cm	16–19cm	18–20cm	19–22cm
Left					
Short	2631-ARS-L	2632-ARS-L	2633-ARS-L	2634-ARS-L	2635-ARS-L
Regular	2631-ARR-L	2632-ARR-L	2633-ARR-L	2634-ARR-L	2635-ARR-L
Long	2631-ARL-L	2632-ARL-L	2633-ARL-L	2634-ARL-L	2635-ARL-L
Right					
Short	2631-ARS-R	2632-ARS-R	2633-ARS-R	2634-ARS-R	2635-ARS-R
Regular	2631-ARR-R	2632-ARR-R	2633-ARR-R	2634-ARR-R	2635-ARR-R
Long	2631-ARL-R	2632-ARL-R	2633-ARL-R	2634-ARL-R	2635-ARL-R