



**Patient Last Name:** \_\_\_\_\_ **Patient First Name:** \_\_\_\_\_  
**Fitter Last Name:** \_\_\_\_\_ **Fitter First Name:** \_\_\_\_\_  
**Fitter Title:** \_\_\_\_\_ **(example: PT/OT/PTA)**  
**Date:** \_\_\_\_\_

# sigvaris

## Compreflex Reduce Foot

### Circumference


Left    Right

I    \_\_\_\_\_    \_\_\_\_\_

J    \_\_\_\_\_    \_\_\_\_\_

### Length

K    \_\_\_\_\_    \_\_\_\_\_



**K:** Measure **medial** length from heel to 1st metatarsal head

Product Information			
Product includes one Compreflex Reduce Foot.			
<input type="checkbox"/> Right Foot	Size:	Item #:	Quantity:
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige			
<input type="checkbox"/> Left Foot	Size:	Item #:	Quantity:
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige			

Compreflex Reduce Foot Size Chart (I-K)		
	Range 1	Range 3
I	30-44cm	39-54cm
J	24-33cm	30-43cm
K	16-25cm	20-27cm
<b>Black</b>	2801-BTL	2803-BTL
<b>Beige</b>	2811-BTL	2813-BTL

**Note:** Compreflex Reduce Calf options available, each sold separately.

## Measuring Instructions

### Compreflex Reduce Foot:

Measure circumference of ankle bend and heel and record on line I. Measure circumference across metatarsal heads and record on line J. Measure length on medial side of foot from heel to 1st metatarsal head and record on line K.