



Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example: PT/OT/PTA)  
 Date: \_\_\_\_\_

# sigvaris

## Compreflex Reduce Knee

### Product Information

Product includes one Compreflex Reduce Knee wrap and one Butterfly Strap

Right Knee      Size: \_\_\_\_\_      Item #: \_\_\_\_\_

Color:  Black    Beige      Length: \_\_\_\_\_      Quantity: \_\_\_\_\_

Left Knee      Size: \_\_\_\_\_      Item #: \_\_\_\_\_

Color:  Black    Beige      Length: \_\_\_\_\_      Quantity: \_\_\_\_\_

### Companion Products (sold separately)

Compreflex Reduce Thigh      Color:  Black    Beige

Size: \_\_\_\_\_      Item #: \_\_\_\_\_      Quantity: \_\_\_\_\_

Size: \_\_\_\_\_      Item #: \_\_\_\_\_      Quantity: \_\_\_\_\_

### Measuring Instructions

#### Compreflex Reduce Knee:

Measure circumferences where indicated, beginning 5cm below the lateral aspect fibular head. Continue 5cm above the lateral aspect top of patella. Record in the corresponding lines C1 and D.

#### Compreflex Reduce Thigh:

Measure circumferences where indicated on the measurement chart, beginning 5cm above the top of patella and continue 23cm up the thigh, just below the gluteal fold. Record in the corresponding lines D and F. Measure length from the lateral aspect top of patella to the gluteal fold, and record in box G1.

### Circumference

Left      Right

F \_\_\_\_\_

E \_\_\_\_\_

D \_\_\_\_\_

C1 \_\_\_\_\_


### G1 Length

Regular: 20–24.9cm

Tall: 25–30cm

*Top of patella to  
gluteal fold*

G1 \_\_\_\_\_



### Compreflex Reduce Knee Size Chart

	Range 1	Range 3
D	<75cm	<115cm
C1	<65cm	<90cm
<b>Black</b>	2401-KP	2403-KP
<b>Beige</b>	2411-KP	2413-KP

### Compreflex Reduce Thigh Size Chart

	Range 1	Range 3
F	<80cm	<125cm
D	<65cm	<110cm

		<b>Black</b>	
Left	Regular	2401-TCR-L	2403-TCR-L
	Tall	2401-TCT-L	2403-TCT-L
Right	Regular	2401-TCR-R	2403-TCR-R
	Tall	2401-TCT-R	2403-TCT-R
		<b>Beige</b>	
Left	Regular	2411-TCR-L	2413-TCR-L
	Tall	2411-TCT-L	2413-TCT-L
Right	Regular	2411-TCR-R	2413-TCR-R
	Tall	2411-TCT-R	2413-TCT-R