



Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____

sigvaris

Compreflex Reduce Thigh

Product Information

Product includes one Compreflex Reduce Thigh wrap, one pair of Cotton Liners, and one Belt and Hip attachment.

<input type="checkbox"/> Right Thigh	Size: _____	Item #: _____
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length: _____	Quantity: _____
<input type="checkbox"/> Left Thigh	Size: _____	Item #: _____
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length: _____	Quantity: _____

Companion Products (sold separately)

<input type="checkbox"/> Compreflex Reduce Knee	Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige
Size: _____	Item #: _____
Quantity: _____	
Size: _____	Item #: _____
Quantity: _____	

Measuring Instructions

Compreflex Reduce Thigh:

Measure circumferences where indicated on the measurement chart, beginning 5cm above the top of patella and continue 23cm up the thigh, just below the gluteal fold. Record in the corresponding lines D and F. Measure length from the lateral aspect top of patella to the gluteal fold, and record in box G1.

Compreflex Reduce Knee:

Measure circumferences where indicated, beginning 5cm below the lateral aspect fibular head. Continue 5cm above the lateral aspect top of patella. Record in the corresponding lines C1 and D.


Circumference

	Left	Right
F	_____	_____
E	_____	_____
D	_____	_____
C1	_____	_____

G1 Length
 Regular: 20–24.9cm
 Tall: 25–30cm

Top of patella to gluteal fold

G1 _____



Compreflex Reduce Thigh Size Chart

	Range 1	Range 3
F	<80cm	<125cm
D	<65cm	<110cm

		Black	
Left	Regular	2401-TCR-L	2403-TCR-L
	Tall	2401-TCT-L	2403-TCT-L
Right	Regular	2401-TCR-R	2403-TCR-R
	Tall	2401-TCT-R	2403-TCT-R

		Beige	
Left	Regular	2411-TCR-L	2413-TCR-L
	Tall	2411-TCT-L	2413-TCT-L
Right	Regular	2411-TCR-R	2413-TCR-R
	Tall	2411-TCT-R	2413-TCT-R

Compreflex Reduce Knee Size Chart

	Range 1	Range 3
D	<75cm	<115cm
C1	<65cm	<90cm

		Black	
		2401-KP	2403-KP
		Beige	
		2411-KP	2413-KP