



Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____

sigvaris

Compreflex Standard Arm

Product Information

Product includes one Compreflex Standard Arm and one pair of Cotton Liners.

<input type="checkbox"/> Right Arm	Size:	Item #:
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:
<input type="checkbox"/> Left Arm	Size:	Item #:
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:

Companion Products (sold separately)

<input type="checkbox"/> Dorsal Pocket Glove	Size:	Item #:
<input type="checkbox"/> Right Hand	Quantity:	
<input type="checkbox"/> Dorsal Pocket Glove	Size:	Item #:
<input type="checkbox"/> Left Hand	Quantity:	

Dorsal Pocket Glove Size Chart (A- E) (Available in Black Only)

	Small	Medium	Large	X Large
E	15-19cm	19-22cm	22-26cm	26-29cm
A	14-20cm	16-22cm	18-24cm	20-26cm
Black				
Right	1201-DGR-R	1202-DGR-R	1203-DGR-R	1204-DGR-R
Left	1201-DGR-L	1202-DGR-L	1203-DGR-L	1204-DGR-L

Circumference

Left Right


C _____
 BI _____
 B _____
 A _____
 E _____

D Length

Short: 38-42.9cm
 Regular: 43-47.9cm
 Long: 48-53cm

Measure length along dorsal aspect

D _____



Compreflex Standard Arm Size Chart (A-C)

	Small	Medium	Large	X Large
C	23-32cm	28-37cm	33-43cm	39-49cm
BI	22-30cm	26-34cm	30-39cm	35-44cm
B	20-27cm	24-31cm	28-35cm	32-39cm
A	14-19cm	15-20cm	16-21cm	18-23cm
Black				
Short	1401-ARS	1402-ARS	1403-ARS	1404-ARS
Regular	1401-ARR	1402-ARR	1403-ARR	1404-ARR
Long	1401-ARL	1402-ARL	1403-ARL	1404-ARL
Beige				
Short	1411-ARS	1412-ARS	1413-ARS	1414-ARS
Regular	1411-ARR	1412-ARR	1413-ARR	1414-ARR
Long	1411-ARL	1412-ARL	1413-ARL	1414-ARL