



**Patient Last Name:** \_\_\_\_\_ **Patient First Name:** \_\_\_\_\_  
**Fitter Last Name:** \_\_\_\_\_ **Fitter First Name:** \_\_\_\_\_  
**Fitter Title:** \_\_\_\_\_ **(example: PT/OT/PTA)**  
**Date:** \_\_\_\_\_

## FEATURES

Highly Breathable Spacer Fabric TR for cool-wearing comfort

EZ-Close Pressure System of eight tabs makes one-handed, lateral donning simple, even for patients with dexterity issues

"Max" sizing for patients with small ankles and large calves

Durable metal clasps, not hook-and-loop closure

Each COOLFLEX NF includes one pair of moisture wicking basic liners to keep legs cool and dry

## MATERIAL CONTENTS

BASIC LINER: 73% Polyester 13% Nylon  
13% Cotton 1% Spandex

COOLFLEX does not contain natural rubber latex

## INDICATIONS

- Chronic Venous Disease
- Venous Leg Ulcers
- Post-thrombotic Syndrome
- Phlebolympheidema
- Lymphedema
- Dependent Edema

## PRECAUTIONS & CONTRAINDICATIONS

Refer to SIGVARIS 2018 catalog or package insert for full product information

## CARE INSTRUCTIONS

- Wash on a gentle cycle using a mild detergent
- DO NOT use bleach or fabric softeners
- Dry on low/no heat or air dry (air drying is recommended to extend garment life)



## SIZING & MEASURING

The COOLFLEX sizing system includes "MAX" sizing (highlighted in green below) to accommodate patients with small ankles and large calves

	SMALL	SMALL MAX	MEDIUM / LARGE	MED / LRG MAX	XLARGE / XXLARGE
C1	28-39cm	35-50cm	34-50cm	45-64cm	45-64cm
C	28-39cm	35-50cm	34-50cm	45-64cm	45-64cm
B	23-33cm	27-38cm	29-43cm	35-49cm	37-54cm
A	16-26cm	16-26cm	23-33cm	23-33cm	30-42cm

### LEFT

REGULAR	2101-NF-BKR-L	2101-NF-BMR-L	2103-NF-BKR-L	2103-NF-BMR-L	2105-NF-BKR-L
TALL	2101-NF-BKT-L	2101-NF-BMT-L	2103-NF-BKT-L	2103-NF-BMT-L	2105-NF-BKT-L

### RIGHT

REGULAR	2101-NF-BKR-R	2101-NF-BMR-R	2103-NF-BKR-R	2103-NF-BMR-R	2105-NF-BKR-R
TALL	2101-NF-BKT-R	2101-NF-BMT-R	2103-NF-BKT-R	2103-NF-BMT-R	2105-NF-BKT-R

