



Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example: PT/OT/PTA)  
 Date: \_\_\_\_\_

# sigvaris

## Coolflex Standard Calf

### Product Information

Product includes one Coolflex Standard Calf and one pair of non-compressive Basic Liners.

<input type="checkbox"/> Right Calf	Size:	Item #:
Color: Black	Length:	Quantity:
<input type="checkbox"/> Left Calf	Size:	Item #:
Color: Black	Length:	Quantity:

**Note:** Liner and foot compression options available, each sold separately.


### Circumference

Left Right

C1 \_\_\_\_\_  
 C \_\_\_\_\_  
 B \_\_\_\_\_  
 A \_\_\_\_\_

**G Length**  
 Regular: 30–35.9cm  
 Tall: 36–42cm

G \_\_\_\_\_



### Coolflex Standard Calf Size Chart (A-C1) (Available in Black Only)

	Small	Small Max	Medium / Large	Medium / Large Max	X Large / XX Large
C1	28–39cm	35–50cm	34–50cm	45–64cm	45–64cm
C	28–39cm	35–50cm	34–50cm	45–64cm	45–64cm
B	23–33cm	27–38cm	29–43cm	35–49cm	37–54cm
A	16–26cm	16–26cm	23–33cm	23–33cm	30–42cm

		Black				
Left	Regular	2101-NF-BKR-L	2101-NF-BMR-L	2103-NF-BKR-L	2103-NF-BMR-L	2105-NF-BKR-L
	Tall	2101-NF-BKT-L	2101-NF-BMT-L	2103-NF-BKT-L	2103-NF-BMT-L	2105-NF-BKT-L
Right	Regular	2101-NF-BKR-R	2101-NF-BMR-R	2103-NF-BKR-R	2103-NF-BMR-R	2105-NF-BKR-R
	Tall	2101-NF-BKT-R	2101-NF-BMT-R	2103-NF-BKT-R	2103-NF-BMT-R	2105-NF-BKT-R