




SIGVARIS

NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASURING FORM & CHART

Patients Name:	Date:
Fitters Name:	Fitters Title: (example PT/OT/PTA)

SEND TO:
Luna Medical, Inc.
 1057 W Grand Ave, Ste. 1
 Chicago, IL 60642
lunamedical.com
 P: 800.380.4339
 F: 888.696.0299

MEDICALLY COMPLEX EDEMA



Series & Features	Indications
CHIPSLEEVE <ul style="list-style-type: none"> • 15–25mmHg of mild, gradient compression from hand to axilla when used with the OverSleeve (included) • Open cell foam chips massage fibrotic tissue, stimulating fluid movement and helping to soften the superficial tissue • Use as a bandage liner or independent compression garment • Four-way stretch fabric contains moisture-wicking, odor control properties • Safe for day or nighttime use • Custom size available (see product catalog) • Easy care - machine wash and dry • Latex-free 	<ul style="list-style-type: none"> • Stage I or II Lymphedema • Mild fibrosis • Post-surgical edema
COMPREFLEX ARM <ul style="list-style-type: none"> • Provides mild/moderate compression • Unique elbow design for 90 degree range of motion • Soft, conforming Breathe-O-Prene® provides a comfortable fit and wicks away moisture • Easy to adjust straps • Available in black or beige • Easy care - machine wash and dry • Latex-free 	<ul style="list-style-type: none"> • Stage I or II Lymphedema • Post-surgical Edema
MEDAFIT ARM <ul style="list-style-type: none"> • 30–40mmHg of firm, gradient compression from hand to axilla • The inner WaveFoam liner creates high/low pressure ridges to promote lymph drainage • Outer shell fabric is very durable, and provides inelastic compression • Replaceable straps and foam liner extend the garment life • Separate MedaHand garment constructed with chipped foam • Unistrap can be wrapped around the hand piece for additional compression • Safe for day or night use • Easy care - machine wash and dry • Latex-free 	<ul style="list-style-type: none"> • Stage II or III lymphedema • Limb fills quickly when compression is removed • Fibrotic tissue is present • Higher compression levels needed

Measuring & Circumference					
	C	B1	B	A	D
RIGHT					
LEFT					

Products	Compression	Includes	Colors
CHIPSLEEVE	15–20 mmHg Graduated compression when used with the OverSleeve (included)	<ul style="list-style-type: none"> • CHIPSLEEVE • Oversleeve (5–10mmHg) • Pair of Cotton Liners 	Chipsleeve is only available in one color. Oversleeve Colors: <input type="checkbox"/> Black <input type="checkbox"/> Leopard <input type="checkbox"/> Navy <input type="checkbox"/> Tie Dye <input type="checkbox"/> Navy <input type="checkbox"/> Blue Horizon
COMPREFLEX ARM*	Provides mild/moderate compression	<ul style="list-style-type: none"> • COMPREFLEX ARM • Pair of Cotton Liners 	<input type="checkbox"/> Black <input type="checkbox"/> Beige
MEDAFIT ARM	30–40 mmHg	<ul style="list-style-type: none"> • MEDAFIT ARM • MEDAHAND • Pair of Cotton Liners • Donning Loop • Unistrap 	Products are only available in one color.

*For all available styles, colors, and compression, please see the SIGVARIS product catalog.

For indications and contraindications, please see prescribing information at sigvarisusa.com

Sigvaris and SIGVARIS GROUP are registered trademarks of SIGVARIS, CH-9014 Svt.Gallen/Switzerland, in many countries worldwide. © Copyright 2019 SIGVARIS, Inc.