



Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____

sigvaris

Foot Compression Wear

Circumference


Left Right

I _____

J _____

Length

K _____



K: Measure **medial** length from heel to 1st metatarsal head

Product Information

<input type="checkbox"/> Compreboot Lite Foot	<input type="checkbox"/> Right	Color: Black	Size:	Item #:	Quantity:
<input type="checkbox"/> Compreboot Lite Foot	<input type="checkbox"/> Left	Color: Black	Size:	Item #:	Quantity:
<input type="checkbox"/> Compreboot Standard Foot	<input type="checkbox"/> Right	Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Size:	Item #:	Quantity:
<input type="checkbox"/> Compreboot Standard Foot	<input type="checkbox"/> Left	Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Size:	Item #:	Quantity:
<input type="checkbox"/> Compreboot Plus Foot	<input type="checkbox"/> Right	Color: Black	Size:	Item #:	Quantity:
<input type="checkbox"/> Compreboot Plus Foot	<input type="checkbox"/> Left	Color: Black	Size:	Item #:	Quantity:
<input type="checkbox"/> Coolflex Standard Foot	<input type="checkbox"/> Right	Color: Black	Size:	Item #:	Quantity:
<input type="checkbox"/> Coolflex Standard Foot	<input type="checkbox"/> Left	Color: Black	Size:	Item #:	Quantity:

Measuring Instructions

Foot Size:

I: Measure circumference of ankle bend and heel
J: Measure circumference across metatarsal heads
K: Measure length from heel to 1st metatarsal head

Compreboot Lite Foot Size Chart (I-K)

(Available in Black Only)

	Small / Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long
I	28-36cm	28-36cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤33cm	≤33cm
K	14-18cm	19-23cm	20-22cm	25-27cm
Black	1803-BLR	1803-BLL	1805-BLR	1805-BLL

Compreboot Standard Foot Size Chart (I-K)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm
Black	1801-BTR	1801-BTL	1803-BTR	1803-BTL	1805-BTR	1805-BTL
Beige	1811-BTR	1811-BTL	1813-BTR	1813-BTL	1815-BTR	1815-BTL

Coolflex Standard Foot Size Chart (I-K)

(Available in Black Only)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-39cm	N/A	30-33cm	N/A	34-49cm	N/A
J	22-26cm	N/A	25-30cm	N/A	28-33cm	N/A
K	14-18cm	19-21cm	18-20cm	20-22cm	20-22cm	23-25cm
Black	2101-BTR	2101-BTL	2103-BTR	2103-BTL	2105-BTR	2105-BTL

Compreboot Plus Foot Size Chart (I-K)

(Available in Black Only)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-36cm	28-34cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-25cm	≤30cm	≤28cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm
Black	1801-BPR	1801-BPL	1803-BPR	1803-BPL	1805-BPR	1805-BPL