



SIGVARIS

NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASURING FORM & CHART

Patients Name:	Date:
Fitters Name:	Fitters Title: (example PT/OT/PTA)


SEND TO:
Luna Medical, Inc.
 1057 W Grand Ave, Ste. 1
 Chicago, IL 60642
lunamedical.com
 P: 800.380.4339
 F: 888.696.0299

MEDICALLY COMPLEX EDEMA

Step One

Measure the circumference at the fullest part of bust.


Bust Size



Step Two

Measure the circumference in inches directly under the bust (where a regular bra sits).

Band Size



Step Three

Determine the **Cup Size**:
 subtract the **Band** from the **Bust**.

Cup Size

2" or less = **A/B**
 2" - 4" = **C/D**
 4" - 6" = **DD/DDD**

Series & Features	Indications
COMPREVEST • 15-25mmHg of mild, supportive compression Caution: elastic band on sizes 2XLRG and above contain Natural Rubber Latex which may cause allergic reactions	• Stage I or II Lymphedema • Pitting edema, non-fibrotic • Post-surgical edema from mastectomy, lumpectomy or reconstructive procedures
• Features anterior hook adjustment for optimal fit • Straps are easily adjusted and feature durable, medical grade material • Interior pockets hold Chip Pad Bra or prosthetic insert • Easy care - machine wash and dry	