




PATIENT INFORMATION

| | | |
|---------|------|-----------------|
| Patient | Date | Number of Pairs |
|---------|------|-----------------|

| Style | | | | | <input type="checkbox"/> Closed Toe | <input type="checkbox"/> Open Toe |
|--------------------------|--------------------------|--------------------------|-------------------------------------------------------|--------------------------|-------------------------------------|-----------------------------------|
| | | | | | | |
| <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> L <input type="checkbox"/> R | <input type="checkbox"/> | | |
| Calf | Thigh | Pantyhose | Thigh w/ waist | Maternity/ Plus Sizes | | |

| Compression Level | Accessories |
|------------------------------------|-------------------------------------------------------------------|
| <input type="checkbox"/> 15–20mmHg | <input type="checkbox"/> Latex-Free Gloves |
| <input type="checkbox"/> 18–25mmHg | <input type="checkbox"/> Rubber Gloves |
| <input type="checkbox"/> 20–30mmHg | <input type="checkbox"/> Cotton Underliner |
| <input type="checkbox"/> 30–40mmHg | <input type="checkbox"/> SIGVARIS Doff N' Donner (Donning Device) |
| <input type="checkbox"/> 40–50mmHg | <input type="checkbox"/> S.O.S. (Slip On SIGVARIS Donning Device) |
| <input type="checkbox"/> 50–60mmHg | <input type="checkbox"/> Cornu-Thenard Extensor (Donning Device) |

LOWER EXTREMITY



CALF STYLE

- 1 Ankle Circumference at narrowest part of ankle, above ankle bone **(b)** _____
- 2 Calf Circumference at fullest part of the calf **(c)** _____
- 3 Calf Length from floor up to the bend of the knee (outside of calf) **(A–D)** _____

PANTYHOSE AND THIGH-HIGH STYLES

- (continue from steps 1–3)
- 4 Thigh Circumference at its widest circumference **(g)** _____
 - 5 Thigh Length from floor up to the gluteal fold **(A–G)** _____
 - 6 Hip Circumference at widest part of the hip **(f)** _____

Send to:

Luna Medical, Inc.

1816 W. Belmont Ave., Ste. 1
Chicago, IL 60657

P: 800.380.4339 F: 888.696.0299
lunamedical.com

WOMEN'S SHOE SIZES:

PETITE LEG 4–8, SHORT LEG 5.5–9.5, LONG LEG 8–12

MEN'S SHOE SIZES: SHORT LEG 7–11, LONG LEG 9.5–14

PLUS SIZES HIP CIRCUMFERENCE GUIDELINES:

SMALL: <56" (142CM), MEDIUM: <62" (157.5CM),
LARGE: <70" (178CM)

SIZING CHART

| ANKLE CIRCUMFERENCE | SMALL 7"–8.5" (18–21.5cm) | | | MEDIUM 8.5"–10" (21.5–25.5cm) | | | LARGE 10"–11.5" (25.5–29cm) | | EXTRA LARGE 11.5"–14" (29–36cm) | |
|----------------------------------|----------------------------------|--------------------------|-------------------------|-----------------------------------|--------------------------|-------------------------|----------------------------------|--------------------------|-------------------------------------------------------------------|------------------------|
| CALF CIRCUMFERENCE | SMALL 11"–15.5" (28–39cm) | | | MEDIUM 12.5"–17.5" (31.5–45cm) | | | LARGE 14"–20" (35.5–51cm) | | EXTRA LARGE 16"–24" (41–61cm) | |
| CALF LENGTH | PETITE* <14" (35.5cm) | SHORT 16" (41cm) | LONG >16" (41cm) | PETITE* <14" (35.5cm) | SHORT 16" (41cm) | LONG >16" (41cm) | SHORT 16" (41cm) | LONG >16" (41cm) | SHORT 16" (41cm) | LONG >16" (41cm) |
| THIGH CIRCUMFERENCE | SMALL 18–26" (45.5–66cm) | | | MEDIUM 19–28.5" (48–73cm) | | | LARGE 20.5–33" (52–84cm) | | *SP & MP SIZES ARE AVAILABLE IN CLOSED TOE ONLY, WITHOUT GRIP-TOP | |
| GUIDELINES FOR HIP CIRCUMFERENCE | SMALL 32.5–44.5" (83–112.5cm) | | | MEDIUM 34.5–50" (87–127cm) | | | LARGE 38.5–59.5" (98–150.5cm) | | | |
| LEG LENGTH: WOMEN & OPEN TOE | PETITE* 25.5" (65cm) | SHORT 28.5" (72cm) | LONG >28.5 (72cm) | PETITE* 25.5" (65cm) | SHORT 28.5" (72cm) | LONG >28.5 (72cm) | SHORT 28.5" (72cm) | LONG >28.5" (72cm) | | |
| LEG LENGTH: MEN | SHORT 30" (76cm) | | LONG >30 (76cm) | SHORT 30" (76cm) | | LONG >30 (76cm) | SHORT 30" (76cm) | LONG >30" (76cm) | | |
| SIZE | SP* | SS | SL | MP* | MS | ML | LS | LL | XS | XL |

For indications and contraindications, please see prescribing information sigvarisusa.com

SIGVARIS is a registered trademark of SIGVARIS AG, CH-9014 Svt. Gallen/Switzerland, in many countries © copyright 2015 SIGVARIS, Inc.