




# SIGVARIS

## NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASURING FORM & CHART


Patients Name:	Date:
Fitters Name:	Fitters Title: (example PT/OT/PTA)

**SEND TO:**  
**Luna Medical, Inc.**  
 1057 W Grand Ave, Ste. 1  
 Chicago, IL 60642  
**lunamedical.com**  
 P: 800.380.4339  
 F: 888.696.0299

### MEDICALLY COMPLEX EDEMA



Series & Features	Indications
<b>MEDAHAND</b> <ul style="list-style-type: none"> <li>Breathable / moisture wicking / anti-odor and soil-release fabric</li> <li>Open-cell foam chips gently massage tissue</li> <li>Conforming design ensures proper fit and gradient compression</li> </ul>	<ul style="list-style-type: none"> <li>Stage I or II Lymphedema</li> <li>Post-operative edema</li> <li>Other swelling issues</li> </ul>
<b>MEDAGLOVE</b> <ul style="list-style-type: none"> <li>Breathable / moisture wicking / anti-odor and soil-release fabric</li> <li>Open-cell foam chips gently massage tissue</li> <li>Conforming design with finger separation ensures proper fit</li> <li>Effective as a compression garment or bandage liner</li> <li>Gradient compression from fingers to wrist</li> </ul>	<ul style="list-style-type: none"> <li>Stage I or II Lymphedema</li> <li>Post-operative edema</li> <li>Other swelling issues</li> </ul>
<b>DORSAL POCKET GLOVE</b> <ul style="list-style-type: none"> <li>Dorsal pocket holds Chip Pad Dorsum or other therapeutic insert</li> <li>Full finger compression with doffing loops for easy use</li> <li>Palm pad provides soft, mild compression</li> <li>Ultra-thin profile for full hand function during use</li> </ul>	<ul style="list-style-type: none"> <li>Stage I or II Lymphedema</li> <li>Post-operative edema</li> <li>Other swelling issues</li> <li>Treats areas that are difficult to compress, including the full length of an individual fingers</li> </ul>



Series & Features	Indications
<b>NECK MANDIBLE PAD</b> <ul style="list-style-type: none"> <li>Provides 15–25mmHg of mild compression for the neck and lower face</li> <li>Conforming design and materials ensure a comfortable fit</li> <li>Straps are easily adjusted and feature durable, medical grade material</li> <li>Easy care - machine wash and dry</li> <li>Four-way stretch fabric contains moisture-wicking, odor control properties</li> <li>Latex-free</li> </ul>	<ul style="list-style-type: none"> <li>Stage I or II Lymphedema</li> <li>Post-surgical Edema</li> </ul>

**Head**  
 Measure circumference starting at the top of the head down to under the chin, back to top of head.

**Neck**  
 Measure circumference at largest part of neck.

For indications and contraindications, please see prescribing information at [sigvarisusa.com](http://sigvarisusa.com)

Sigvaris and SIGVARIS GROUP are registered trademarks of SIGVARIS, CH-9014 Svt.Gallen/Switzerland, in many countries worldwide. © Copyright 2019 SIGVARIS, Inc.

Rev.04.19