



Patient Last Name: _____ **Patient First Name:** _____
Fitter Last Name: _____ **Fitter First Name:** _____
Fitter Title: _____ **(example: PT/OT/PTA)**
Date: _____

sigvaris

Neck and Mandible Pad

Product Information

Product includes one Neck and Mandible Pad.

Neck & Mandible Color: Black
 Size: _____ Quantity: _____ Item #: _____

Measuring Instructions

Head:

Measure the circumference by encircling the head; beginning at the top, moving below the chin; and finishing at the top of the head.


Neck:

Measure circumference at largest part of neck.

Circumference

Head _____

Neck _____



Neck & Mandible Pad Sizing Chart (Available in Black Only)

	Small	Medium	Large	X Large
Head	46-56cm	52-62cm	56-66cm	58-68cm
Neck	28-33cm	33-38cm	38-43cm	43-58cm
Black	2231-NM	2232-NM	2233-NM	2234-NM