




# SIGVARIS

## NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASURING FORM & CHART

Patients Name:	Date:
Fitters Name:	Fitters Title: (example PT/OT/PTA)

**SEND TO:**  
**Luna Medical, Inc.**  
 1816 W. Belmont Ave., Ste. 1  
 Chicago, IL 60657  
**lunamedical.com**  
 P: 800.380.4339  
 F: 888.696.0299

### LOWER EXTREMITY



#### CALF STYLE

- 1 Ankle Circumference at narrowest part of ankle, above ankle bone (**b**) \_\_\_\_\_
- 2 Calf Circumference at fullest part of the calf (**c**) \_\_\_\_\_
- 3 Calf Length from floor up to the bend of the knee (outside of calf) (**A-D**) \_\_\_\_\_

#### PANTYHOSE AND THIGH-HIGH STYLES

(continue from steps 1-3)

- 4 Thigh Circumference at its widest circumference (**g**) \_\_\_\_\_
- 5 Thigh Length from floor up to the gluteal fold (**A-G**) \_\_\_\_\_
- 6 Hip Circumference at widest part of the hip (**f**) \_\_\_\_\_

#### Donning Accessories

- |   |  |
|---|--|
| <input type="checkbox"/> SIGVARIS ARION Sim-Slide (Open Toe Donning Device) | <input type="checkbox"/> SIGVARIS ARION Magnide On/Off (Closed Toe Donning Device) |
| <input type="checkbox"/> Latex-Free Donning Gloves                          | <input type="checkbox"/> SIMON (Donning Device)                                    |
| <input type="checkbox"/> Rubber Donning Gloves                              | <input type="checkbox"/> SIGVARIS Doff N' Donner (Donning Device)                  |
| <input type="checkbox"/> Cotton Underliner                                  | <input type="checkbox"/> S.O.S. (Slip On SIGVARIS Donning Device)                  |

Series	Features	Best Suited For
SELECT COMFORT  64% Nylon 36% Spandex  (contents may vary by style)  Latex-Free	From petite to plus, SELECT COMFORT offers a complete line of sizes designed to fit more than 99% of the population  Durable construction combined with a smooth, opaque look  Comfortable, stretchable material constructed from double covered yarn easily glides on and off  Plus-size pantyhose is a great alternative for maternity  Cotton-lined gusset on pantyhose style	Post-vein treatment (stripping, ligation, phlebectomy, sclero, EVLT, RF ablation)
COTTON  50% Nylon 25% Spandex 25% Supima Cotton  (contents may vary by style)  Latex-Free	Only Supima Cotton touches the skin to ensure comfort  Breathable and durable  Classic ribbed look  Ideal for everyday wear and for sensitive skin	Containment of venous edema and early stages of lymphedema  All-day cotton comfort in all types of climates  Ideal for dress or casual shoes  Ideal for long-distance travel

Products	Compression	Style	Colors
860 SELECT COMFORT*	<input type="checkbox"/> 20-30 mmHg <input type="checkbox"/> 30-40 mmHg	<b>Closed Toe:</b> <input type="checkbox"/> Calf <input type="checkbox"/> Petite calf <input type="checkbox"/> Calf plus <input type="checkbox"/> Calf w/grip-top <input type="checkbox"/> Thigh w/grip-top <input type="checkbox"/> Petite thigh-high w/grip-top  <input type="checkbox"/> Extra-large thigh w/grip-top <input type="checkbox"/> Pantyhose <input type="checkbox"/> Petite pantyhose <input type="checkbox"/> Plus-size Pantyhose  <b>Open Toe:</b> <input type="checkbox"/> Calf <input type="checkbox"/> Calf plus <input type="checkbox"/> Calf w/grip-top <input type="checkbox"/> Thigh w/grip-top <input type="checkbox"/> Extra-large thigh w/grip-top <input type="checkbox"/> Thigh-high w/ waist attachment	<input type="checkbox"/> Natural (33) <input type="checkbox"/> Suntan (36) <input type="checkbox"/> Crispa (66) <input type="checkbox"/> Black (99)
230 COTTON*	<input type="checkbox"/> 20-30 mmHg <input type="checkbox"/> 30-40 mmHg	<b>Closed Toe:</b> <input type="checkbox"/> Calf	<input type="checkbox"/> White (00) <input type="checkbox"/> Navy (10) <input type="checkbox"/> Black Mist (14)  <input type="checkbox"/> Crispa (66) <input type="checkbox"/> Chocolate (88) <input type="checkbox"/> Black (99)

\*For all available styles, colors, and compression, please see the SIGVARIS product catalog.

For indications and contraindications, please see prescribing information at [sigvarisusa.com](http://sigvarisusa.com)