




SIGVARIS

NON-CUSTOM CIRCULAR KNIT LEG STOCKING MEASURING FORM & CHART

Patients Name:	Date:
Fitters Name:	Fitters Title: (example PT/OT/PTA)

SEND TO:
Luna Medical, Inc.
 1816 W. Belmont Ave., Ste. 1
 Chicago, IL 60657
lunamedical.com
 P: 800.380.4339
 F: 888.696.0299

LOWER EXTREMITY



CALF STYLE

- 1 Ankle Circumference at narrowest part of ankle, above ankle bone **(b)** _____
- 2 Calf Circumference at fullest part of the calf **(c)** _____
- 3 Calf Length from floor up to the bend of the knee (outside of calf) **(A-D)** _____

PANTYHOSE AND THIGH-HIGH STYLES (continue from steps 1-3)

- 4 Thigh Circumference at its widest circumference **(g)** _____
- 5 Thigh Length from floor up to the gluteal fold **(A-G)** _____
- 6 Hip Circumference at widest part of the hip **(f)** _____

Donning Accessories

- | | |
|---|--|
| <input type="checkbox"/> SIGVARIS ARION Sim-Slide (Open Toe Donning Device) | <input type="checkbox"/> SIGVARIS ARION Magnide On/Off (Closed Toe Donning Device) |
| <input type="checkbox"/> Latex-Free Donning Gloves | <input type="checkbox"/> SIMON (Donning Device) |
| <input type="checkbox"/> Rubber Donning Gloves | <input type="checkbox"/> SIGVARIS Doff N' Donner (Donning Device) |
| <input type="checkbox"/> Cotton Underliner | <input type="checkbox"/> S.O.S. (Slip On SIGVARIS Donning Device) |

Series	Features	Best Suited For
SOFT OPAQUE	Comfortable tights with lasting microfibre softness	Patients who prefer fashionable options
79% Nylon 21% Spandex (content can vary slightly by style) Latex-Free	Advanced moisture management that keeps legs cool and dry Comfortable wide top band prevents knee-highs from rolling down Cotton-lined gusset on pantyhose style	Pre- and post-vein treatment with minimal edema Pregnancy Long-distance travel
EVERSHEER	Fashionable, sheer medical graduated compression hosiery	Patients who prefer fashionable options
64% Nylon 36% Spandex (content may vary by style) Latex-Free	Constructed from double-covered inlay yarns for softness against the skin Excellent breathability and moisture management for added comfort Cotton-lined gusset on pantyhose style	Pre- and post-vein treatment with minimal edema Ideal for business or dress

Products	Compression	Style	Colors
780 EVERSHEER* (for women only)	<input type="checkbox"/> 15-20 mmHg <input type="checkbox"/> 20-30 mmHg <input type="checkbox"/> 30-40 mmHg	Closed Toe: <input type="checkbox"/> Calf <input type="checkbox"/> Thigh-high w/grip-top <input type="checkbox"/> Pantyhose <input type="checkbox"/> Maternity Pantyhose Open Toe: <input type="checkbox"/> Calf <input type="checkbox"/> Thigh-high w/grip-top <input type="checkbox"/> Pantyhose <input type="checkbox"/> Maternity Pantyhose	<input type="checkbox"/> Dark Navy (08) <input type="checkbox"/> Suntan (36) <input type="checkbox"/> Natural (33) <input type="checkbox"/> Café (73) <input type="checkbox"/> Warm Sand (34) <input type="checkbox"/> Mocha (85) <input type="checkbox"/> Toasted Almond (37) <input type="checkbox"/> Nightshade (34) <input type="checkbox"/> Black (99)
840 SOFT OPAQUE* (for women only)	<input type="checkbox"/> 15-20 mmHg <input type="checkbox"/> 20-30 mmHg <input type="checkbox"/> 30-40 mmHg	Closed Toe: <input type="checkbox"/> Calf <input type="checkbox"/> Calf w/grip-top <input type="checkbox"/> Thigh-high w/grip-top <input type="checkbox"/> Pantyhose <input type="checkbox"/> Maternity Pantyhose Open Toe: <input type="checkbox"/> Calf <input type="checkbox"/> Thigh-high w/grip-top <input type="checkbox"/> Pantyhose	<input type="checkbox"/> Graphite (91) <input type="checkbox"/> Black (99)

*For all available styles, colors, and compression, please see the SIGVARIS product catalog.

For indications and contraindications, please see prescribing information at sigvarisusa.com

SIGVARIS is a registered trademark of SIGVARIS, CH-9014 Svt. Gallen/Switzerland, in many countries worldwide. © Copyright 2018 SIGVARIS, Inc. Rev.12.18