




Account Number: 4057807

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example PT/OT/PTA)  
 Date: \_\_\_\_\_

## FarrowHybrid™ Foot Compression Order Form

Business: \_\_\_\_\_ Patient: \_\_\_\_\_ PO: \_\_\_\_\_  
 Person Measuring: \_\_\_\_\_ Phone/Fax: \_\_\_\_\_ / \_\_\_\_\_ Date: \_\_\_\_\_



**FarrowHybrid™ AD Foot Compression** is half compression stocking (over the foot and ankle) and half sock (over the leg) and may be used in place of our footpieces with our legpiece. Easy to apply and remove, they are comfortable and readily fit in patients' shoes.

FarrowHybrid™ AD		SKU	Small	Medium	Large	XLarge
Male Shoe Size			3 – 6	7 – 9	10 – 12	13 – 15
Female Shoe Size			4 – 7	8 – 10	11 – 13	14 – 16
<b>FarrowHybrid™ ADI</b>	FLFH-ADI	Qty		Qty		Qty
<b>FarrowHybrid™ ADI Wide</b>	FLFH-ADIW	Qty		Qty		Qty
<b>FarrowHybrid™ ADIIS</b>	FLFH-ADIIS	Qty		Qty		Qty
<b>FarrowHybrid™ ADIIS Wide</b>	FLFH-ADIISW	Qty		Qty		Qty

FarrowHybrid™ Foot Compression Size Selector		
	20 – 30mmHg	30 – 40mmHg
A@Circumference 20 – 27cm	FLFH-ADI	FLFH-ADIIS
A@Circumference 28 – 35cm	FLFH-ADIW	FLFH-ADIISW

FarrowHybrids cover a wide range of circumferences in standard and wide sizes. The maximum leg circumference must be 60 cm or less. FarrowHybrid™ ADI and ADIW provide Class I (20-30mmHg) compression on most swollen feet and FarrowHybrid™ ADIIS and ADIISW provide Class II (30-40mmHg) compression. The FarrowHybrid™ ADIIS foot compression CONTAINS SILVER and should not be used on patients allergic to silver.

### FarrowHybrid™ AD Foot Compression

The **FarrowHybrid™ AD Foot Compression** garments are easily stretchable in the legging, which makes them easy to apply and remove, while being low profile enough to fit into the patient's current shoes.

### FarrowHybrid™ Measuring Instructions

Measure the circumference at the midfoot and widest calf. Midfoot circumference ranges are between 20 – 27cm for **standard** sizes and 28 – 35cm for the **wide** sizes. The widest calf circumference must be less than 60 cm and shoe size must be between 3.5-16 US. Indicate in the blank white box the number of FarrowHybrid's™ desired. FarrowHybrid™ AD are sold as pairs and are designed for use with Farrow brand leg garments.




Account Number: 4057807

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example PT/OT/PTA)  
 Date: \_\_\_\_\_

## FarrowWrap® Liners Order Form

Business: \_\_\_\_\_ Patient: \_\_\_\_\_ PO: \_\_\_\_\_  
 Person Measuring: \_\_\_\_\_ Phone/Fax: \_\_\_\_\_ / \_\_\_\_\_ Date: \_\_\_\_\_



Farrow Silver Liners		Sizes			
		Small (Flesh Color)		Large (Silver Color)	
Calf (C) or Arm Circumference (F@)		15 – 58 cm		52 – 70 cm	
Silver AD Liner	FLS_-AD_	Qty		Qty	
Silver Arm Liner	FLSF-A2	Qty			
Silver Arm Hand Liner	FLSF-AH2	Qty			

TG® Soft is used on patients who are allergic or adversely react to silver, do not fit into a silver liner, who need a little more cushioning than our silver liners, or who desire an open toe.

TG® Soft Liner		Sizes					
		Small		Medium		Large	
Extremity Circumference		< 40 cm		40 – 70 cm		70 – 120 cm	
Arm (F@)	FLTG-ARM	Qty		Qty		Qty	
AD Knee High (C)	FLTG-AD	Qty		Qty		Qty	
Are the AD liners a SUBSTITUTE to our standard liners? (Only if you circle "YES" here will you not be charged.)						YES / NO	
AG Thigh High +G,	FLTG-AG	Qty		Qty		Qty	