




Account Number: 1014233

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Fitter Title: \_\_\_\_\_ (example PT/OT/PTA)  
 Date: \_\_\_\_\_

## Size Charts For Dynamic Max Compression Stockings 20-30 mmHg and 30-40 mmHg



### 20-30, 30-40 mmHg Dynamic Max Compression Stockings 3511 / 3512

|           | III                         | IV                          | V                           |  |
|-----------|-----------------------------|-----------------------------|-----------------------------|--|
| <b>cG</b> | 60-82 cm<br>23 3/4"-32 1/4" | 62-87 cm<br>24 1/2"-34 1/4" | 77-94 cm<br>30 1/4"-37"     | length A-D<br>A-D Regular - 15 3/4"-18 1/4" / 40-46 cm<br>A-D Short - 13-15 3/4" / 34-40 cm      |
|           | 49-55 cm<br>19 1/4"-21 3/4" | 53-60 cm<br>21 1/4"-23 1/2" | 55-66 cm<br>21 3/4"-26"     | length A-G<br>A-G Regular - 28 1/4"-32 3/4" / 72-83 cm<br>A-G Short - 24 3/4"-28 1/4" / 63-72 cm |
| <b>cC</b> | 24-27 cm<br>9 1/2"-10 3/4"  | 27-32 cm<br>10 3/4"-12 3/4" | 32-37 cm<br>12 3/4"-14 3/4" |  |
|           | <b>cB</b>                   |                             |                             |  |

PCSZ-01-47