



Account Number: 4057807



Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example PT/OT/PTA)
 Date: _____

Upper Extremity Order Form

| Quantity/Class | CCL 1 18.5-21.5 mmHg* | CCL 2 25.1-31.5 mmHg* | CCL 3 36.1-45.5 mmHg* | CCL 3F 36.1-45.5 mmHg* | CCL 4 58-70 mmHg* | CCL 4S 60-90 mmHg* |
|---------------------------------------|--------------------------|--------------------------|--------------------------|---------------------------|----------------------|-----------------------|
| Left | | | | | | |
| Right | | | | | | |
| Pressure panty Elvarex (Body Bandage) | | | | | | |

Basic styles

- AD Knee
- AF Mid-Thigh
- AG Thigh
- AG-T (AG with chapstyle)
 - Piece
 - Pair
- AG-HT 1½ Leg Pantyhose
- AT Pantyhose

Options

- Profile/Ankle pad
- T-Heel (Class 2-3 forte only)

Comments

* Design Pressure

CAUTION: This product contains natural rubber latex which may cause allergic reactions.

Variations

- B*G
- BG
- FG (leg extension)
- GT Biker Shorts
- B1G-T Capri

Color

- Beige
- Black
- Cocoa
- Navy
- Grey


Special options

- Knee Comfort Zone (ccl 2 & ccl 3)
- Silicone dotted band 2.5cm:
 - Top
 - Inside
 - Inside ¾
- Silicone dotted band 5cm:
 - Top
 - Inside
 - Inside ¾
- Zipper B to D only:
 - Inside (medial)
 - Outside (Lateral)
- Zipper E to G only:
 - Inside (medial)
 - Outside (Lateral)
- Pocket
 - Instep
 - Back of Knee
- Fly for Men
- Open Pubis
- Adjustable Waistband
- Top Comfort Zone

Profile measuring points


Point P1: is the horizontal circumference from the inner to the outer malleolus.

Point P2: is the "stirrup" or vertical circumference from the inner to the outer malleolus.



Circumference

Left Waist Right



| Circum. (c) | Length (ℓ) | Length (ℓ) |
|-------------|------------|------------|
| cT | K2-T | ℓT |
| cH | K1-T | ℓH |

_____ cG _____ cG _____ ℓG _____ ℓG _____

_____ cF _____ cF _____ ℓF _____ ℓF _____

_____ cE _____ cE _____ ℓE _____ ℓE _____

_____ cD _____ cD _____ ℓD _____ ℓD _____

_____ cC _____ cC _____ ℓC _____ ℓC _____

_____ cB' _____ cB' _____ ℓB' _____ ℓB' _____

_____ cB _____ cB _____ ℓB _____ ℓB _____

Straight Open Toe
Lateral Length _____ cm

Straight Closed Toe
Total Foot Length _____ cm

Slant Open Toe
Medial Length _____ cm
Lateral Length _____ cm

Slant Closed Toe
Medial Length _____ cm
Lateral Length _____ cm
Total Foot Length _____ cm

Profile measurement points

Left leg P1 _____ P2 _____

Right leg P1 _____ P2 _____

