



Comfort, Health and Style!

A photograph of a woman and a young girl sitting on a light-colored rug in a living room. The woman, on the left, has long dark hair and is wearing a blue top and grey pants. She is wearing a white, textured compression sleeve on her right arm. She is holding an open book and pointing to a page. The girl, on the right, has long dark hair and is wearing a floral dress. She is also looking at the book. In the background, there is a grey sofa with a stuffed rabbit on it and a wooden coffee table with some books on it.

JOBST® RELAX
REDEFINING LYMPHEDEMA NIGHT CARE

THERAPIES. HAND IN HAND.



JOBST® RELAX – DESIGNED FOR MORE THAN JUST SLEEPING

JOBST® Relax is an innovative, new product designed to complement daytime compression garments and offers an alternative to nighttime compression solutions such as compression bandages. Night care is an often neglected, but critical, part of lymphedema management, and **JOBST® Relax** was specifically designed to address patients' complaints with current nighttime garments. **JOBST® Relax** provides the edema control patients need, combined with the look and comfort that patients want. **JOBST® Relax** is intended for mild to moderate lymphedema patients in the maintenance phase, who are looking for a garment that will better suit their lifestyle. It is ideal for those who need containment, combined with comfort and style.



Monofilament spacer layer

- Flexible, lightweight, and slim-line
- Provides cushioning and promotes airflow
 - ▶ Comfortable wearing experience

Soft textured surface

- Luxurious feel, inside and out
- Creates a micro-massaging effect
 - ▶ Stimulation of lymph flow

Customized, flat-knitting technology

- Made to patient's unique dimensions
- Provides consistent gradient compression
- Anatomically shaped
 - ▶ Stays in place and contains swelling



The **integrated zipper** option with stay-down puller facilitates donning and doffing. The fabric underneath the zipper protects the skin and promotes uniform compression.

COOLMAX®

Coolmax® yarns and advanced knitting technique

- Moisture wicking properties
- Increased breathability
 - ▶ Cooling effect

THE SLIM LINE, FLEXIBLE, AND LIGHTWEIGHT DESIGN OF JOBST® RELAX:



Adjusts to almost any sleeping position.







Ensures the freedom of movement that patients want. This makes **JOBST® Relax** ideal for relaxing before bedtime.



Means the garment can be folded to fit comfortably into a purse or travel bag, making **JOBST® Relax** the perfect travel companion.

Greater comfort and freedom encourage patient compliance.

JOBST® Relax Compression classes / Styles	AG1 	CG1 	AD 	AG 
CCL1 (15–20 mmHg)	●	●	●	●
CCL2 (20–30 mmHg)			●	●
Colors				
Beige	○	○	○	○
Rosé	○	○	○	○
<ul style="list-style-type: none"> • The leg sleeves have an integrated durable anti-slip velour pad on the sole to prevent from slipping and sliding • Exclusively delivered in its very own satin bag to protect and store the garment. Perfect for traveling. • Zipper option available on all styles 				