



THUASNE

When it comes to insurance or private pay, we've got you covered!

National in-network insurance contracts for compression products

Elastic and in-elastic support are covered by insurance



We're so swell – you don't have to be.
Compression wear it counts.

www.lunamedical.com
[800] 380-4339

TAKE CARE DIFFERENTLY

We believe in challenging the status quo, in everything we do. Through our unique approach, we are able to manufacture products that are efficient, easy to use, comfortable and effective. We are changing the industry by taking care, differently.


A French Medical Company in the service of health since 1847

For over a 150 years Thuasne has provided medical solutions that empower patients to take an active role in their own health.

Care – Designing and distributing customized solutions built for today's lifestyle

Empowerment – Providing patients with the products and resources to help manage their own health

Wellness – Working with patients to optimize mobility and quality of life.



A leader in industrial skills, off-the-shelf solutions and custom devices

Expertise built over decades in the manufacture of yarn to the production of finished goods

Use of diverse technologies including wrapping, weaving, knitting, coating, welding, assembly, silicone, molding and shaping of specialty metals

A capacity to produce off-the-shelf, made-to-measure and custom-to-cast products with a high level of precision

MOBIDERM

A unique technology that utilizes pressure difference to improve the flow of lymphatic fluid and optimize drainage.

EFFECTIVE

- rapid reduction of the volume of the edema⁽²⁾

PRACTICAL


- the light and flexible design allows the patient to maintain freedom of movement⁽⁴⁾

COMFORTABLE

- soft skin contact
- conforms easily to the patient's anatomy

THE TECHNOLOGY

A pattern of foam cubes layered between two non-woven bandages. When applied to the patient, Mobiderm creates a pressure difference between the supported area of the blocks and the non-supported area surrounding. This pressure difference facilitates the flow of lymphatic fluid and encourages drainage efficacy⁽¹⁾.



Two sizes of foam blocks to adapt to the various needs of the affected location



SMALL BLOCKS

- More effective in highly induced areas
- Better for small morphologies or body areas: pediatric patients, hands/foot, breast, and more.




BIG BLOCKS

- Better adapted along the limb

MOBIDERM autofit garments are recommend for patients with all stages of lymphedema, dysmorphic or regular morphology and progressive edema.

Reduction Phase



MOBIDERM
Pads



MOBIDERM
Bandages



MOBIDERM
Kit

Maintenance Phase



MOBIDERM autofit
Standard Garments



MOBIDERM autofit
Adjustable Garments

(1) Todd - Mobiderm Autofit: an adjustable sleeve that enables patients to self-manage lymphoedema - Chronic oedema - April 2018 - Page 6
 (2) Quéré - Prospective multicenter observational study of lymphoedema therapy: POLIT study - Vascular Disease Journal - 2014 - Pages 4 and 5
 (3) Only in following countries: Algeria, Austria, Belgium, Canada, China, Czech Republic, Denmark, Finland, France, Germany, Hungary, India, Morocco, Italy, Netherlands, Poland, Portugal, Romania, Russian Federation, Slovakia, Slovenia, Spain, Sweden, Switzerland, Tunisia, United-Kingdom
 (4) Mestre - Interest of an auto-adjustable night-time compression sleeve (MOBIDERM Autofit) in maintenance phase of upper limb lymphoedema: the MARILYN pilot RCT - Support Care Cancer (25:2455-2462) - 2017 - Page 7

Mobiderm pads & bandages - Indications: Drainage of edema (venous / lymphatic / post-traumatic). Contraindications: Do not use on a skin surface that is infected or very swollen.
 Mobiderm Autofit - Indications: Maintenance of volume reduction in lymphedema. Contraindications: Do not use on broken skin or on a skin surface that is infected or very swollen. Do not use in the following cases: known allergies to the materials used, advanced diabetic microangiopathy, stage III or IV arterial disease of the limbs, decompensated heart failure, exudative skin condition on the limb treated. Mobiderm standard garments - Indications: Maintenance of volume reduction in lymphoedema. Contraindications: Do not use on a skin surface that is infected or very swollen.
 Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region. All the medical devices mentioned on this document are CE marked according to the European council directive 93/42/EEC and its relatives or the Regulation 2017/745 on medical devices, unless specifically identified as "NOT CE MARKED". Products mentioned in this document are CE class I devices. Please contact Thuasne should you need any additional information on devices classification.
 Please read carefully the instructions for use, indications and contraindications of the product.

Your Thuasne contact:

Ryan Hansen
Business Development
Manager - Compression



Tel: 800-432-3466
Mail: rhansen@thuasneusa.com