



## Seamless Options

Models:

### 2300 Seamless Glove & Gauntlets



- Seamless
- Smooth appearance

Compression	(2301) 20-30 (2302) 30-40 mmHg
Styles	AC Gauntlet (2301/02) ACFS Glove (2301)
Regular Sizes	S, M, L, XL
Colors	Beige, Black

### Soft Seamless Glove & Gauntlets



- Added durability in the palm
- Channeled fingers for greater ease in donning

Compression	(2000) 15-20, (2001) 20-30, (2002) 30-40 mmHg
Styles	AC Gauntlet (2000/01/02) ACFS Glove (2000/01)
Regular Sizes	S, M, L, XL specify right or left
Colors	Beige, Black, Chocolate, Cinnamon, Chestnut, White Pink, Trend Colors & Prints


### 2300 Seamless Foot Portion



- Anatomical shape, providing a precise fit
- 3D comfort zones in sensitive areas, providing high wearing comfort.

Compression	(2300) 15-20, (2301) 20-30 mmHg
Styles	FP (2300/01)
Regular Sizes	S, M, L, XL specify right or left
Colors	Beige, Black

**New**  
2300 Foot Portion  
Soft Seamless  
Glove / Gauntlet



## Size chart

### Juzo 2000 & 2300 Seamless Glove or Gauntlet

Item # 2000/01/02 & 2301/02 AC & ACFS		S	M	L	XL
	A	17.5 – 19.5 cm 6 7/8" - 7 5/8"	19.5 – 21.0 cm 7 5/8" - 8 1/4"	21.0 – 23.0 cm 8 1/4" - 9 1/8"	23.0 – 25.5 cm 9 1/8" - 10 1/8"
	C	14.5 – 16.5 cm 5 5/8" - 6 1/2"	16.5 – 18.0 cm 6 1/2" - 7"	18.0 – 20.0 cm 7" - 7 7/8"	20.0 – 22.0 cm 7 7/8" - 8 5/8"
PCSZ-01-62					

### Juzo 2300 Seamless Foot Portion (available soon)

Item # 2300, 2301 FP		S	M	L	XL
	A'	19.5 – 22.5 cm 7 5/8" - 8 7/8"	22.5 – 25.0 cm 8 7/8" - 9 3/4"	25.0 – 27.5 cm 9 3/4" - 10 7/8"	27.5 – 30.0 cm 10 7/8" - 11 7/8"
	A	19.5 – 22.5 cm 7 5/8" - 8 7/8"	22.5 – 25.0 cm 8 7/8" - 9 3/4"	25.0 – 27.5 cm 9 3/4" - 10 7/8"	27.5 – 30.0 cm 10 7/8" - 11 7/8"
PCSZ-01-61					

**Juzo Seamless Compression for Lymphology**  
Seamless Glove, Gauntlet & Foot Portion



## Juzo Seamless

### Seamless compression for lymphology

Juzo Seamless Products are the ideal solution for early stages of Lymphedema therapy.

#### Innovative fabric

The innovative knitting technology enables seamless manufacturing of compression garments. The soft 3-dimensional fabric conforms perfectly to anatomical features and produces well-defined gradient compression. The innovative seamless fabric is breathable and is comfortable against the skin. Thanks to its soft characteristics, the fabric offers high wearing comfort and makes donning and doffing easy.

#### Optimal fit

The lengthwise knitted structure guarantees easy donning while allowing efficient compression modeling and additional stability. The top borders are especially soft and pressure-minimized. This prevents the borders from rolling up and avoiding any constrictions or pressure points.

#### Features

- light weight material
- breathable and comfortable against the skin
- smooth and seamless
- easy donning and doffing

#### Indications

- early stages of edema and swelling
- mild to moderate lymphedema



seamless knitting



standard sizes

## New Juzo Soft Seamless Glove & Gauntlet

The anatomical shape of fingers and thumb, along with the enclosure of the thumb webbing ensure a precise and comfortable fit. The 3D comfort zone in the sensitive area between the thumb and index finger prevents the fabric from creasing and protects the skin from irritation. It is designed not to restrict the thumb movement in any way.



## New Juzo 2300 Seamless Foot Portion (available soon)

The knitted length of the toes ensures additional stability and easy donning of the toe portion. The 3D comfort zone in the sensitive area between the big toe and the ball of the foot has a precise anatomical shape. This allows a precise fit and prevents the fabric from creasing, protecting the skin from irritation and providing a high wearing comfort.

