




Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____



Mobiderm Autofit - Lower Extremity

SOCK							
Size	Ankle circumference (cB) <i>(at the narrowest point)</i>	Calf circumference (cC) <i>(at the widest point)</i>	Below knee circumference (cD) <i>(on the head of the fibula)</i>	Reference (right/left)			
				Short ℓD = 38 - 40 cm	Normal ℓD = 40 - 42 cm	Long ℓD = 42 - 44 cm	Xtra Long ℓD = 44 - 46 cm
1	16 - 24 cm	24 - 35 cm	25 - 36 cm	37340120100231	37340120100232	37340120100233	37340120100234
2	20 - 28 cm	32 - 43 cm	32 - 43 cm	37340120200231	37340120200232	37340120200233	37340120200234
3	24 - 32 cm	40 - 51 cm	39 - 50 cm	37340120300231	37340120300232	37340120300233	37340120300234
4	28 - 36 cm	48 - 59 cm	46 - 57 cm	37340120400231	37340120400232	37340120400233	37340120400234

STOCKING							
Size	Ankle circumference (cB) <i>(at the narrowest point)</i>	Calf circumference (cC) <i>(at the widest point)</i>	Thigh circumference (cG) <i>(at the widest point)</i>	Reference (right/left)			
				Short ℓG = 67 - 72 cm	Normal ℓG = 72 - 77 cm	Long ℓG = 77 - 82 cm	Xtra Long ℓG = 82 - 87 cm
1	16 - 24 cm	24 - 35 cm	25 - 36 cm	37350120100231	37350120100232	37350120100233	37350120100234
2	20 - 28 cm	32 - 43 cm	32 - 43 cm	37350120200231	37350120200232	37350120200233	37350120200234
3	24 - 32 cm	40 - 51 cm	39 - 50 cm	37350120300231	37350120300232	37350120300233	37350120300234
4	28 - 36 cm	48 - 59 cm	46 - 57 cm	37350120400231	37350120400232	37350120400233	37350120400234

ℓD: floor-to-knee height (measurement taken below the bend of the knee, back side)
 ℓG: floor-to-thigh height (measurement at the gluteal fold)