




Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Fitter Title: _____ (example: PT/OT/PTA)
 Date: _____



Mobiderm
 Standard-
 Upper
 Extemity

ARMSLEEVE					
Size	Wrist circumference (C)	Elbow circumference (E)	Axillary arm circumference (G)	Reference (right / left)	
				Long (45 - 50 cm)	Normal (40 - 45 cm)
1	14 - 16 cm	22 - 25 cm	25 - 29 cm	37300220100233	37300220100232
1+		24 - 27 cm	28 - 32 cm	37300240100233	37300240100232
2	16 - 18 cm	25 - 28 cm	29 - 33 cm	37300220200233	37300220200232
2+		27 - 30 cm	31 - 35 cm	37300240200233	37300240200232
3	18 - 20 cm	28 - 31 cm	32 - 36 cm	37300220300233	37300220300232
3+		30 - 33 cm	35 - 39 cm	37300240300233	37300240300232
4	20 - 22 cm	31 - 34 cm	35 - 39 cm	37300220400233	37300220400232
4+		33 - 36 cm	38 - 42 cm	37300240400233	37300240400232
5	22 - 24 cm	34 - 37 cm	38 - 42 cm	37300220500233	37300220500232
5+		36 - 39 cm	41 - 45 cm	37300240500233	37300240500232