



Compression for
those who love to
**Move
Freely**

Kimberly,
Secondary lymphedema patient

JOBST® Confidence

The Next Generation of Made-to-Measure
Flat-Knit Compression Garments

**NEW
& UNIQUE**


Introducing a New Era of Comfort

JOBST® Confidence helps your patients feel unrestricted, with a one-of-a-kind fit garment that's as unique as they are. The new generation of flat-knit technology used in JOBST® Confidence allows for a precise fit to the real anatomical body shapes. This personally refined fit not only provides comfort, but also supports your patients with a greater degree of movement.

The new Contour Fit technology

JOBST® Confidence has been expertly developed with a new knitting technology that, unlike traditional flat knit, allows the fit of the garment to adapt in more positions.

This supports a true personalization of garment fit. Specially designed knitting machines were created for the JOBST® Confidence line to deliver 3D contour fit.



JOBST® Confidence is designed to:

- provide more comfort
- support easy donning
- be more skin-friendly
- deliver improved climate control
- provide better movement
- offer a range of colors and design options

JOBST® Confidence



Design options available in **Deco line** and **patient initials**




3D innovative knitting method allows to **conform to individual body shapes**



Soft & Firm

Conforming fabric that also offers **firm support**



Moisture Wicking

Moisture wicking with double-layer construction for **optimal comfort**



Easy Donning

Flexible, conforming fabric to make putting on and taking off easier



Skin-friendly fabric that meets the needs for **sensitive skin, pediatric, elderly, and lipedema** patients


Garment shown includes knee comfort, Deco line and patient initials.



Next-generation technology


Unique Moisture Wicking Yarns for breathability and eliminates heat build-up

Only JOBST® Confidence utilizes a high-performance double-layer fabric to wick away moisture from the skin and towards the garment's outer layer. This is designed to keep the skin dry and leave patients feeling comfortable all day long.




Softer yarns to support a comfortable feeling

JOBST® Confidence is created with supple yarns that have been used with a double-layered knitting construction. This provides the flexibility your patients want when wearing compression.



Firm support that adapts to each individual

The innovation in JOBST® Confidence garments means that for the first time ever, the properties of the yarn adapt to the real size of any patient. This always delivers the optimal compression level according to the RAL standard. The RAL standard is used as a basic requirement for certifying medical compression stockings and guarantees the highest quality by an independent group.




JOBST® Confidence - Innovative 3D Technology for the freedom of flexible movement

Colors



Styles




Compression Classes

Upper Limb	Lower Limb
CCL 1 ▶ 15 - 21 mmHg	CCL 1 ▶ 18 - 21 mmHg
CCL 2 ▶ 23 - 32 mmHg	CCL 2 ▶ 23 - 32 mmHg
	CCL 3 ▶ 34 - 46 mmHg



with SoftFit technology

New Design Options



Deco Line: The deco line is positioned on the front of the garment.
Patient Initials: Initials are located under the top band of the garment.

Comfort Options



Elbow Functional Zone, Slant Toe, Ankle Functional Zone

