

JOBST® Elvarex® Portfolio – Offers a complete

JOBST® Elvarex® – with Natural Rubber Core

JOBST® Elvarex Key Features

- Flat-knit structure bridges skin folds instead of working into them
- Latex core provides excellent containment and allows garment to return to original shape for extended wear
- Excellent breathability for comfort and patient compliance
- Seven compression classes for optimal therapeutic outcomes
- Vast product assortment covers a diverse range of specifications and requirements of patients
- SPF 80 UV protection



Basic Styles



Compression classes



Extremities / Colors

Upper	1-2	Beige	Black	Navy	Honey	Cranberry	Cherry	Caramel		
	2-Forte	Beige	Black	Navy	Honey	Cranberry	Cherry			
Lower	1-4 Super	Beige	Black	Navy	Cocoa	Cranberry	Grey	Henna**	Graphite**	Denim**

* Not available in Toe Cap


JOBST® Elvarex® Plus®

Seamless 3D Knitting Technology – Pairs with Confidence Line


JOBST® Elvarex® Plus® Key Features

- 3D technology to closely conform to hands and feet
- The seamless helps eliminate potential irritation
- Elvarex Plus styles pairs with the Confidence line

Basic Styles



Compression classes



Type / Colors

Hand	Beige	Black	Navy	Honey	Cranberry	Caramel	Cherry	Red
Foot	Beige	Black	Navy	Cocoa	Cranberry	Grey	Red	

te range that's as individual as your patients

JOBST® Elvarex® Soft

Softer, Latex-Free Garment

**LATEX
FREE
FORMULA**

JOBST® Elvarex® Soft Key Features

- Specialized seam yarns for a smoothhand gentle seam against the skin
- Well suited for pediatric or elderly patients and patients with lipedema or sensitive skin
- Excellent anatomical fit
- Latex-free garment



Basic Styles



Compression classes



Extremities / Colors

Upper CCL 1-2	● Beige	● Black	● Navy	● Cocoa	● Grey	● Cranberry	● Cherry	● Honey
Lower CCL 1	● Beige	● Black	● Navy	● Cocoa	● Grey	● Cranberry	● Cherry	

JOBST® Elvarex® Soft Seamless

Soft and Seamless for Hands and Feet


**LATEX
FREE
FORMULA**

JOBST® Elvarex® Soft Seamless Key Features

- Soft, durable, and breathable flat-knitted fabric
- Flat edge finish at fingers and toes decreases bulk and prevents curling
- Easy to don, comfortable to wear
- Latex-free garment to pair with Confidence line



Basic Styles



Compression classes



Extremities / Colors

Upper CCL 1-2	● Beige	● Black	● Navy	● Cocoa	● Grey	● Cranberry	● Cherry	● Honey
Lower CCL 1-3	● Beige	● Black	● Navy	● Cocoa	● Grey	● Cranberry	● Cherry	