

Power Socks

Models:

Power Comfort



- Men's Casual Sock
- Stylish, broad ribbed design
- Padded sole

Compression	(2600) 15-20 mmHg (2601) 20-30 mmHg
Styles	AD
Regular Sizes	2600 - M, L, XL 2601 - S, M, L, XL, 2XL
Colors	Black, White, Heather Gray Retro Stripes - Black, White, Heather Gray

Power Lite



- Light weight
- Moisture wicking & fast drying using *DriRelease*
- Easy to put on



Compression	(2620) 15-20 mmHg
Styles	AD
Regular Sizes	S, M, L, XL
Colors	Black, White

Power Vibe



- Trendy patterns
 - Moisture wicking
 - Easy to put on
- | | |
|---------------|---|
| Compression | (2650) 15-20 mmHg |
| Styles | AD |
| Regular Sizes | S, M, L, XL |
| Colors | Deep Black, Groovy Zebra,
Cool Dot, Super Stripe |

Better Sizing & Fit

- Easy to use 15-20 mmHg size chart
- 70% fewer skus (15-20 mmHg)
- Retail friendly
- Power Lite 15-20 mmHg available Summer 2021
- Power Lite & Power Vibe 20-30 mmHg available Fall 2021

Juzo USA
3690 Zorn Drive
Cuyahoga Falls, OH 44223
www.juzousa.com
support@juzousa.com

Juzo Canada Ltd.
1100 Burloak Dr. Ste. 300
Burlington, ON L7L 6B2
www.juzo.ca
support@juzo.ca



POWER Comfort



POWER Lite



POWER Vibe



POWER SOCKS
Socks That Fit Your Lifestyle

POWER COMFORT

Men’s Comfort Compression Sock

Juzo **Power Comfort Socks** for men designed with thick padded sole, providing all day cushioning.



Size Chart for Juzo 2601 Power Comfort 20-30 mmHg						
	S	M	L	XL	2XL	Length
Calf	29 - 38 cm 11 1/2" - 15"	34 - 43 cm 13 1/4" - 17"	37 - 49 cm 14 1/2" - 19 1/4"	41 - 53 cm 16 1/4" - 21"	46 - 58 cm 18" - 23"	AD Standard 13-15 3/4" / 34-40 cm
Ankle	18 - 21 cm 7" - 8 1/4"	21 - 24 cm 8 1/4" - 9 1/2"	24 - 27 cm 9 1/2" - 10 3/4"	27 - 31 cm 10 3/4" - 12 1/4"	31 - 35 cm 12 1/4" - 13 3/4"	AD Tall 15 3/4" - 18 1/4" / 40-46 cm

PCSZ-01-63

Size Chart for Juzo 2600 Power Comfort 15-20 mmHg				
	M	L	XL	
Men’s Shoe Size	6.5 - 8	8.5 - 10.5	11 - 14	Choose the appropriate size based off your shoe size.
Women’s Shoe Size	7.5 - 9.5	10 - 12	12 +	Verify that your ankle & calf measurements fall within the given ranges for that size.
Calf Circumference	32 - 44 cm 12 1/2" - 17 1/4"	34 - 46 cm 13 1/4" - 18"	37 - 49 cm 14 1/2" - 19 1/4"	
Ankle Circumference	19.5 - 23.5 cm 7 3/4" - 9 1/4"	20.5 - 24.5 cm 8" - 9 3/4"	21.5 - 25.5 cm 8 1/2" - 10"	Use shoe size with ankle/calf measurements for best fit.

PCSZ-01-64

POWER LITE

Unisex, Lightweight Compression Sock

Juzo **POWER Lite Socks** designed with specially blended *DriRelease®* fabric, keeps sensitive skin dry, refreshed and provides a gentle touch.



Power Lite 20-30 mmHg available Fall 2021

Size Chart for 2620 Juzo Power Lite 15-20 mmHg				
	S	M	L	XL
Men’s Shoe Size	6	6.5 - 8	8.5 - 10.5	11 - 14
Women’s Shoe Size	5 - 7	7.5 - 9.5	10 - 12	12 +
Calf Circumference	30 - 41 cm 11 3/4" - 16"	32 - 44 cm 12 1/2" - 17 1/4"	34 - 46 cm 13 1/4" - 18"	37 - 49 cm 14 1/2" - 19 1/4"
Ankle Circumference	18.5 - 22.5 cm 7 1/4" - 8 3/4"	19.5 - 23.5 cm 7 3/4" - 9 1/4"	20.5 - 24.5 cm 8" - 9 3/4"	21.5 - 25.5 cm 8 1/2" - 10"

PCSZ-01-65

POWER VIBE

Unisex, Stylish Compression Sock

Juzo **POWER Vibe Socks** designed with colorful patterns for everyday fashion.



Power Vibe 20-30 mmHg available Fall 2021

Size Chart for 2650 Juzo Power Vibe 15-20 mmHg				
	S	M	L	XL
Men’s Shoe Size	6	6.5 - 8	8.5 - 10.5	11 - 14
Women’s Shoe Size	5 - 7	7.5 - 9.5	10 - 12	12 +
Calf Circumference	30 - 41 cm 11 3/4" - 16"	32 - 44 cm 12 1/2" - 17 1/4"	34 - 46 cm 13 1/4" - 18"	37 - 49 cm 14 1/2" - 19 1/4"
Ankle Circumference	18.5 - 22.5 cm 7 1/4" - 8 3/4"	19.5 - 23.5 cm 7 3/4" - 9 1/4"	20.5 - 24.5 cm 8" - 9 3/4"	21.5 - 25.5 cm 8 1/2" - 10"

PCSZ-01-66