

SIGVARIS
GROUP

2021 Compression Solutions



Graduated Compression
Socks & Hosiery

Previous	Series Number	Series	Segment	
COTTON	230	Cotton	Essential	
SELECT COMFORT	860	Opaque		
SEA ISLAND COTTON	220	Sea Island Cotton	Style	
ALL-SEASON MERINO WOOL	240	Merino Wool		
ALLURE	710	Patterns		
MIDSHEER	750	Medium Sheer		
EVERSHEER	780	Sheer		
MIDTOWN MICROFIBER	820	Microfiber		
MICROFIBER SHADES	830	Microfiber Patterns		
SOFT OPAQUE	840	Soft Opaque		
CUSHIONED COTTON	360	Cushioned Cotton		Motion
TRAVERSE SOCKS	412	High Tech		
MERINO WOOL	422	Thermoregulating Wool		
SWISS MOUNTAIN SOCK	462	Swiss Mountain Sock		
CANADA MOUNTAIN SOCK	472	Canada Mountain Sock		
DAILY COMFORT	850	Comfort		

The logo for SIGVARIS GROUP, featuring the company name in white and blue text on a dark blue square background.

Established in 1864, SIGVARIS GROUP is, and has always been, a family-owned company. Founded and headquartered in Switzerland, our passion for helping people feel their best is why we've been a global market leader for 155 years and counting.

Our compression garment portfolio includes Medical, Consumer Health Care and products for Medically Complex Extremities that are specifically designed to restore, assist and enhance people's lives. We are renowned for making exceptional products that combine advanced technology and quality craftsmanship with diverse styles that fit our modern customers' lives.

When heritage, competence, and dedication come together, there's no limit to where we'll go.

Live your life. Be yourself.

For more information, please visit:
sigvaris.com

TO ORDER:

800.363.4999 (phone)
800.263.8736 (fax)
canadaorders@sigvaris.com

Online Dealer Portal:

sigvaris-online.com

	SERIES	PAGE	COMPRESSION LEVELS	STYLES
Essential	Cotton (230)	8	20–30mmHg 30–40mmHg	closed toe: calf, calf w/grip-top, thigh-high w/grip-top open toe: calf, thigh-high w/grip-top
	Opaque (860)	10	20–30mmHg 30–40mmHg	closed toe: calf, petite calf, calf plus, calf w/grip-top, thigh w/grip-top, petite thigh-high w/grip-top, extra-large thigh w/grip-top, pantyhose, petite pantyhose, plus-sized pantyhose open toe: calf, calf plus, calf w/grip-top, thigh w/grip-top, extra-large thigh w/grip-top, thigh-high w/waist attachment
Style	Sheer (780) <i>(for women only)</i>	14	15–20mmHg 20–30mmHg 30–40mmHg	closed toe: calf, thigh-high w/grip-top, pantyhose, maternity pantyhose open toe: calf, thigh-high w/grip-top, pantyhose
	Medium Sheer (750) <i>(for women only)</i>	16	20–30mmHg	closed toe: calf, thigh-high w/grip-top, pantyhose
	Soft Opaque (840) <i>(for women only)</i>	18	15–20mmHg 20–30mmHg 30–40mmHg	closed toe: calf, calf w/grip-top, thigh-high w/grip-top, pantyhose, maternity pantyhose open toe: calf, thigh-high w/grip-top, pantyhose
	Patterns (710) <i>(for women only)</i>	20	15–20mmHg 20–30mmHg	closed toe: thigh-high w/grip-top, pantyhose
	Sea Island Cotton (220)	ZÜRICH COLLECTION 22	20–30mmHg	closed toe: calf
	Merino Wool (240)	ZÜRICH COLLECTION 24	20–30mmHg	closed toe: calf
	Microfiber (820) <i>(for men only)</i>	26	15–20mmHg 20–30mmHg 30–40mmHg	closed toe: calf, calf w/grip-top, thigh w/grip-top, extended foot sizes
	Microfiber Patterns (830)	28	20–30mmHg	closed toe: calf
	Motion	Cushioned Cotton (360)	32	20–30mmHg
High Tech (412)		34	20–30mmHg	closed toe: calf
Thermoregulating Wool (422)		36	20–30mmHg	closed toe: calf
Swiss Mountain Sock (462) Canada Mountain Sock (472)		38	20–30mmHg	closed toe: calf
Comfort (850)		40	20–30mmHg	closed toe: calf
Specialty	Natural Rubber (500)	44	30–40mmHg 40–50mmHg 50–60mmHg	open toe: calf, thigh-high, thigh w/grip-top, thigh-high w/waist attachment, pantyhose
	Advance Armsleeve (910)	46	15–20mmHg 20–30mmHg 30–40mmHg	armsleeve w/ grip-top w/ gauntlet, armsleeve w/ grip-top w/out gauntlet

SIGVARIS CONSUMER HEALTH

SERIES	PAGE	COMPRESSION LEVEL	STYLES
CASUAL COTTON (146/186)	50	15–20mmHg	closed toe: calf
TRAVENO TRAVEL SOCKS (146/T / 186/T)	52	15–20mmHg	closed toe: calf
CUSHIONED COTTON (142/182) <small>ZURICH COLLECTION</small>	54	15–20mmHg	closed toe: calf
SEA ISLAND COTTON (151/191) <small>ZURICH COLLECTION</small>	56	15–20mmHg	closed toe: calf
ALL-SEASON MERINO WOOL (152/192)	58	15–20mmHg	closed toe: calf
MICROFIBER SHADES (143/183)	60	15–20mmHg	closed toe: calf
SHEER FASHION (120) <i>(for women only)</i>	62	15–20mmHg	closed toe: calf, thigh-high w/grip-top, pantyhose, maternity pantyhose open toe: calf
SOFT SILHOUETTE LEGGINGS (170) <i>(for women only)</i>	64	15–20mmHg	open foot: leggings
ATHLETIC RECOVERY SOCKS (401)	66	15–20mmHg	closed toe: calf
PERFORMANCE SLEEVES (412)	68	20–30mmHg	open foot: calf sleeve
MERINO OUTDOOR SOCKS (421)	70	15–20mmHg	closed toe: calf
EVERSOFT DIABETIC SOCKS (160)	72	8–15mmHg	closed toe: calf

Medically Complex Edema & Inelastic Wraps	76
Accessories	86
CoverOne™ Face Masks	90

sigvaris

We are all Natural

SIGVARIS GROUP products have always been made for everyone.

With that in mind, SIGVARIS GROUP is changing some colour names to better reflect the natural range of skin tones in the world.

Same Colours, New Names

Previous Colour	#	New Colour
Natural	33	Honey
Suntan	36	Golden
Nude	35	Chai
Crispa	66	Light Beige

Find your Natural

See individual product pages for more information.



Be yourself.

ESSENTIAL

Support Yourself

Sigvaris Essential solutions always meet your customers' needs. These straightforward stockings are an everyday go-to for people who seek the comfort, support, and reliability of traditional therapeutic socks and stockings.

Why choose **Essential**?

- Reduce swelling and discomfort
- Long-distance travel
- Everyday comfort and reliability

STYLE

Express Yourself

Choice is empowering, and fashion can be as unique as you are. With Sigvaris Style, personality will shine as your customers select from distinctive hosiery and socks.

Why choose **Style**?

- Fashionable, modern patterns and designs
- Special occasion and formal wear
- Our most luxurious feel

MOTION

Push Yourself

When daily comfort for hard-working legs is required, put your customers in motion. Sigvaris Motion socks offer all day strong support for people with active lives.

Why choose **Motion**?

- Cushioned sole for added comfort
- Durable thermoregulating wool
- Moisture-wicking yarns

sigvaris



Essential

ESSENTIAL

Support Yourself

Why choose **Essential**?

- Reduce swelling and discomfort
- Long-distance travel
- Everyday comfort and reliability

Essential Products

Previous	Series Number	New	Page
COTTON	230	Cotton	8
SELECT COMFORT	860	Opaque	10

230

Essential COTTON

for women & men

Essential Cotton combines all-day cotton comfort and a classic ribbed design that is ideal for sensitive skin.

Compression Level

20-30mmHg

30-40mmHg

Features and Benefits



**Sensitive
Skin**



Cotton



Breathable

- Supima® cotton for highest level of comfort and breathability
- Ideal for sensitive skin and allergies
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Soft to the skin with a silicone dotted grip-top band to prevent slipping

Fibre Content


- 50% Nylon, 25% Spandex, 25% Supima Cotton (content can vary slightly by style)
- Latex free

Best Suited for

- Containment of various types of edema, including venous edema, phlebolympheidema, and early stages of lymphedema
- All-day cotton comfort in all types of climates
- Ideal for dress or casual shoes
- Ideal for long-distance travel



SUPIMA®
WORLD'S FINEST COTTONS



Silicone
dotted grip-top
band in calf
and thigh

20-30mmHg
Closed Toe Calf w/ Grip-top
Medium Long in Black

232CMLM99/S



20-30mmHg
Closed Toe Calf Medium Short in Black Mist

232CMSW14

Sizing Chart

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")		29–36cm (11.5"–14")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")		41–61cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Thigh Circumference	45.5–66cm (18"–26")		48–73cm (19"–28.5")		52–84cm (20.5"–33")		NA	
Leg Length: Women & Open Toe	Short	Long	Short	Long	Short	Long		
	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")		
Leg Length: Men	Short	Long	Short	Long	Short	Long		
	≤76cm (30")	>76cm (30")	≤76cm (30")	>76cm (30")	≤76cm (30")	>76cm (30")		
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

Styles and Colours

Closed Toe = ● Open Toe = ○					
	WHITE 00	NAVY 10	BLACK MIST 14	LIGHT BEIGE 66	BLACK 99
Calf (C)	○	●	●	●○	●
Calf w/Grip-Top (C/S) Silicone dotted grip-top band				●*	●**
Thigh w/Grip-Top (N) Silicone dotted grip-top band				●○	●

* 20–30mmHg only

** 30–40mmHg calf with grip-top only available in black, 30–40mmHg extra large calf only available in black with grip-top

Product Code Example: **232NMLW66**

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
23	2	N	ML	W	66

Product Code Example: **232CLSM66/S**

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour	Grip-top
23	2	C	LS	M	66	/S

Essential

OPAQUE

for women & men

Essential Opaque offers durability you can trust in our largest range of sizes, making it your go-to choice for everyday hosiery.

Compression Level

20-30mmHg

30-40mmHg

Features and Benefits



Durable



**Largest
Size Range**



Opaque


- Largest size range
- Durable construction combined with a smooth, opaque look
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Soft to the skin with a silicone dotted grip-top band to prevent slipping
- Cotton-lined gusset for pantyhose styles
- Plus size pantyhose great alternative for maternity

Fibre Content

- 64% Nylon, 36% Spandex
(content can vary slightly by style)
- Latex free

Best Suited for

- Containment of various types of edema including venous edema, phlebolympheidema, and early stages of lymphedema
- Post-vein treatment including thermal and non-thermal techniques
- People looking for petite and plus-sized options



Silicone
dotted grip-top
band in thigh
and calf
w/ Grip-top

20-30mmHg
Closed Toe Thigh Medium Long in Golden

862NMLW36



20-30mmHg
Closed Toe Calf Medium Short in Black

862CMSM99

Sizing Chart

	Small			Medium			Large		X Large	
Ankle Circumference	18–21.5cm (7"–8.5")			21.5–25.5cm (8.5"–10")			25.5–29cm (10"–11.5")		29–36cm (11.5"–14")	
Calf Circumference	28–39cm (11"–15.5")			31.5–45cm (12.5"–17.5")			35.5–51cm (14"–20")		41–61cm (16"–24")	
Calf Circumference: Calf Plus w/Grip-Top	NA			39.5–53cm (15.5"–21")			45.5–61cm (18"–24")		51–71cm (20"–28")	
Calf Length	Petite	Short	Long	Petite	Short	Long	Short	Long	Short	Long
	≤35.5cm (14")	≤41cm (16")	>41cm (16")	≤35.5cm (14")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Thigh Circumference	45.5–66cm (18"–26")			48–73cm (19"–28.5")			52–84cm (20.5"–33")		64–95cm* (25"–37.5")	
Leg Length: Women & Open Toe	Petite	Short	Long	Petite	Short	Long	Short	Long	Short	Long
	≤65cm (25.5")	≤72cm (28.5")	>72cm (28.5")	≤65cm (25.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")
Hip Circumference Guidelines	83–112.5cm (32.5"–44.5")			87–127cm (34.5"–50")			98–150.5cm (38.5"–59.5")		NA	
Size	SP	SS	SL	MP	MS	ML	LS	LL	XS	XL
Plus (Maternity) Size Hip Circumference Guidelines	≤142cm (56")			≤157.5cm (62")			≤178cm (70")		NA	
Calf Plus w/Grip-Top					MS/P	ML/P	LS/P	LL/P	XS/P	XL/P

Women's Shoe Sizes: Petite Leg 4–8, Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

Styles and Colours

Closed Toe = ●
Open Toe = ○

HONEY 33	GOLDEN 36	LIGHT BEIGE 66	BLACK 99
-------------	--------------	----------------------	-------------

Women's Calf (C)	●	●	○	●●
Women's Calf Petite Sizes (C) Only available in 20–30mmHg, not available with dotted grip-top band	●	●		●
Women's Calf w/Grip-Top (C/S) Silicone dotted grip-top band			●●	●
Women's Calf Plus w/Grip-Top (C/P) Only available with silicone dotted grip-top band			●●	
Women's Thigh w/Grip-Top (N) Silicone dotted grip-top band**	●	●	○	●
Women's Thigh Petite Sizes (N) Only available in 20–30mmHg	●	●		●
Women's Pantyhose (P)	●	●		●
Women's Pantyhose Plus (P/P)***			●	●
Men's Calf (C) Only available in 20–30mmHg			●	●
Unisex Open Toe Thigh w/Grip-Top (N) Silicone dotted grip-top band**			○	
Unisex Thigh w/Waist Attachment (W) Indicate leg LEFT (/L) or RIGHT (/R)*			○	

* Not available in XL with waist attachment

** Women's Thigh w/Grip-Top (N) and Unisex Open Toe Thigh w/Grip-Top (N) sizes XS and XL are only available in Light Beige

*** Pantyhose plus is an alternative for maternity

Product Code Example: 863PMSW33

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
86	3	P	MS	W	33

Product Code Example: 862CLLW66/P

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour	Grip-top with Plus size
86	2	C	LL	W	66	/P

sigvaris



STYLE

Express Yourself

Why choose **Style**?

- Fashionable, modern patterns and designs
- Special occasion and formal wear
- Our most luxurious feel

Style Products

Previous	Series Number	New	Page
EVERSHEER	780	Sheer	14
MIDSHEER	750	Medium Sheer	16
SOFT OPAQUE	840	Soft Opaque	18
ALLURE	710	Patterns	20
SEA ISLAND COTTON	220	Sea Island Cotton	22
ALL-SEASON MERINO WOOL	240	Merino Wool	24
MIDTOWN MICROFIBER	820	Microfiber	26
MICROFIBER SHADES	830	Microfiber Patterns	28

780

Style

SHEER

for women

Style Sheer hosiery offers fashionable sheerness you can feel comfortable in, making it ideal for dress and professional wear.

Compression Level

15-20mmHg

20-30mmHg

30-40mmHg

Features and Benefits



Sheer



Elegant



Breathable


- Fashionable for professional and dress wear
- Consistent sheerness throughout the leg
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Comfortable Sensinnov® grip-top band prevents slipping
- Cotton-lined gusset for pantyhose styles

Fibre Content

- 15-20, 20-30mmHg: 69% Nylon, 31% Spandex (content can vary slightly by style)
- 30-40mmHg: 60% Nylon, 40% Spandex (content can vary slightly by style)
- Latex free

Best Suited for

- Management of Chronic Venous Disease (CVD)
- Pre- and post-vein treatment with minimal edema
- Pregnancy
- People who prefer fashionable options
- Ideal for business or dress



Sensinnov
grip-top band
in thigh

15-20mmHg
Closed Toe Thigh Medium Long in Dark Navy

781NMLW08



20-30mmHg
Closed Toe Calf Medium Long in Mocha

782CMLW85

Sizing Chart

	Small		Medium		Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Thigh Circumference	45.5–66cm (18"–26")		48–73cm (19"–28.5")		52–84cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")
Hip Circumference Guidelines	83–112.5cm (32.5"–44.5")		87–127cm (34.5"–50")		98–150.5cm (38.5"–59.5")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes Short 5.5–9.5, Long 8–12

Maternity Hip Circumference Guidelines Small ≤142cm (56"), Medium ≤157.5cm (62"), Large ≤178cm (70")

Styles and Colours

	DARK NAVY 08	HONEY 33	WARM SAND 34*	GOLDEN 36	TOASTED ALMOND 37*	CAFÉ 73*	MOCHA 85*	NIGHT-SHADE 94	BLACK 99
Calf (C)	●	●●	●●	●●	●●	●●	●●	●	●●
Thigh w/Grip-Top (N) <small>Sensinnov grip-top band</small>	●	●●	●●	●●	●●	●●	●●	●	●●
Pantyhose (P)	●	●●	●●	●●	●●	●●	●●	●	●●
Maternity Pantyhose (M)*		●							●

Closed Toe = ●
Open Toe = ○

* 20–30mmHg ONLY

Product Code Example: **781NMLW36**

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
78	1	N	ML	W	36

Product Code Example: **783PSSO33**

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
78	3	P	SS	O	33

750

Style

MEDIUM SHEER

for women

Thanks to the durability of Style Medium Sheer, you will have the confidence to wear sheer hosiery every day.

Compression Level

20-30mmHg

Features and Benefits



Medium Sheer



Durable



Breathable


- Fashionable, durable, medium sheer hosiery for everyday wear
- Excellent breathability and moisture management
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Soft to the skin with a silicone dotted grip-top band to prevent slipping
- Comfortable and large 2" top band on calf style
- Cotton-lined gusset for pantyhose styles

Fibre Content

- 68% Nylon, 32% Spandex
(Content can vary slightly by style)
- Latex free

Best Suited for

- Management of Chronic Venous Disease (CVD)
- Pre- and post-vein treatment with minimal edema
- Daily wear for business or casual
- Active professionals needing a durable sheer



Silicone
dotted grip-top
band in thigh

20-30mmHg
Closed Toe Thigh Medium Long in Black

752NMLW99



20-30mmHg
Closed Toe Calf Small Long in Nightshade

752CSLW94


Sizing chart

	Small		Medium		Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Thigh Circumference	45.5–66cm (18"–26")		48–73cm (19"–28.5")		52–84cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")
Hip Circumference Guidelines	83–112.5cm (32.5"–44.5")		87–127cm (34.5"–50")		98–150.5cm (38.5"–59.5")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes Short leg 5.5–9.5, Long leg 8–12.

Styles and Colours

Closed Toe = ●



	DARK NAVY 08	HONEY 33	NIGHT-SHADE 94	BLACK 99
Calf (C)	●	●	●	●
Thigh w/Grip-Top (N) <small>Silicone dotted grip-top band</small>	●	●	●	●
Pantyhose (P)	●	●	●	●

Product Code Example: **752NMLW33**

Series	mmHg Level	Style	Size	Gender	Colour
75	2	N	ML	W	33

Style

SOFT OPAQUE

for women

Style Soft Opaque hosiery blends a fashionable opaque look with long-lasting microfiber softness for your professional and everyday wear.

Compression Level

15-20mmHg


20-30mmHg

30-40mmHg

Features and Benefits



Opaque



Microfiber



Soft


- Fashionable opaque hosiery for professional and everyday wear
- Lasting microfiber softness and excellent moisture management
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Comfortable Sensinnov grip-top band on thigh style prevents slipping
- Cotton-lined gusset for pantyhose styles

Fibre Content

- 79% Nylon, 21% Spandex
(Content can vary slightly by style)
- Latex free

Best Suited for

- Management of Chronic Venous Disease (CVD)
- Pre- and post-vein procedure with minimal edema
- People who prefer fashionable options
- Pregnancy
- Long distance travel



Sensinnov
grip-top band
in thigh

20-30mmHg
Closed Toe Thigh Small Long in Midnight Blue

842NSLW09



Silicone
dotted grip-top
band in calf

20-30mmHg
Closed Toe Calf Medium Long in Black

842CMLW99/S

Sizing Chart

	Small		Medium		Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Thigh Circumference	45.5–66cm (18"–26")		48–73cm (19"–28.5")		52–84cm (20.5"–33")	
Hip Circumference Guidelines	83–112.5cm (32.5"–44.5")		87–127cm (34.5"–50")		98–150.5cm (38.5"–59.5")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes Short leg 5.5–9.5, Long leg 8–12 Maternity hip circumference guidelines Small: <142cm (<56"), Medium: <157.5cm (<62"), Large: <178cm (<70")

Styles and Colours

	MIDNIGHT BLUE 09	CHAI 35	PECAN 40	MULBERRY 60*	BLUE IRIS 68*	CARAMEL 83*	ESPRESSO 89*	GRAPHITE 91	BLACK 99
Calf (C)	●	●○	●	●	●	●	●	●	●○
Calf w/Grip-Top (C/S)*									●
Thigh w/Grip-Top (N) <small>Sensinnov grip-top band</small>	●	●○	●	●	●	●	●	●	●○
Pantyhose (P)	●	●○	●	●	●	●	●	●	●○
Maternity Pantyhose (M)*									●

* 20–30mmHg ONLY

Product Code Example: 842CMSW68

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
84	2	C	MS	W	68

Style PATTERNS

for women

Style Patterns hosiery is available in fashionable patterns, making it a staple for your professional and everyday wear.

Compression Level

15-20mmHg

20-30mmHg

Features and Benefits



Sheer



Elegant



Breathable


- Stylish diamond-dot and stripe patterns for professional and dress wear
- Excellent breathability and moisture management
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Comfortable Sensinnov top band on thigh style prevents slipping
- Cotton-lined gusset for pantyhose styles

Fibre Content


- 64% Nylon, 36% Spandex
(Content can vary slightly by style)
- Latex free

Best Suited for

- Management of Chronic Venous Disease (CVD)
- Pre- and post-vein treatment with minimal edema
- People who prefer fashionable options
- Ideal for business or dress



Scalloped Sensinnov® Thigh band on Stripe patterns



Sensinnov thigh band on Diamond Dot patterns

15-20mmHg
Closed Toe Thigh Medium Long in Graphite

711NMLW91



20-30mmHg
Closed Toe Pantyhose Small Long in Onyx Stripe

712PSLW19

Sizing chart

	Small		Medium		Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")	
Thigh Circumference	45.5–66cm (18"–26")		48–73cm (19"–28.5")		52–84cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")	≤72cm (28.5")	>72cm (28.5")
Hip Circumference Guidelines	83–112.5cm (32.5"–44.5")		87–127cm (34.5"–50")		98–150.5cm (38.5"–59.5")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes Short leg 5.5–9.5, Long leg 8–12

Styles and Colours

	GRAPHITE STRIPE 17*	ONYX STRIPE 19*	GRAPHITE 91	BLACK 99
Closed Toe = ●	●	●	●	●
Thigh w/Grip-Top (N) Sensinnov grip-top band	●	●	●	●
Pantyhose (P)	●	●	●	●

* Available in 20–30mmHg only

Product Code Example: **712NLSW17**

Series	mmHg Level	Style	Size	Gender	Colour
71	2	N	LS	W	17

220

Style

SEA ISLAND COTTON

for women & men

Style Sea Island Cotton socks are made with the world's rarest cotton, ensuring a soft and luxurious feel.

Compression Level

20-30mmHg

Also available in 15-20mmHg, see page 56

Features and Benefits



Soft



Cotton



Breathable

- Knitted with the world's rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

Fibre Content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Latex free

Best Suited for

- Management of Chronic Venous Disease (CVD)
- People who prefer natural fibres for everyday wear
- Long-distance travel



20-30mmHg
Closed Toe Calf Small Long in Brown

222CSLW11



20-30mmHg
Closed Toe Calf Medium Long in Black

222CMLM99


Sizing chart

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")		29–36cm (11.5"–14")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")		41–61cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

Styles and Colours

Closed Toe = ●



Calf (C)	●	●	●
----------	---	---	---

Product Code Example: **222CLLM10**

Series	mmHg Level	Style	Size	Gender	Colour
22	2	C	LL	M	10

240

Style

MERINO WOOL

for women & men

Style Merino Wool socks combine extra-fine Merino wool with a stylish design, making it the ultimate sock for year-round wear.

Compression Level

20-30mmHg

Also available in 15-20mmHg, see page 58

Features and Benefits



Wool



Thermoregulating



Stylish

- Made with the highest quality extra-fine Merino wool
- Thermoregulating for year-round wear
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

Fibre Content

- 53% Merino Wool, 38% Nylon, 9% Spandex
- Latex free

Best Suited for

- Management of Chronic Venous Disease (CVD)
- People who prefer fashionable and premium fabric options
- Long-distance travel



20-30mmHg
Closed Toe Calf Medium Long in Brown

242CMLM11



20-30mmHg
Closed Toe Calf Medium Long in Navy

242CMLW10


Sizing chart

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")		29–36cm (11.5"–14")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")		41–61cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

Styles and Colours

Closed Toe = ●



Calf (C)	●	●	●
-----------------	---	---	---

Product Code Example: **242CMSW99**

Series	mmHg Level	Style	Size	Gender	Colour
24	2	C	MS	W	99

Style

MICROFIBER

for men

Style Microfiber socks offer microfiber softness for your professional and everyday wear.

Compression Level

15-20mmHg

20-30mmHg

30-40mmHg

Features and Benefits



Microfiber



Soft



Stylish

- 100% of microfiber sits next to skin for a soft feel
- Excellent breathability and moisture management for added comfort
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Sensinno grip-top band on thigh style

Fibre Content

- 87% Nylon , 13% Spandex
(content can vary slightly by style)
- Latex free

Best Suited for


- Management of Chronic Venous Disease (CVD)
- People who prefer fashionable options
- Men with larger feet
- Long-distance travel



Silicone
dotted grip-top
band in calf

20-30mmHg
Closed Toe Calf Medium Long in Black

822CMLM99/S



Sensinno
grip-top band
in thigh

30-40mmHg
Closed Toe Thigh Medium Long in Tan Khaki

823NMLM32

Sizing Chart


	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")		29–36cm (11.5"–14")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")		41–61cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤ 41cm (16")	> 41cm (16")	≤ 41cm (16")	> 41cm (16")	≤ 41cm (16")	> 41cm (16")	≤ 41cm (16")	> 41cm (16")
Thigh Circumference	45.5–66cm (18"–26")		48–73cm (19"–28.5")		52–84cm (20.5"–33")		NA	
Leg Length	Short	Long	Short	Long	Short	Long		
	≤ 76cm (30")	> 76cm (30")	≤ 76cm (30")	> 76cm (30")	≤ 76cm (30")	> 76cm (30")		
Size Shoe size	SS	SL	MS	ML	LS	LL	XS	XL
	7–11	9.5–14	7–11	9.5–14	7–11	9.5–14	7–11	9.5–14
				MX* 12.5–17.5		LX* 12.5–17.5		XX* 12.5–17.5

Men's shoe sizes: Short leg 7–11, long leg 9.5–14 (For extended foot sizes: long leg 12.5–17.5)

* MX, LX, and XX extended foot sizes only available in calf style in Black & Tan Khaki in 15–20 & 20–30mmHg only. Extended foot sizes not available with grip top

Styles and Colours

Closed Toe = ●



	NAVY 10***	TAN KHAKI 32	STEEL GREY 93***	BLACK 99
Calf (C)	●	●	●	●
Calf w/Grip-Top (C/S)** Silicone dotted grip-top band		●		●
Thigh w/Grip-Top (N) Sensinnov grip-top band		●		●

** Calf w/ Grip-Top not available in 15–20mmHg

*** Navy and Steel Grey colours only available in 15–20mmHg and 20–30mmHg

Product Code Example: **822NLSM99**

Series	mmHg Level	Style	Size	Gender	Colour
82	2	N	LS	M	99

Product Code Example: **822CMLM32/S**

Series	mmHg Level	Style	Size	Gender	Colour	Grip-top
82	2	C	ML	M	32	/S

Style MICROFIBER PATTERNS

for women & men

Style Microfiber Patterns socks combine ultimate comfort with a variety of fashionable colours and designs to spice up your wardrobe.

Compression Level

20-30mmHg

Also available in 15-20mmHg, see page 60

Features and Benefits



Microfiber



Stylish



Patterns

- 100% of microfiber sits next to skin for a soft feel
- Fashionable, durable socks for everyday wear
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

Fibre Content

- 87% Nylon, 13% Spandex
(content can vary slightly by style)
- Latex free


Best Suited for

- Management of Chronic Venous Disease (CVD)
- People who prefer fashionable options
- Long-distance travel



20-30mmHg
Closed Toe Calf Small Long in Royal Blue Argyle

832CSLW48



20-30mmHg
Closed Toe Calf Large Short in Pink Stripe

832CLSW15

Sizing Chart

	Small		Medium		Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Size	SS	SL	MS	ML	LS	LL

Men's Shoe sizes: Short leg 7–11, Long leg 9.5–14 Women's Shoe Sizes: Short leg 5.5–9.5, Long leg 8–12

Styles and Colours

Closed Toe = ●	GRAPHITE MINI-STRIPE 07	DARK NAVY STRIPE 13	PINK STRIPE 15*	GRAPHITE STRIPE 17	ONYX STRIPE 19	GRAPHITE HEATHER 27
	●	●	●	●	●	●
Calf (C)	PINK ARGYLE 45	PURPLE ARGYLE 46	GRAPHITE ARGYLE 47	ROYAL BLUE ARGYLE 48	ONYX ARGYLE 49	
	●	●	●	●	●	

* Available for women only

Product Code Example: **832CMLM17**

Series	mmHg Level	Style	Size	Gender	Colour
83	2	C	ML	M	17

sigvaris



MOTION

Push Yourself

Why choose **Motion**?

- Cushioned sole for added comfort
- Durable thermoregulating wool (422)
- Moisture-wicking yarns

Motion Products

Previous	Series Number	New	Page
CUSHIONED COTTON	360	Cushioned Cotton	32
TRAVERSE SOCKS	412	High Tech	34
MERINO WOOL	422	Thermoregulating Wool	36
SWISS MOUNTAIN SOCK	462	Swiss Mountain Sock	38
CANADA MOUNTAIN SOCK	472	Canada Mountain Sock	38
DAILY COMFORT	850	Comfort	40

360

Motion CUSHIONED COTTON

for women & men

Motion Cushioned Cotton is soft and breathable with an extra-padded sole for optimal comfort, and is ideal for sensitive skin.

Compression Level

20-30mmHg

Also available in 15-20mmHg, see page 54

Features and Benefits



**Cushioned
Sole**



Cotton



**Sensitive
Skin**

- Supima cotton for highest level of comfort and breathability
- Extra-padded sole, heel, and toe to protect feet
- Flat, low-friction toe seam
- Diabetic friendly

Fibre Content

- 64% Supima Cotton, 28% Nylon, 8% Spandex
- Latex free

Best Suited for

- Management of Chronic Venous Disease (CVD)
- People who stand for prolonged periods of time
- People with sensitive skin and toes
- Work shoes and boots
- Casual shoes and sports shoes



SUPIMA®
WORLD'S FINEST COTTONS



20-30mmHg
Closed Toe Calf Medium Long in Black

362CMLM99



20-30mmHg
Closed Toe Calf Small Long in White

362CSLW00

Sizing Chart

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")		29–36cm (11.5"–14")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")		41–61cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes Short 5.5–9.5 Long 8–12 Men's Shoe Sizes Short 7–11 Long 9.5–14

Styles and Colours

Closed Toe = ●

	 WHITE 00	 BLACK 99
Calf (C)	●	●

Product Code Example: **362CMLM00**

Series	mmHg Level	Style	Size	Gender	Colour
36	2	C	ML	M	00

Motion HIGH TECH

for women & men

Motion High Tech delivers optimal comfort where it counts with an anatomically contoured shape and a cushioned sole.

Compression Level

20-30mmHg

Features and Benefits



Moisture Wicking



Anatomically Contoured



Colours

- High-tech fibres for improved thermoregulation
- Cushioned footbed, toe, and Achilles tendon zones for optimal comfort
- Additional compression in calf improves circulation

Fibre Content

- 64% Nylon, 17% Polyester, 19% Spandex
- Latex free


Best Suited for

- Management of Chronic Venous Disease (CVD)
- People who stand for prolonged periods of time
- Walking and everyday activities
- Casual and occupational footwear
- Active individuals



20-30mmHg
Closed Toe Calf Small Medium in Limeade

412CSM64



20-30mmHg
Closed Toe Calf Medium Large in Pink

412CML56

Sizing Chart

	Small			Medium			Large			X Large		
Ankle Circumference	18–21.5cm (7"–8.5")			21.5–25.5cm (8.5"–10")			25.5–29cm (10"–11.5")			29–33cm (11.5"–13")		
Calf Circumference	28–39cm (11"–15.5")			31.5–45cm (12.5"–17.5")			35.5–51cm (14"–20")			40.5–58cm (16"–23")		
Women's Shoe Size	4.5-7.5	8-10.5	11-13.5	4.5-7.5	8-10.5	11-13.5	8-10.5	11-13.5	>14	8-10.5	11-13.5	>14
Men's Shoe Size	≤6.5	7-9.5	10-12.5	≤6.5	7-9.5	10-12.5	7-9.5	10-12.5	13-15.5	7-9.5	10-12.5	13-15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL	XS	XM	XL

Styles and Colours

	WHITE 00	BLUE 50	LIME 54	PINK 56	RED 59	STEEL BLUE 62	BLOOD ORANGE 63	LIMEADE 64	BLACK 99
Closed Toe = ●									
Calf (C)	●	●	●	●	●	●	●	●	●

Product Code Example: **412CML62**

Series	mmHg Level	Style	Size	Colour
41	2	C	ML	62

Motion THERMOREGULATING WOOL

for women & men

Motion Thermoregulating Wool combines the benefits of natural Merino wool with cushioned soles, making it ideal for wherever the road takes you!

Compression Level

20-30mmHg

Also available in 15-20mmHg, see page 70

Features and Benefits



Wool



Cushioned Thermoregulating
Sole



- Natural properties of wool make these socks odour resistant
- Cushioned, padded soles help prevent blisters
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

Fibre Content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Latex free


Best Suited for

- Management of Chronic Venous Disease (CVD)
- People who stand for prolonged periods of time
- Work boots, steel-toed shoes, and most occupational footwear
- Outdoor work environments



20-30mmHg
Closed Toe Calf Large Long in Charcoal

422CLL12



20-30mmHg
Closed Toe Calf Small Long in Olive


422CSL31

Sizing Chart

	Small			Medium			Large		
Ankle Circumference	18-21.5cm (7"-8.5")			21.5-25.5cm (8.5"-10")			25.5-29cm (10"-11.5")		
Calf Circumference	28-39cm (11"-15.5")			31.5-45cm (12.5"-17.5")			35.5-51cm (14"-20")		
Women's Shoe Size	4.5-7.5	8-10.5	11-13.5	4.5-7.5	8-10.5	11-13.5	8-10.5	11-13.5	>14
Men's Shoe Size	≤6.5	7-9.5	10-12.5	≤6.5	7-9.5	10-12.5	7-9.5	10-12.5	13-15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL

Styles and Colours

Closed Toe = ●



Calf (C)	●	●
----------	---	---

Product Code Example: **422CML12**

Series	mmHg Level	Style	Size	Colour
42	2	C	ML	12

462/472

Motion

SWISS MOUNTAIN SOCK CANADA MOUNTAIN SOCK

for women & men

These socks are moisture-wicking and anatomically contoured, with a padded sole and Achilles tendon zone for ultimate comfort on the go.

Compression Level

20-30mmHg

Features and Benefits



**Cushioned
Sole**



**Anatomically
Contoured**



**Moisture
Wicking**

- Thermal and odour control fibres to promote moisture management
- Footbed, toe and Achilles tendon cushion help reduce blisters
- Swiss precision graduated compression

Fibre Content

- 52% Nylon, 25% Polypropylene, 12% Polyester, 11% Spandex
- Latex free

Best Suited for

- People who stand for prolonged periods of time
- Work boots, steel-toed shoes, and most occupational footwear
- Walking and outdoor activities



20-30mmHg Swiss Mountain Sock
Closed Toe Calf Large Long in Red

462CLL59



20-30mmHg Canada Mountain Sock
Closed Toe Calf Medium Short in Red

472CMS59

Sizing Chart

	Small			Medium			Large			X Large		
Ankle Circumference	18–21.5cm (7"–8.5")			21.5–25.5cm (8.5"–10")			25.5–29cm (10"–11.5")			29–33cm (11.5"–13")		
Calf Circumference	28–39cm (11"–15.5")			31.5–45cm (12.5"–17.5")			35.5–51cm (14"–20")			40.5–58cm (16"–23")		
Women's Shoe Size	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14	8–10.5	11–13.5	>14
Men's Shoe Size	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5	7–9.5	10–12.5	13–15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL	XS	XM	XL

Styles and Colour

Closed Toe = ●



Calf (C) ●

Swiss Mountain Sock
 Product Code Example: **462CML59**

Series	mmHg Level	Style	Size	Colour
46	2	C	ML	59

Canada Mountain Sock
 Product Code Example: **472CSM59**

Series	mmHg Level	Style	Size	Colour
47	2	C	SM	59

850

Motion COMFORT

for women & men

Motion Comfort offers cushioned soles and moisture wicking to keep your legs feeling energized all day long.

Compression Level

20-30mmHg

Features and Benefits



**Cushioned
Sole**



Durable



**Moisture
Wicking**

- NILIT® Aquarius fibres provide moisture wicking to keep feet dry and comfortable
- Cushioned soles to absorb impact and help prevent blisters
- Smooth, flat toe seam

Fibre Content

- 86% Nylon, 14% Spandex
- Latex free

Best Suited for

- Management of Chronic Venous Disease (CVD)
- Healthcare professionals
- Uniforms requiring black or solid colour
- Business travelers
- Individuals working in physically-demanding occupations



20-30mmHg
Closed Toe Calf Large Long in Black

852CLLM99



20-30mmHg
Closed Toe Calf Medium Long in Black

852CMLW99


Sizing Chart

	Small		Medium		Large	
Ankle Circumference	18–21.5cm (7"–8.5")		21.5–25.5cm (8.5"–10")		25.5–29cm (10"–11.5")	
Calf Circumference	28–39cm (11"–15.5")		31.5–45cm (12.5"–17.5")		35.5–51cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")	≤41cm (16")	>41cm (16")
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short leg 5.5–9.5, Long leg 8–12 Men's Shoe sizes: Short leg 7–11, Long leg 9.5–14

Styles and Colours

Closed Toe = ●



Calf (C) ●



Product Code Example: **852CMLM99**

Series	mmHg Level	Style	Size	Gender	Colour
85	2	C	ML	M	99

sigvaris



SPECIALTY

Take care of Yourself

Why choose **Specialty**?

- Specific medical needs
- Need low fatigue fibers
- Natural rubber (500)

Specialty Products

Product	Series Number	Page
Natural Rubber	500	44
Advance Armsleeve	910	46

500

Specialty NATURAL RUBBER

for women & men

Natural Rubber offers patients a durable, long-term containment solution for compression therapy with minimal fibre fatigue throughout the day.

Compression Level

30-40mmHg

40-50mmHg

50-60mmHg


Features and Benefits



Natural Rubber



Durable



Wide Size Range


- Ideal for patients with more serious indications
- Comfortable silicone dotted grip-top band on thigh styles to prevent slipping
- Extended sizing for a custom-like fit

Fibre Content

- 70% Nylon, 30% Natural Latex Rubber (Calf & Thigh-high styles)
- 75% Nylon, 20% Natural Latex Rubber, 5% other fibres (Pantyhose style)
- **Caution:** contains natural rubber latex that may cause allergic reaction

Best Suited for

- Severe varicosities
- Chronic edema post DVT
- Severe venous stasis
- Prevention of ulcer recurrence
- Early stages of lymphedema
- Active venous ulcers
- Ideally paired with the Cotton Underliner (page 87)



Silicone
dotted grip-top
band in thigh

30-40mmHg
Open Toe Thigh Small Average Long in Beige

503NS2077



40-50mmHg
Open Toe Calf Medium Full Long in Beige

504CM4077

Sizing Chart

	Small				Medium				Large				X Large			
Ankle Circumference	19-23cm (7.5"-9")				23-30cm (9"-12")				30-40cm (12"-15.5")				40-46cm (15.5"-18")			
Calf Circumference	Average		Full		Average		Full		Average		Full		Average		Full	
	28-38cm (11"-15")		≤43cm (17")		32-43cm (12.5"-17")		≤48cm (19")		36-48cm (14.5"-19")		≤57cm (22.5")		45-57cm (18"-22.5")		≤68cm (27")	
Calf Circumference for 505C (Full only)	505C: full calf only		32-43cm (12.5"-17")		505C: full calf only		36-48cm (14.5"-19")		505C: full calf only		45-57cm (18"-22.5")		505C: not available in XL			
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	≤40cm (≤15.5")	>40cm (>15.5")	≤40cm (≤15.5")	>40cm (>15.5")	≤40cm (≤15.5")	>40cm (>15.5")	≤40cm (≤15.5")	>40cm (>15.5")	≤40cm (≤15.5")	>40cm (>15.5")	≤40cm (≤15.5")	>40cm (>15.5")	≤40cm (≤15.5")	>40cm (>15.5")	≤40cm (≤15.5")	>40cm (>15.5")
Thigh Circumference	Average		Full		Average		Full		Average		Full		N/A			
	41-54cm (16"-21")		≤61cm (24")		48-61cm (19"-24")		≤71cm (28")		58-71cm (23"-28")		≤81cm (32")					
Thigh Circumference for 505T (Full only)	505T: full thigh only		≤61cm (24")		505T: full thigh only		≤71cm (28")		505T: full thigh only		≤81cm (32")					
Leg Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long				
	≤71cm (28")	>71cm (28")	≤71cm (28")	>71cm (28")	≤71cm (28")	>71cm (28")	≤71cm (28")	>71cm (28")	≤71cm (28")	>71cm (28")	≤71cm (28")	>71cm (28")				
Hip Circumference Guidelines	≤122cm (48")				≤137cm (54")				≤155cm (61")							
Size	S1	S2	S3	S4	M1	M2	M3	M4	L1	L2	L3	L4	X1	X2	X3	X4

Styles and Colours

Open Toe = O



Caution: this product contains natural rubber latex which may cause allergic reactions. For all indications, cautions, and contraindications, please see page 89 for details. Sigvaris accepts no liability for non-observance of contraindications or cautions.

Calf (C)	O	Calf only available in S, M, and L full circumference for 50-60mmHg
Thigh w/out Grip-Top (T)	O	
Thigh w/Grip-Top (N) Silicone dotted grip-top band	O	Not available in 50-60mmHg and non-available in XL
Panty hose (P)	O	Panty hose only available in 30-40mmHg and only in full thigh circumference
Thigh w/Waist Attachment (W) Indicate leg LEFT (/L) or RIGHT (/R)	O	Thigh w/ waist attachment is available in average and full thigh circumferences in 30-40mmHg and in full thigh circumferences only in 40-50mmHg

Product Code Example: **504CS2077**

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
50	4	C	S2	O	77

Product Code Example: **503WL1077/L**

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour	Leg
50	3	W	L1	O	77	/L

Specialty ADVANCE ARMSLEEVE

for women & men

The Advance Armsleeve combines ultimate softness with silver ion properties and a Sensinnov grip-top band that's ideal for all-day comfort and effectiveness.

Compression Level

15-20mmHg

20-30mmHg

30-40mmHg

Features and Benefits



**Moisture
Wicking**



**Sensinnov
Grip-top
Band**



Breathable


- Ideal for patients with arm lymphedema or post-mastectomy care
- Sensinnov grip-top band ensures armsleeve stays in place all day long
- Provides moisture wicking properties to keep skin dry and comfortable
- Breathable to provide comfortable, year-round wear

Fibre Content

- 50% Nylon, 28% Rayon, 22% Spandex (Content can vary slightly by style)
- Latex free

Best Suited for


- Early stages of lymphedema
- People with arm lymphedema
- Post-mastectomy care



**Sensinnov
grip-top
band**

15-20mmHg
Armsleeve Without Gauntlet
Small Long in Beige

911ASLO77/S



**With
Gauntlet**

20-30mmHg
Armsleeve With Gauntlet
Small Long in Beige


912ASLG77/S

Sizing Chart

	Small		Medium				Large				X Large			
Palm	17-21cm (6.5"-8")		18-22cm (7"-8.5")				19-23cm (7.5"-9")				20-24cm (8"-9.5")			
Wrist	14-16cm (5.5"-6.5")		16-18cm (6.5"-7")				18-20cm (7"-8")				20-23cm (8"-9")			
Upper Mid-Arm	Regular		Regular		Plus		Regular		Plus		Regular		Plus	
	26-34cm 10"-13.5"		29-37cm 11.5"-14.5"		37-43cm 14.5"-17"		32-40cm 12.5"-15.5"		40-46cm 15.5"-18"		34-42cm 13.5"-16.5"		42-48cm 16.5"-19"	
Arm Length	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long
	≤43cm (<17")	>43cm (>17")	≤43cm (<17")	>43cm (>17")	≤43cm (<17")	>43cm (>17")	≤43cm (<17")	>43cm (>17")	≤43cm (<17")	>43cm (>17")	≤43cm (<17")	>43cm (>17")	≤43cm (<17")	>43cm (>17")
Size	SR	SL	MR	ML	MR/P	ML/P	LR	LL	LR/P	LL/P	XR	XL	XR/P	XL/P

Styles and Colours

Sleeve = ●	BEIGE 77
Armsleeve w/Gauntlet (G) Sensinnov grip-top band	●
Armsleeve w/out Gauntlet (O) Sensinnov grip-top band	●



Note: all armsleeves come with Grip-top (/S)

15-20mmHg not available with gauntlet and not available in plus sizing

Armsleeve without gauntlet
Product Code Example: **912ASLO77/S**

Series	mmHg Level	Style	Size	W/out Gauntlet	Colour	Grip-Top
91	2	A	SL	O	77	/S

Armsleeve with gauntlet
Product Code Example: **912ALRG77/PS**

Series	mmHg Level	Style	Size	W/ Gauntlet	Colour	Plus Size	Grip-Top
91	2	A	LR	G	77	/P	S

NEW!



MARINER STRIPE

CONSUMER HEALTH

Care for Yourself

Why choose **Consumer Health**?

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease (CVD) Prevention of swelling and varicose veins
- Simplified sizing

Consumer Health Products

Product	Series Number	Page
CASUAL COTTON	146 / 186	50
TRAVENO TRAVEL SOCKS	146/T / 186/T	52
CUSHIONED COTTON	142 / 182	54
SEA ISLAND COTTON	151 / 191	56
ALL-SEASON MERINO WOOL	152 / 192	58
MICROFIBER SHADES	143 / 183	60
SHEER FASHION	120	62
SOFT SILHOUETTE LEGGINGS	170	64
ATHLETIC RECOVERY SOCKS	401	66
PERFORMANCE SLEEVES	412	68
MERINO OUTDOOR SOCKS	421	70
EVERSOFT DIABETIC SOCKS	160	72

146/186

CASUAL COTTON

for women & men

Compression Level

15-20mmHg

Features and Benefits

- » Made with premium Supima Cotton for optimum softness and durability
- » Flat, low-friction toe seam to prevent irritation
- » Soft, stretchable, knit-in band prevents pinching

Fibre Content

- » 66% Supima Cotton, 26% Nylon, 8% Spandex
- » Latex free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease (CVD)








SUPIMA[®]
WORLD'S FINEST COTTONS

Sizing Chart

146	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23cm (7"-9")	19-25cm (7.5"-10")	20-27cm (8"-10.5")
186	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	19-25cm (7.5"-10")	20-27cm (8"-10.5")	21-29cm (8.5"-11.5")
Size		A	B	C

Style and Colours

Closed Toe = ●

					
	WHITE 00	NAVY 10	BROWN 11	KHAKI 30	BLACK 99
Calf (C)	●	●	●	●	●

Product Code Example: **146CB10**

Series	Style	Size	Colour
146	C	B	10



146
Women's Calf
Style

15-20mmHg
Closed Toe Calf Size B in Khaki

146CB30



186
Men's Calf
Style

15-20mmHg
Closed Toe Calf Size C in Black

186CC99

146/T - 186/T

TRAVENO TRAVEL SOCKS

for women & men

Compression Level

15-20mmHg

Features and Benefits

- » Made with premium Supima Cotton for optimum softness and durability
- » Flat, low-friction toe seam to prevent irritation
- » Soft, stretchable knit-in band prevents pinching

Fibre Content

- » 66% Supima Cotton, 26% Nylon, 8% Spandex
- » Latex free

Best Suited for

- » Travelers
- » Tired, heavy, restless, aching legs
- » Occupational leg health risks, expecting mothers, family history of chronic venous disorders



SUPIMA®
WORLD'S FINEST COTTONS

Style and Colours

Closed Toe = ●

BLACK
99

Calf (C)



Sizing Chart

146/T	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23cm (7"-9")	19-25cm (7.5"-10")	20-27cm (8"-10.5")
186/T	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	19-25cm (7.5"-10")	20-27cm (8"-10.5")	21-29cm (8.5"-11.5")
Size		A	B	C

IMPROVE YOUR
CIRCULATION
IN THE AIR
OR ON
THE GROUND

Traveno 
TRAVEL SOCKS



Travel Socks for Women & Men
Available in 15–20mmHg in black

Countertop displays available. Contact your Territory
Manager or call **1.800.363.4999** for more information

CUSHIONED COTTON

for women & men

Compression Level

15–20mmHg

Also available in 20–30mmHg, see page 32

Features and Benefits

- » Cushioned sole, heel, and toe area to protect feet
- » Flat, low-friction toe seam
- » Soft, stretchable knit-in band prevents pinching
- » Ideal for everyday wear
- » Suitable for patients with diabetes
- » No compression in heel and instep for ease of donning

Fibre Content

- » 68% Supima Cotton, 25% Nylon, 7% Spandex
- » Latex free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease (CVD)
- » Sensitive skin




SUPIMA[®]
WORLD'S FINEST COTTONS

Sizing Chart

142	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23cm (7"-9")	19-25cm (7.5"-10")	20-27cm (8"-10.5")
182	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	19-25cm (7.5"-10")	20-27cm (8"-10.5")	21-29cm (8.5"-11.5")
Size		A	B	C

Style and Colours

Closed Toe = ●



Calf (C)



Product Code Example: 142CB00


Series	Style	Size	Colour
142	C	B	00



142
Women's
Calf Style

15-20mmHg
Closed Toe Calf Size B in White

142CB00



182
Men's Calf
Style

15-20mmHg
Closed Toe Calf Size C in Black

182CC99

SEA ISLAND COTTON

for women & men



Compression Level

15-20mmHg

Also available in 20-30mmHg, see page 22

Features and Benefits

- » Knitted from the world's rarest cotton
- » Sea Island Cotton yarns are lightweight and have great strength, a silky luster, and an incredibly soft feel
- » Breathable and comfortable for sensitive skin
- » Ideal for travel and everyday wear

Fibre Content

- » 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- » Latex free

Best Suited for


- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease (CVD)

Sizing Chart

151	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23cm (7"-9")	19-25cm (7.5"-10")	20-27cm (8"-10.5")
191	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	19-25cm (7.5"-10")	20-27cm (8"-10.5")	21-29cm (8.5"-11.5")
	Size	A	B	C

Style and Colours

Closed Toe = ●



Calf (C)	●	●	●
----------	---	---	---


Product Code Example: **151CA10**

Series	Style	Size	Colour
151	C	A	10



151 Women's Calf Style
15-20mmHg Closed Toe Calf Size A in Navy

151CA10



191 Men's Calf Style
15-20mmHg Closed Toe Calf Size B in Brown

191CB11

ALL-SEASON MERINO WOOL

for women & men

Compression Level

15-20mmHg

Also available in 20-30mmHg, see page 24

Features and Benefits

- » Soft and comfortable
- » Natural properties of wool make these socks odour resistant
- » Thermoregulating for year-round wear and keeps feet dry

Fibre Content

- » 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- » Latex free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease (CVD)




ZURICH
COLLECTION

Sizing Chart

152	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23cm (7"-9")	19-25cm (7.5"-10")	20-27cm (8"-10.5")
192	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	19-25cm (7.5"-10")	20-27cm (8"-10.5")	21-29cm (8.5"-11.5")
	Size	A	B	C

Style and Colours


Closed Toe = ●



Calf (C)	●	●	●
----------	---	---	---

Product Code Example: **192CB11**

Series	Style	Size	Colour
192	C	B	11



152
Women's Calf
Style

15-20mmHg
Closed Toe Calf Size A in Black

152CA99



192
Men's Calf
Style

15-20mmHg
Closed Toe Calf Size C in Navy

192CC10

MICROFIBER SHADES

for women & men

Compression Level

15-20mmHg

Also available in 20-30mmHg, see page 28

Features and Benefits

- » Advanced microfiber technology for ultimate comfort
- » Glides on and off with ease
- » Excellent breathability and moisture management
- » Soft and durable sock for everyday wear

Fibre Content

- » 87% Nylon, 13% Spandex
- » Latex free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease (CVD)



Sizing Chart

143	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23cm (7"-9")	19-25cm (7.5"-10")	20-27cm (8"-10.5")
183	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	19-25cm (7.5"-10")	20-27cm (8"-10.5")	21-29cm (8.5"-11.5")
Size		A	B	C

Style and Colours

Closed Toe = ● <i>Actual colour may vary</i>							
	Calf (C)	●	●	●	●	●	●
Closed Toe = ● <i>Actual colour may vary</i>							
	Calf (C)	●	●	●	●	●	●


Mariner Stripe, Grey Plaid, and Black Plaid available in 15-20mmHg only

* Available for women only

** Available for men only

Product Code Example: 143CB02


Series	Style	Size	Colour
143	C	B	02



143
Women's
Calf Style

15-20mmHg
Closed Toe Calf Size B in Pink Argyle

143CB45



183
Men's
Calf Style

15-20mmHg
Closed Toe Calf Size C in Black Plaid

183CC42

NEW!

SHEER FASHION

for women

Compression Level

15-20mmHg

Features and Benefits

- » Fabric feels soft against the skin
- » Incredible sheer design in a variety of fashionable colours
- » Sheer toe and heel area
- » Comfortable stretch is excellent for everyday wear
- » Cotton-lined gusset on pantyhose styles

Fibre Content

- » 79% Nylon, 21% Spandex
- » Content can vary slightly by style
- » Latex free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease (CVD)



15-20mmHg
Closed Toe Calf Size C in Charcoal

120CC12



15-20mmHg
Closed Toe Pantyhose Size C in Black

120PC99

Pantyhose and Maternity Sizing Chart

Ft/Cm	Kg/Lbs	45 / 100	48 / 105	50 / 110	52 / 115	54 / 120	57 / 125	59 / 130	61 / 135	64 / 140	66 / 145	68 / 150	70 / 155	73 / 160	75 / 165	77 / 170	80 / 175	82 / 180	84 / 185	86 / 190	89 / 195	91 / 200	93 / 205	95 / 210	98 / 215	100 / 220	102 / 225	105 / 230		
145cm / 4'9"																														
147cm / 4'10"			SIZE A																											
150cm / 4'11"																														
152cm / 5'0"			ANKLE CIRCUMFERENCE 19-23cm (7.5"-9")																											
155cm / 5'1"																														
157cm / 5'2"																														
160cm / 5'3"																														
162cm / 5'4"			SIZE B																											
165cm / 5'5"																														
167cm / 5'6"			ANKLE CIRCUMFERENCE 19-23cm (7.5"-9")																											
170cm / 5'7"																														
173cm / 5'8"																														
175cm / 5'9"																														
178cm / 5'10"																														
180cm / 5'11"																														
183cm / 6'0"																														

Note: please use pre-pregnancy weight for maternity pantyhose. * Maternity pantyhose not available in size F.

Calf / Thigh Sizing Chart

120	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	18-23cm (7"-9")	19-25cm (7.5"-10")	20-27cm (8"-10.5")
	Size	A	B	C

Weight Guidelines For Thigh Style: Size A: <135Lbs., Size B: <180Lbs., Size C: <230Lbs.

Styles and Colours

	CHARCOAL 12	TAUPE 29	HONEY 33	GOLDEN 36	CAFÉ 73	BLACK 99
Calf (C)	●	●●	●●	●●	●	●
Thigh (N)	●	●	●	●	●	●
Pantyhose (P)	●	●	●	●	●	●
Maternity Pantyhose (M)			●			●

Product Code Example: 120PE73

Series	Style	Size	Colour
120	P	E	73

Product Code Example: 120CCO36

Series	Style	Size	Open Toe	Colour
120	C	C	O	36

170

SOFT SILHOUETTE LEGGINGS

for women

Compression Level

15-20mmHg

Features and Benefits


- » Stylish, subtle ribbed design looks chic and pairs perfectly with dresses, blouses, tunics, and shirts
- » Foldable, comfort-stretch waist band provides a seamless look for all-day wearing comfort
- » Advanced soft stretch microfiber fabric shapes legs and smooths skin

Fibre Content

- » 83% Nylon, 17% Spandex
- » Latex free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Occupational leg health risks, family history of Chronic Venous Disease (CVD)



Non-binding ankle cuff
stays comfortably in place
for all day wear

15-20mmHg
Legging Size B in Black

170LB99

Sizing Chart


cm / ft	Kg / Lbs	45 / 100	48 / 105	50 / 110	52 / 115	54 / 120	57 / 125	59 / 130	61 / 135	64 / 140	66 / 145	68 / 150	70 / 155	73 / 160	75 / 165	77 / 170	80 / 175	82 / 180	84 / 185	86 / 190	89 / 195	91 / 200	93 / 205	95 / 210	98 / 215	100 / 220	102 / 225	105 / 230	
147cm / 4'10"																													
150cm / 4'11"																													
152cm / 5'0"																													
155cm / 5'1"																													
157cm / 5'2"																													
160cm / 5'3"																													
162cm / 5'4"																													
165cm / 5'5"																													
167cm / 5'6"																													
170cm / 5'7"																													
173cm / 5'8"																													
175cm / 5'9"																													
178cm / 5'10"																													
180cm / 5'11"																													
183cm / 6'0"																													

Guidelines for Ankle Circumference

Size A: 18–23cm (7"–9"), **Size B & C:** 20–26cm (8"–10"), **Size D:** 22–29cm (8.5"–11.5")

Style and Colours

Open Foot = ●



Leggings (L)



Product Code Example: 170LB99

Series	Style	Size	Colour
170	L	B	99

ATHLETIC RECOVERY SOCKS

for women & men

Compression Level

15-20mmHg

Features and Benefits

- » Helps reduce lactic acid buildup
- » Comfortable and easy to put on and take off
- » *drirelease* yarns provide moisture wicking and odour control
- » No compression in heel and instep for ease of donning

Fibre Content

- » 67% *drirelease* (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- » Latex free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Post-workout recovery
- » Mild athletic activity, walking




Sizing Chart

401	Women's Shoe	<7.5	7.5-9.5	10-12	12.5-15	≥15.5
	Men's Shoe	<6.5	6.5-8.5	9-11	11.5-14	14.5-17
	Ankle Circumference	18-24cm (7"-9.5")	19-26cm (7.5"-10")	20-27cm (8"-10.5")	21-29cm (8.5"-11.5")	21-29cm (8.5"-11.5")
	Size	S	M	L	X	2

Style and Colours

Closed Toe = ●



Calf (C)



Product Code Example: **401C200**

Series	Style	Size	Colour
401	C	2	00



15-20mmHg
Closed Toe Calf Size X-Large in Black

401CX99



15-20mmHg
Closed Toe Calf Size Medium in White

401CM00

412

PERFORMANCE SLEEVES

for women & men

Compression Level

20-30mmHg

Features and Benefits

- » Fabric resists deterioration in salt, chlorine, or fresh water
- » Stabilizes and supports active muscles
- » Quick drying
- » High-visibility reflective transfer

Fibre Content

- » 75% Nylon, 25% Spandex
- » Latex free

Best Suited for

- » Running, triathlons, endurance sports, and gym workouts
- » Team sports



Sizing Chart					
412	Ankle Circumference	18–21.5cm (7"–8.5")	21.5–25.5cm (8.5"–10")	25.5–29cm (10"–11.5")	29–33cm (11.5"–13")
	Calf Circumference	28–39cm (11"–15.5")	31.5–45cm (12.5"–17.5")	35.5–51cm (14"–20")	40.5–58cm (16"–23")
	Size	S	M	L	X

Style and Colours

Calf Sleeve = ●	WHITE 00	BLUE 50	PINK 56*	PURPLE 58*	STEEL BLUE 62	BLACK 99
Calf Sleeve (V)	○	●	○	●	●	●

* Pink and purple not available in X-Large

Product Code Example: 412VL56

Series	Style	Size	Colour
412	V	L	56



20–30mmHg
Calf Sleeve Size Small in Purple

412VS58



20–30mmHg
Calf Sleeve Size Medium in Steel Blue

412VM62

MERINO OUTDOOR SOCKS

for women & men

Compression Level

15–20mmHg

Also available in 20–30mmHg, see page 36

Features and Benefits

- » Soft and comfortable
- » Natural properties of wool make these socks odour resistant
- » Mesh flex zones for comfort and ventilation
- » Cushioned, padded soles protect against blisters
- » No compression in heel and instep for ease of donning
- » Thermoregulating for year-round comfort

Fibre Content

- » 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- » Latex free

Best Suited for

- » Tired, heavy, restless, aching legs
- » Outdoor activities




Sizing Chart

421	Women's Shoe	<7.5	7.5–9.5	10–12	12.5–15
	Men's Shoe	<6.5	6.5–8.5	9–11	11.5–14
	Ankle Circumference	18–24cm (7"–9.5")	19–26cm (7.5"–10")	20–27cm (8"–10.5")	21–29cm (8.5"–11.5")
	Size	S	M	L	X

Style and Colours

Closed Toe = ●



Calf (C)	●	●
----------	---	---

Product Code Example: 421CS31

Series	Style	Size	Colour
421	C	S	31



15-20mmHg
Closed Toe Calf Size Large in Olive

421CL31



15-20mmHg
Closed Toe Calf Size Medium in Charcoal

421CM12

EVERSOFT DIABETIC SOCKS

for women & men

Compression Level

8-15mmHg

Features and Benefits

- » Thick, padded soles, heels, and toes for comfort and protection
- » Flat, low friction toe seam
- » Soft, stretchable knit-in band prevents pinching
- » Compression from ankle to knee, compression-free foot

Fibre Content

- » 67% *drirelease* (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- » Latex free

Best Suited for


- » Individuals with diabetes
- » Individuals with limited dexterity or lower tolerance for compression
- » Individuals with sensitive skin



Sizing Chart					
160	Women's Shoe	<7.5	7.5–9.5	10–12	12.5–15
	Men's Shoe	<6.5	6.5–8.5	9–11	11.5–14
	Ankle Circumference	18–24cm (7"–9.5")	19–26cm (7.5"–10")	20–27cm (8"–10.5")	21–29cm (8.5"–11.5")
	Size	S	M	L	X

Style and Colours

Closed Toe = ●



Calf (C)



Product Code Example: 160CS00

Series	Style	Size	Colour
160	C	S	00



8–15mmHg
Closed Toe Calf Size Small in White

160CS00



8–15mmHg
Closed Toe Calf Size Large in White

160CL00

sigvaris



Medically Complex Edema & Inelastic Wraps

MEDICALLY COMPLEX EDEMA & INELASTIC WRAPS

The SIGVARIS GROUP line of products for medically complex edema (MCE) are for the treatment and management of conditions impacting the lymphatic and venous circulatory systems.

The product portfolio contains ready-to-wear garments for the active treatment and maintenance of lymphedema and for patients managing chronic venous disorders in the lower extremities.

Sigvaris is looking towards a brighter future with new therapies for complex venous and lymphatic diseases. Our inelastic wraps deliver maximum therapeutic results with increased patient compliance.

Why Choose Inelastic Wraps?

- Adjustable pressure
- Large variety of solutions

Medically Complex Edema & Inelastic Wraps


Product	Page
Compreflex Standard Calf & Foot	76
Compreflex Standard Calf	76
Compreflex Transition Calf	78
Compreflex Complete Calf	78
Compreflex Strap Extender	80
Assorted Liners	81
Assorted Foot Garments	82

Compreflex® Standard CALF & FOOT

20-50mmHg



Colours



When to Use

- Adjustable ranges of resting compression is desired (20-50mmHg)
- Non-attached foot piece is desired
- No compression in the liner is desired
- Patient requires a calf and foot garment

Conditions

- Lymphedema
- Edema
- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Lipedema



Compreboot® Standard

Variable compression
(see page 82 for details)




Basic Liner

Non-Compressive
(see page 81 for details)
One pair of Liners
included with each leg



Accutab® Devices

Accutab Devices
(see page 80 for details)



Compreflex Standard CALF

20-50mmHg



Colours



When to Use

- Adjustable ranges of resting compression is desired (20-50mmHg)
- No compression in the liner is desired
- No Compreboot needed
- Patient requires compression in the calf only or a different foot option
- Ability to wear standard footwear over liner is desired

Conditions

- Lymphedema
- Edema
- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Lipedema



Basic Liner

Non-Compressive
(see page 81 for details)
One pair of Liners
included with each leg



Accutab Devices

Accutab Devices
(see page 80 for details)




Medically Complex Edema & Inelastic Wraps

COMPREFLEX WRAPS



Variable compression Compreboot for added foot and ankle compression (see page 82 for details)

REQUIRED MEASUREMENTS



Measuring Instructions

Foot Size:

- I:** Measure circumference of ankle bend and heel
- J:** Measure circumference across metatarsal heads
- K:** Measure **medial** length from heel to 1st metatarsal head

Available Styles and Colours

	Black	Beige	Heathered Silver
	(99)	(77)	(N/A)
Compreflex Standard – Calf & Foot	●	●	
Compreflex Standard – Calf	●	●	
Basic Liner			●

Compreflex Standard Calf & Foot Sizing Chart

	Small	Medium	Large	X Large	XX Large	
C1	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm	
C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm	
B	24–34cm	29–39cm	34–44cm	39–49cm	44–55cm	
A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm	
Black	Regular	1401-BKR	1402-BKR	1403-BKR	1404-BKR	1405-BKR
	Tall	1401-BKT	1402-BKT	1403-BKT	1404-BKT	1405-BKT
Beige	Regular	1411-BKR	1412-BKR	1413-BKR	1414-BKR	1415-BKR
	Tall	1411-BKT	1412-BKT	1413-BKT	1414-BKT	1415-BKT

Product includes one Compreflex Standard Calf wrap, one Compreboot Standard Foot wrap, and 1 pair of Basic Liners.

Compreflex Standard – Calf Sizing Chart (A-C1)

	Small	Medium	Large	X Large	XX Large	
C1	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm	
C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm	
B	24–34cm	29–39cm	34–44cm	39–49cm	44–55cm	
A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm	
Black	Regular	1401-NF-BKR	1402-NF-BKR	1403-NF-BKR	1404-NF-BKR	1405-NF-BKR
	Tall	1401-NF-BKT	1402-NF-BKT	1403-NF-BKT	1404-NF-BKT	1405-NF-BKT
Beige	Regular	1411-NF-BKR	1412-NF-BKR	1413-NF-BKR	1414-NF-BKR	1415-NF-BKR
	Tall	1411-NF-BKT	1412-NF-BKT	1413-NF-BKT	1414-NF-BKT	1415-NF-BKT

Compreflex Transition

CALF

20-50mmHg



Colours



Black



Beige

When to Use

- Adjustable ranges of resting compression is desired (20-50mmHg)
- Patient requires a compressive liner; the Transition Liner provides 10-15mmHg in the foot and ankle
- Mild foot compression in liner is desired
- Ability to wear standard footwear over liner is desired

Conditions

- Lymphedema
- Edema
- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Lipedema



Transition Liner

10-15 mmHg of mild compression in the foot and ankle
(see page 81 for details)
One pair of Liners included with each leg



Accutab Devices

Accutab Devices
(see page 80 for details)

Compreflex Complete

CALF

20-50mmHg



Colours



Black



Beige

When to Use

- Adjustable ranges of resting compression is desired (20-50mmHg)
- Patient requires a compressive liner; the Complete Liner provides 10-15mmHg in the foot, ankle, and calf
- Mild foot and calf compression in liner is desired
- Ability to wear standard footwear over liner is desired


Conditions

- Lymphedema
- Edema
- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Lipedema



Complete Liner


Mild compression in foot, ankle, and calf
(see page 81 for details)
One pair of Liners included with each leg



Accutab Devices

Accutab Devices
(see page 80 for details)

REQUIRED MEASUREMENTS



Compreflex Complete Calf
20-50mmHg
Medium Regular in Beige
1412-FC-BKR

Available Styles and Colours

	Black	Beige	Heathered Silver
	(99)	(77)	(N/A)
Compreflex Transition – Calf	●	●	
Compreflex Complete – Calf	●	●	
Basic Liner			●

Liner included will match the colour of the wrap

Compreflex Transition – Calf Sizing Chart (A-C1)

	Small	Medium	Large	X Large	XX Large
C1	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
C	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
B	24-34cm	29-39cm	34-44cm	39-49cm	44-55cm
A	19-26cm	21-30cm	26-36cm	31-41cm	36-46cm

Black	Regular	1401-UC-BKR	1402-UC-BKR	1403-UC-BKR	1404-UC-BKR	1405-UC-BKR
	Tall	1401-UC-BKT	1402-UC-BKT	1403-UC-BKT	1404-UC-BKT	1405-UC-BKT
Beige	Regular	1411-UC-BKR	1412-UC-BKR	1413-UC-BKR	1414-UC-BKR	1415-UC-BKR
	Tall	1411-UC-BKT	1412-UC-BKT	1413-UC-BKT	1414-UC-BKT	1415-UC-BKT

Compreflex Complete – Calf Sizing Chart (A-C1)

	Small	Medium	Large	X Large	XX Large
C1	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
C	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
B	24-34cm	29-39cm	34-44cm	39-49cm	44-55cm
A	19-26cm	21-30cm	26-36cm	31-41cm	36-46cm

Black	Regular	1401-FC-BKR	1402-FC-BKR	1403-FC-BKR	1404-FC-BKR	1405-FC-BKR
	Tall	1401-FC-BKT	1402-FC-BKT	1403-FC-BKT	1404-FC-BKT	1405-FC-BKT
Beige	Regular	1411-FC-BKR	1412-FC-BKR	1413-FC-BKR	1414-FC-BKR	1415-FC-BKR
	Tall	1411-FC-BKT	1412-FC-BKT	1413-FC-BKT	1414-FC-BKT	1415-FC-BKT


Accutab Devices



Easy To Use

- Gently lay the Compreflex straps in place without stretching the fabric
- Attach one Accutab at the end of each strap
- Pull each strap to the prescribed dosage indicated by the lines on the Accutab: 20-30, 30-40, or 40-50mmHg

5 Accutabs are shipped for Regular Length and 6 Accutabs are shipped for Tall Length products



Compreflex Strap Extender

- Provides up to 10cm of additional circumferential size range to Compreflex calf wraps
- Detaches easily
- Potential solution for patients outside of standard size ranges of Compreflex calf wraps or needing custom size

	1 Unit	10 Unit Pack
Black	1401-SE1	1401-SEB
Beige	1411-SE1	1411-SEB



Basic Liner

Non-compressive

Features

- Non-compressive liner for the protection of the skin
- Easy to apply knee-high liners
- *drive/lease* moisture wicking and odour control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Heathered Silver colour only

Size: One size fits all	
Pair	10 Pair Pack
1700 - NK1	1700 - NKC
Included With Compreflex Standard Calf & Foot Compreflex Standard Calf	



Transition Liner

10-15mmHg

Features

- Mild compression in foot and ankle only
- *drive/lease* moisture wicking and odour control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Black and Beige

Transition Liner Sizing Chart (Included with Compreflex Transition)

SIZE	Ankle Circ. (A)	Calf Circ. (C)	Women's Shoe	Men's Shoe	Pair (Black)	Pair (Beige)
Small	19-26cm	29-39cm	<7½	<6½	1601-CK	1611-CK
Medium	21-30cm	34-44cm	7.5-9.5	6.5-8.5	1602-CK	1612-CK
Large	26-36cm	40-50cm	10-12	9-11	1603-CK	1613-CK
X Large	31-41cm	46-56cm	12.5-15	11.5-14	1604-CK	1614-CK
XX Large	36-46cm	54-64cm	12.5-15	11.5-14	1605-CK	1615-CK



Complete Liner

10-15mmHg

Features

- Mild compression in foot and ankle with graduated compression in calf
- Comfortable knee-high liners
- *drive/lease* moisture wicking and odour control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Black and Beige

Complete Liner Sizing Chart (Included with Compreflex Complete)

SIZE	Ankle Circ. (A)	Calf Circ. (C)	Women's Shoe	Men's Shoe	Pair (Black)	Pair (Beige)
Small	19-26cm	29-39cm	<7½	<6½	1601-CK-F1	1611-CK-F1
Medium	21-30cm	34-44cm	7.5-9.5	6.5-8.5	1602-CK-F1	1612-CK-F1
Large	26-36cm	40-50cm	10-12	9-11	1603-CK-F1	1613-CK-F1
X Large	31-41cm	46-56cm	12.5-15	11.5-14	1604-CK-F1	1614-CK-F1
XX Large	36-46cm	54-64cm	12.5-15	11.5-14	1605-CK-F1	1615-CK-F1



Compreboot Lite Foot

Variable Compression



Features

- Ultra-thin profile that easily fits into most footwear
- Highly adjustable design may be trimmed to fit
- Provides mild, therapeutic compression

Colour



Black



Compreboot Standard Foot

Variable Compression



Features

- Four-way stretch for firm compression around ankle
- Ideal for mild to moderate chronic edema
- Low stretch material conforms to foot shape

Colours



Black



Beige



Compreboot Plus Foot

Variable Compression



Features

- Chip foam lining creates high / low pressure ridges
- Removable gripper patches provide traction
- Provides moderate, therapeutic compression

Colour



Black



Circumference


Left Right

I _____

J _____

Length

K _____



Measuring Instructions

Foot Size:

I: Measure circumference of ankle bend and heel

J: Measure circumference across metatarsal heads

K: Measure **medial** length from heel to 1st metatarsal head

Compreboot Lite Foot Size Chart (I, J, K)

(Available in Black Only)

	Small / Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long
I	28–36cm	28–36cm	≤44cm	≤44cm
J	22–26cm	22–26cm	≤33cm	≤33cm
K	14–18cm	19–23cm	20–22cm	25–27cm
Black	1803-BLR	1803-BLL	1805-BLR	1805-BLL

Compreboot Standard Foot Size Chart (I, J, K)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28–36cm	28–36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22–26cm	22–26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14–18cm	19–23cm	18–20cm	23–25cm	20–22cm	25–27cm
Black	1801-BTR	1801-BTL	1803-BTR	1803-BTL	1805-BTR	1805-BTL
Beige	1811-BTR	1811-BTL	1813-BTR	1813-BTL	1815-BTR	1815-BTL

Compreboot Plus Foot Size Chart (I, J, K)

(Available in Black Only)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28–36cm	28–34cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22–26cm	22–25cm	≤30cm	≤28cm	≤33cm	≤33cm
K	14–18cm	19–23cm	18–20cm	23–25cm	20–22cm	25–27cm
Black	1801-BPR	1801-BPL	1803-BPR	1803-BPL	1805-BPR	1805-BPL

SIGVARIS
GROUP



ACCESSORIES

Help Yourself

Why choose **Accessories**?

- Aid in putting on and removing stockings
- Stocking cleaning to maximize usage
- Extra liners for added comfort

Accessories

Product	Page
Gloves	86
Sim-Slide	86
Cotton Underliner	86
Magnide On/Off	87
Easy Slide Arm	87
Washing Solution	87
Foot Slip	88
Sigvaris Doff N' Donner	88
Sigvaris Doff N' Donner Cone	88
MELANY	88

Gloves

Rubber Gloves* – Ridged
The practical and protective way to apply and remove your compression stockings.



Package: 12 per Pack

Size	Product Code
Small	591R400S
Medium	591R400M
Large	591R400L
X-Large	591R400X
ASST'D 4x S/M/L	591R400A

Latex-Free Gloves
Gripping donning gloves recommended for use in a latex-free environment.



Package: 12 per Pack

Size	Product Code
Small	592R400S
Medium	592R400M
Large	592R400L
X-Large	592R400X



Textile Gloves* – These high-quality gloves have a latex surface to provide a good grip when putting on and taking off compression stockings.

Package: 12 per Pack

Size	Product Code
X-Small	593R400XS
Small	593R400S
Medium	593R400M
Large	593R400L
X-Large	593R400X
ASST'D 4x XS/S/M	593R400A



Cotton Underliner

10–15mmHg stocking. Can be used as an underlayer to protect fragile skin, secure wound dressings, add 10–15mmHg of compression, or create an easy glide surface to aid in putting on and taking off compression stockings. For sizing, please refer to the Natural Rubber Calf sizing on page 45.

Package: 2 per pack

Size	Product Code	Natural Rubber Sizing
Small / Short	588USS	S1/S3
Small / Long	588USL	S2/S4
Medium / Short	588UMS	M1/M3
Medium / Long	588UML	M2/M4
Large-XL / Short	588ULS	L1/L3 – X1/X3
Large-XL / Long	588ULL	L2/L4 – X2/X4

Rubber Glove and Latex-Free Glove Sizing Chart			
Measure Circumference of the Hand at the Palm			
16.5–18cm (6.5"–7")	19–20.5cm (7.5"–8")	21.5–23cm (8.5"–9")	24–25.5cm (9.5"–10")
Small	Medium	Large	X-Large

Textile Glove Sizing Chart					
Measure Length and Width of the Hand					
Length	20.5–21cm	21.5–22cm	22.5–23cm	23.5–24cm	24.5–25cm
Width	7–7.5cm	7.5–8cm	8–8.5cm	9–9.5cm	10–10.5cm
	X-Small	Small	Medium	Large	X-Large

* **Caution:** Rubber Gloves and Textile Gloves contain natural rubber latex that may cause allergic reaction



For more **SIGVARIS** videos, including donning tips and product demonstrations, visit our YouTube channel:
SIGVARIS GROUP Canada

Donning Aids



Sim-Slide

Size	Product Code
Small	587SLIDESM
Medium	587SLIDEMD
Large	587SLIDELG
X-Large	587SLIDEXL

Sizing Chart

For use with Open toe styles

Women's Shoe Size	≤4-6	6.5-8	8.5-9.5	≥10
Men's Shoe Size	≤3-4.5	5-6.5	7-8	≥8.5
Size	S	M	L	XL



Magnide On/Off

Size	Product Code
Medium	587MAGMD
Large	587MAGLG
X-Large	587MAGXL

Sizing Chart

For use with Closed toe styles

Women's Shoe Size	≤8.5	9-12	≥12.5
Men's Shoe Size	≤7.5	8-11	≥11.5
Size	M	L	XL*

* XL suitable for extra-wide feet



Easy Slide Arm

Application aid for armsleeves.

Size	Product Code
Small/Medium	587E01408
Large	587E01425

Washing Solution



Washing Solution

Gentle washing solution specifically formulated to help maintain the quality of SIGVARIS GROUP compression socks and stockings. Also great for washing delicate garments!

Package Options

Contains	Product Code
Six 16 oz. Bottles (473ml)	586W900
One 16 oz. Bottle (473ml)	586W950
Twelve 2 oz. Bottles (60ml)	586W975



Foot Slip

Comfortable applicator for use with open toe stockings. Each pair of Open Toe stockings ships with one complimentary Foot Slip.

Package: 6 per pack

Size	Product Code
One Size	580F200

Donning Aids



SIGVARIS Doff N' Donner

A soft device that allows you to quickly guide a compression sock or stocking up and down your leg or arm, making it simple to put on and remove all styles of graduated compression. It is perfect for both self-donning and caregiver use. Latex free.

Size	Product Code
One Size	583D100



SIGVARIS Doff N' Donner Cone

This loading device makes it easy to roll stockings onto the SIGVARIS Doff N' Donner. The bottom of the base firmly attaches to most clean, flat surfaces via suction.

Size	Product Code
One Size	583D200



MELANY

The new MELANY from Sigvaris makes both donning and doffing compression garments simple. Its innovative one-piece design with adjustable handles was made especially for users with limited range of motion and finger dexterity.

Size	Calf Circumference	Product code
X-Small	25–30cm	579SXS
Small	30–38cm	579SSM
Medium	38–46cm	579SMD
Large	46–56cm	579SLG



Symptoms & Indications

Symptoms / Indications	Compression Level at the Ankle
<ul style="list-style-type: none"> ▪ Tired, achy legs ▪ Occupational or evening edema ▪ Leg discomfort from long hours of standing or sitting ▪ Prevent varicosities & swelling during pregnancy ▪ Preventative for legs predisposed to risk ▪ Reduce swelling during travel ▪ Prevent DVT during travel 	15–20mmHg
<ul style="list-style-type: none"> ▪ Heavy, fatigued, aching legs ▪ Mild edema in lower extremities ▪ Mild varicosities with minimal edema ▪ Mild varicosities during pregnancy ▪ Lymphedema upper extremities ▪ Post-procedure of small veins ▪ Prevention or management of DVT <p><i>Note: A prescription is recommended for stockings with the indications of 20mmHg or higher</i></p>	20–30mmHg
<ul style="list-style-type: none"> ▪ Moderate varicosities with mild to moderate edema (also during pregnancy) ▪ Post-procedure of larger veins to maintain treatment success ▪ CEAP C3 (edema), C4 (skin changes without ulceration), C5 (skin changes with healed ulceration), C6 (skin changes with active ulcerations) ▪ Management and treatment of DVT or Post-Thrombotic Syndrome ▪ Superficial phlebitis ▪ Lymphedema after decongestant therapy to maintain reduction ▪ Orthostatic / postural hypotension ▪ Post-fracture, post-traumatic edema 	30–40mmHg
<ul style="list-style-type: none"> ▪ Lymphedema after decongestant therapy to maintain reduction ▪ Severe Post-Thrombotic Syndrome ▪ Severe varicosities and/or edema 	40mmHg or higher
<ul style="list-style-type: none"> ▪ Severe arterial insufficiency, intermittent claudication, ischemia ▪ Uncontrolled congestive heart failure 	Contraindications
<ul style="list-style-type: none"> ▪ Skin sensitivities or allergies ▪ Advanced neuropathy with sensory loss in leg or foot ▪ Diabetes with arterial compromise ▪ Confinement to bed or non-ambulatory use unless otherwise prescribed by the physician ▪ Eczema, hypodermatitis, dermatosclerosis, dermatitis 	Cautions
<p>The therapeutic effect of graduated compression stockings is based on the compression level. With the proper CEAP classification, the correct compression level at the ankle can be chosen. For more severe symptoms, consult with a vascular specialist.</p>	
<p><i>No liability accepted for non-observance of contraindications and cautions.</i> Fronck HS, Bergan JJ, et al. The Fundamentals of Phlebology: Venous Disease for Clinicians. 2004. pg 151. Blaettner, et al. Leg Symptoms of Healthy People and Their Treatment with Compression Hosiery. Phlebology 2008. Compression therapy of the extremities. Partsch, Rabe, Stemmer. ISBN: 2.85480.770.7.</p>	

CoverOne

By
SIGVARIS
GROUP

REUSABLE FACE MASK

Comfortable • Durable • Breathable

Double-layered with earloops & pocket
Made with soft washable microfiber


Helping people feel their best is our priority

Sigvaris has been a trusted healthcare solutions partner for over 150 years. While we continue to make medical compression socks and hosiery, we have rapidly adapted our production facility in Peachtree City, Georgia to also produce cloth face masks. We're pleased to be able to help people follow healthcare guidelines with CoverOne, a high-quality cloth face mask.



Features and Benefits

- Constructed of two layers of soft microfiber
- Earloops made with microfiber yarns fit comfortably on the ears
- Features a pocket to accommodate an optional disposable filter (available to order)
- Adding filter makes CoverOne a three-layer mask
- Knit-in nose zone minimizes mask sliding
- Text embossed on the mask to indicate which side goes against the face and top vs. bottom orientation
- Hand or machine washable (Sigvaris Washing Solution may be used). May tumble dry on low heat or delicate setting. Do not bleach or iron
- Unisex, One size fits most adults
- Available in Black, Blush Pink, Coyote Brown, Denim, and Desert Yellow
- Fibre Contents: 89% Nylon, 11% Spandex. Latex free
- Manufactured in an ISO 13485 manufacturing facility
- Patent pending




Series	Description	Size	Colours
M901	Mask with Pocket	1 - Adult Regular	Black (199) Blush Pink (202) Coyote Brown (186) Denim (201) Desert Yellow (203)

Product Code Example: **M901199**

Series	Size	Colour	Quantity
M901	1	199	One mask

This product is non-returnable and cannot be exchanged due to health and hygiene reasons



5 Layers PM2.5 Activated Carbon Filter

Series	Description
M905100	PM2.5 Activated Carbon Filter Single 10-pack

Disclaimers / Cautions

- This is a non-medical grade cloth face mask and is not intended for medical use.
- Cloth face coverings should not be placed on young children under age 2, anyone who has trouble breathing, or is unconscious, incapacitated or otherwise unable to remove the mask without assistance.
- The manufacturer is not liable for damage/injuries which have been caused by improper handling or misuse.

Guidance for Canada / Recommandations pour le Canada

The mask must not be shared with other individuals and it does not provide complete protection against contraction of the virus. This mask may not be effective in blocking virus particles that may be transmitted by coughing, sneezing or certain medical procedures.

This product is non-returnable and cannot be exchanged due to health and hygiene reasons. As customer satisfaction is a high priority, please contact us with any questions about your order.

Term	Definition	Sigvaris Usage
CHRONIC VENOUS DISEASE (CVD)	A disease state comprised of a wide range of disorders including tired, achy legs, varicose veins, Chronic Venous Disease (CVD), skin issues, and venous leg ulcers (open and healed).	SIGVARIS GROUP manufactures MEDICAL and CONSUMER HEALTH CARE products to help reduce the risk of, improve symptoms of, and provide treatment for CVD.
DOUBLE-COVERED YARN	Spandex fibres are covered twice with yarn to prevent bare spandex from touching the skin for improved wearing comfort.	Used in nearly all SIGVARIS GROUP graduated compression socks and hosiery to create products that are easy to put on and take off, while also providing greater comfort.
DEEP VEIN THROMBOSIS (DVT)	A blood clot in a deep vein, most commonly in the lower leg or thigh.	SIGVARIS GROUP cares about DVT prevention. Clinical evidence supports the use of graduated compression to help prevent DVTs and the more serious complication of PE (Pulmonary Embolism), which can be life-threatening.
GRADUATED COMPRESSION VS. ANTI-EMBOLISM	Graduated compression is for ambulatory patients, whereas anti-embolism stockings are for non-ambulatory patients.	SIGVARIS GROUP offers true graduated compression, meaning products are tightest at the ankle and decrease in pressure going up the leg to help improve circulation.
HYPOALLERGENIC	Designed to minimize the chance of an allergic reaction.	SIGVARIS GROUP offers Latex free products and cotton fabrics for those with sensitive skin and allergies. Also, many SIGVARIS GROUP accessories are designed for consumers with sensitive skin.
ISO CERTIFIED	SIGVARIS GROUP is ISO (International Organization for Standardization) certified, indicating we have the required quality systems in place.	SIGVARIS GROUP follows ISO Standard 13485:2016 for our products manufactured in North and South America.
MICROFIBER	An extremely fine synthetic fibre knitted into textiles with the texture and drape of natural-fibre fabric. This fabric has excellent breathability.	Used to create products with a soft and comfortable feel. Find this fibre in 840 Soft Opaque, 830, 143 & 183 Microfiber Patterns, and 820 Microfiber.
mmHg (MILLIMETERS OF MERCURY)	Defined as the extra pressure generated by a column of mercury one millimeter high; the same unit of measurement used for blood pressure.	SIGVARIS GROUP products are offered in several ranges of compression levels based on this measurement.
MERINO WOOL	Merino is the softest and finest wool available in the world and is inherently itch-free.	SIGVARIS GROUP uses this premium wool yarn to create soft products that are thermoregulating for use in all seasons. Used in 422 Thermoregulating Wool, 421 MERINO OUTDOOR SOCKS, 240 Merino Wool, and 152 & 192 ALL-SEASON MERINO WOOL.
SEA ISLAND COTTON	Sea Island Cotton is the rarest cotton in the world and represents the finest in cotton softness.	SIGVARIS GROUP uses the highest-quality yarns in the industry to create products that are best in class. Look for this fine quality cotton in our 220 Sea Island and 151 & 191 SEA ISLAND COTTON SOCKS.
SENSINNOV	A patented SIGVARIS silicone top band that gently grips the skin and keeps the band in place all day.	This top band is featured in 780 Sheer, 840 Soft Opaque, 710 Patterns, and 820 Microfiber to help keep thigh-high stockings in place all day long.
SUPIMA COTTON SUPIMA WORLD'S FINEST COTTONS	Supima accounts for only about three percent of annual cotton production in the United States. Its fineness and longer staple length make it a premium cotton fibre.	SIGVARIS GROUP uses this premium cotton yarn in the 230 Cotton, 220 Sea Island and 151 & 191 SEA ISLAND COTTON SOCKS, 146 & 186 CASUAL COTTON, and 360, 142 & 182 Cushioned Cotton products, making them comfortable and ideal for those with sensitive skin.



**Sensitive
Skin**



Thermoregulating



Breathable



**Moisture
Wicking**



Soft




Durable



Cotton




Wool




Microfiber



**Natural
Rubber**



**Largest Size
Range**



**Wide Size
Range**



Sheer



**Medium
Sheer**



Opaque



Stylish



Stylish



Elegant



Colours



Patterns



**Sensinno
Grip-top
Band**



**Anatomically
Contoured**



**Cushioned
Sole**

Please follow these instructions carefully to determine the correct size for different product styles.

Measurements should always be taken in the following order, and as early in the morning as possible when swelling is minimal.

CALF STYLE

Open/Closed Toe

1	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM, LARGE or X-LARGE. Note: This should be your first measurement.
2	Measure the calf circumference at the fullest part of the calf.
3	Determine the calf length by measuring from the floor up to the bend of the knee. This measurement is taken on the outside of the leg.

PANTYHOSE AND THIGH-HIGH STYLE

Open/Closed Toe

1	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM or LARGE. Note: This should be your first measurement.
2	Measure the calf circumference at the fullest part of the calf.
3	Skip Step 3 (but refer to General Guidelines #6 below).
4	Measure the thigh at its widest circumference.
5	Determine the leg length by measuring from the floor up to the gluteal fold.
6	Determine the hip circumference by measuring the widest part of the hips.

GENERAL GUIDELINES


1	Measurements should always be taken as early in the morning as possible, when swelling is minimal.
2	Measure against bare skin, not over clothing (Steps 1-5).
3	The measuring system is based on having the foot flat against the floor and the calf or leg at a right angle to the foot (see images to right for examples).
4	Measurements should be taken by pulling the tape around the leg so that it is snug but not constricting.
5	Length measurements should be taken from the ground up.
6	Take all measurements, even if the patient only requires calf-length products, as the patient might decide at a later date to purchase additional styles.



1 Ankle circumference directly above ankle bone




2 Calf circumference at fullest part of the calf



3 Calf length from the floor to the right-angle bend of the knee



4 Widest circumference of the thigh



5 Leg length from the floor up to the gluteal fold



6 Circumference at the hip

How to Order Medical and Specialty Products

Series Number 23 = Cotton (See individual pages for series number)	Compression Level 1 = 15–20mmHg 2 = 20–30mmHg 3 = 30–40mmHg 4 = 40–50mmHg* 5 = 50–60mmHg*	Product Style C = Knee-high or Calf (See individual pages for available styles)	Size ML = Medium Long (See individual pages for sizing)	Gender M = Men W = Women O = Open toe	Colour 99 = Black (See individual pages for available colours)
23	2	C	ML	M	99

* 40–50mmHg and 50–60mmHg available in Natural Rubber only

To order Plus Calf with Grip-top or Plus Size Pantyhose (860P for women only):
add "/P" after the colour in the ordering code

To order Calf Style with Grip-top:
add "/S" after the colour in the ordering code (see product pages for availability)

To order Thigh with Waist Attachment:
add "/L" to indicate the Left leg or "/R" to indicate the Right leg after the colour in the ordering code (860W and 500W only)

To order Motion High Tech, Motion Thermoregulating Wool, Swiss Mountain Sock, and Canada Mountain Sock:
The letter for Gender is not used in the ordering code.
Example: 412CSM59 = Motion High Tech, 20–30mmHg, Closed toe Calf, Size Small Medium, Colour Red

How to Order CONSUMER HEALTH Products

Series Number 120 = SHEER FASHION (See individual pages for series number)	Product Style N = Thigh-high (See individual pages for available styles)	Size A = Small (See individual pages for sizing)	Colour 33 = Honey (See individual pages for available colours)
120	N	A	33

To order ATHLETIC RECOVERY SOCKS, PERFORMANCE SLEEVES, MERINO OUTDOOR SOCKS, and EVERSOFIT DIABETIC SOCKS:

Size is only one letter and the letter for Gender is not used in the ordering code.

Example: 421CL31 = MERINO OUTDOOR SOCK, 15–20mmHg, Closed toe Calf, Size Large, Colour Olive

To Order:

1.800.363-4999

canadaorders@sigvaris.com

Online Dealer Portal:


sigvaris-online.com

WHY SIGVARIS GROUP?

We produce the best compression garments to enhance our customers' health and quality of life.

We are dedicated to helping legs look and feel their best. Using premium yarns, we've spun ourselves into being the worldwide leader in graduated compression socks and hosiery by embracing our core beliefs of **QUALITY**, **COMFORT**, **FIT** and **INNOVATION**. Staying true to these values, we are continuing our efforts by expanding our business to help customers with medically-complex needs.

More than 20 people handle our products from start to finish, ensuring every garment is of the highest quality.



We offer the most comfortable fabric options available, including the finest cottons, wools and microfibers. We double-cover the yarns for ease of donning and doffing.

We offer premium products designed for all shapes and sizes, from petite to plus, including compression for a broad spectrum of edemas. All of our products are designed using the most recent human data measurement studies. A proper fit is key to adherence.



We use the finest textiles and materials to create modern styles in our circular-knit products, while maintaining medical efficacy. We continue to invest in medical research to help customers.



The Worldwide Leader
in Compression Therapy

Sigvaris, SIGVARIS GROUP, PFS², Sensinnov, Compreflex, Compreboot, and Accutab are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. © Copyright 2021 SIGVARIS, Inc.

Coolflex, Spacer Fabric TR, and CoverOne are trademarks of SIGVARIS AG, CH-9014, St. Gallen/Switzerland, in many countries worldwide.

Lycra is a registered trademark of INVISTA. Supima is a registered trademark of the Supima Association of America Corporation. NILIT is a registered trademark of Nilit. All rights reserved. *drilease* is copyright ©2021 *drilease*. *drilease* is an Optimer Brands company. All rights reserved. Sim-Slide and Magnide are registered trademarks of Arion Holding B.V.

Sigvaris Dealer Portal

SIGVARIS
GROUP

The SIGVARIS GROUP Dealer Portal is your online location for securely ordering SIGVARIS GROUP products and managing your account. Place and track your orders, review order history, contact Customer Care, and much more at sigvaris-online.com.

DEALER PORTAL BENEFITS

24/7 access – order any time!

Three ways to securely place an order online:
Catalog, Quick Entry or Order Upload

Access to order history and tracking

Email canadaorders@sigvaris.com
to get your account activated!



SIGVARIS GROUP Education



SIGVARIS GROUP offers 100% online educational seminars to help increase your leg health knowledge.

Become A SIGVARIS Certified Fitter!

Contact your Territory Manager

or call Customer Care at **1-800-363-4999**

or email canadaseminars@sigvaris.com

SIGVARIS
GROUP



ORDERING

Online: sigvaris-online.com

Call: 1-800-363.4999

Fax: 1-800-263.8736

Customer Service:

Monday–Friday, 8am–7pm (EST)



sigvarisgroup.canada



[sigvarisgroup_canada](https://www.instagram.com/sigvarisgroup_canada)



SIGVARIS GROUP Canada