

The logo features the word "alta" in a bold, dark blue, lowercase sans-serif font with a slight drop shadow. Below it, the word "MEDICAL" is written in a smaller, green, uppercase sans-serif font. The entire logo is centered within a light cream-colored circle. The background is a dark blue field with a repeating pattern of the "alta" logo in a lighter shade.

alta
MEDICAL

Why Alta Medical?



2024

Alta Medical & SiennaWrap Launch

2023

Patient and Therapist Boards

Designing, changing, and feedback process started

2022

Recommendations

Family friend asked Kathy and Kyle for a recommendation and they were surprised to see lots of familiarities.

2014


Sold Solaris

SiennaWrap

A fresh look on velcro Compression Wraps


Hard to don

Finger loops for easier use



Hard to contain

Bowed straps for less gaps



Options for Clinicians

Solutions for different presentations



Our Calf, Contoured

is designed for those who have a more pronounced calf with a smaller ankle.




Our Calf, Seamless

is designed for those who have a more cylindrical shape calf, with little to no difference in size of the ankle.


Hard to wear

Vapour-release limestone is more **breathable** & **durable** than petroleum.




Hard to fit real legs

Designed for limbs not mannequins



Hard to fit real legs

Reducing the risk of tourniquet



Hard to fit real feet



More options to fit unique needs

Foot



● Available in black and beige

Calf - Seamless



● Available in black and beige
Available in regular and tall height

Calf - Contoured



● Available in black and beige
Available in regular and tall height

Knee



● Available in beige

Thigh



● Available in beige



SiennaWrap
With Assist Sleeve



Full Leg
Separate Components



Hybrid Liner
Compression Hybrid Liner
8-15mmHg



Knee-high Liner
Non-Compressive



Thigh-high Liner
Non-Compressive



theofficialdesiretanner 21w
 Alta Medical launches in June!

I'm glad I was able to get my hands their velcro wrap before the launch.

I absolutely love how they added fingers loops which makes it easier patients to put on.

<https://altamedical.com>

♥ - velcro wraps for a day/ night tim garment.
 ♥ - velcro wraps because they may easier for patients to put on vs compression socks

6 likes
 May 30

alta Add a comment...



elastic compression
 AKA Velcro wraps for lymphedema/ lipedema

Reply to Leah's 3 comments
 You are my new favorite person.

TikTok
 @kathleenblissom

so let's take a look at inelastic compression or what we also call Velcro wraps

New Compression wrap has hit the market!
 @alta.medical

New to the market but not new to the compression world. Experienced compression family (m... See more

Original audio

stressreductionforlymphedema 19w
 Let's look at inelastic compression / compression wraps for #lymphedema and #lipedema - good options for clients that refuse or cannot tolerate stockings or decongest fast and need to adjust compression often. Not everyone can go directly from CDT to stockings. As always, consult your #lymphedematherapist to see what is right for you. These are #siennawraps from @alta.medical #lipedematreatment #lymphedematreatment

alariagodwin 19w
 Also helpful for those with lower extremity paralysis. Lymphedema

1 like Reply

30 likes
 June 17

alta Add a comment...



So forgive the amateur video but I just received the sample of the SiennaWrap and wanted to see how easy it is to use. So here you go.



“I ABSOLUTELY LOVE THE WRAP”

“REALLY LOVE THE DESIGN”

“A PATIENT CHOSE THEM OVER ANY OF THE OTHERS”

“LOVE THE OPEN FINGER HOLES”

“I WILL BE RECOMMENDING FOR MY FUTURE CLIENTS”

Adapt-a-Wrap

Affordable, Customizable Compression Garments

What is Adapt-a-Wrap?

Lower Cost

To follow-up on SiennaWrap, we focused on a developing a lower-cost garments for under-insured or cash pay patients.

Affordable & Adjustable

In the end, we developed a product that is not only more affordable, but also is adjustable and can better serve those in the active phase of treatment.

Customizable

However, as we designed the garment, we discovered that that we could make it customizable at no additional cost.



Built on Several of SiennaWrap's Favorite Features

Curved Design

To better conform to typical calf and promote fluid uptake

Billowed Straps

To reduce gaps and tissue pinching, improving comfort



Material

From same US-based manufacturer

Fingerloops

Eases donning and doffing process

USA Made

Proud to be American made

50% Thicker than SiennaWrap

Thicker


SiennaWrap foam is 2mm thick;
Adapt-a-Wrap foam is 3mm thick

Rigid

Provides a more rigid garment for
active phase reduction

Assistive

Additional foam may assist with fibrotic
tissue management



Easy to Trim-to-Fit

Straps

Strap length can be cut down as patient reduces volume

Billowing

The billowing can be reduced as well

Cut Markings

Easier to use 2cm markings



16 Unique Sizes

Variety of Sizing


Because we limited capital costs to manufacture, we can offer the widest variety of sizing of any wrap on the market

Variants

Four heights and four girths – and strap length is still customizable

ADAPT-A-WRAP CALF			
Cirth Measurements	D Circumference Two Finger Widths Below Crease at Back of Knee	C Calf Circumference at Widest Part	B Ankle Circumference at Narrowest Part
<i>Small</i>	29-41cm	29-41cm	19-26cm
<i>Medium</i>	34-46cm	34-46cm	21-31cm
<i>Large</i>	39-51cm	39-51cm	26-36cm
<i>Extra Large</i>	44-61cm	44-61cm	31-41cm

ADAPT-A-WRAP CALF	
Height Measurements	BD Two Finger Widths Below Knee Crease to Ankle Bony Prominence
<i>30</i>	30cm
<i>33</i>	33cm
<i>36</i>	36cm
<i>39</i>	39cm



Value

Less to Manufacture

The design is less for us to manufacture so we are passing that savings on to you and your patients

30% More Affordable

Development is not complete, but we think Adapt-a-Wrap will be ~30% more affordable than SiennaWrap




SiennaWrap

Flexible Fit System

Premium Trim-to-Fit Compression Wraps

What is the SiennaWrap FFS?

It is a 3 part modular garment that can be **customised** to fit a variety of presentation and can **reduce** as the limb decongests and can be **renewed**.




Built on the SiennaWrap Foundation?

Curved Spine

Finger Loops

Billowed
Straps



Breathable,
wickleable,
comfortable
fabric

Made in the
USA

What's Better?

Fully Independent Spine Changes Everything

- Easy to use
- Bilateral adjustment of the velcro to better contour to limb
- Adjustable from 100cm to 30cm
- All the cost is in the velcro straps, so the spine is replaceable at 20% the cost of the garment
- No-no velcro straps
- Spine is modular so you can pick the right curve for presentation

Indigo Nighttime

Foam Based - Therapeutic Garments

**COMING
IN 2025**